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IT

User manual CINCO MAS
10G Self loading stretcher with 5 heights

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EN

1. MODELS	15
2. INTENDED USE	15
3. REFERENCE STANDARD	15
4. INTRODUCTION	15
4.1 USE OF THE MANUAL	15
4.2 LABELLING AND TRACKING CONTROL OF THE DEVICE	15
4.3 SYMBOLS	16
4.4 WARRANTY AND SUPPORT	16
5. WARNINGS	16
6. SPECIFIC WARNINGS	17
6.1 PHYSICAL REQUIREMENTS OF OPERATORS	18
7. RESIDUAL RISK	18
8. TECHNICAL DATA AND COMPONENTS	19
9. INSTALLATION AND START-UP	20
10. FUNCTIONAL CHARACTERISTICS	20
11. INSTRUCTIONS FOR USE	21
11.1 MEDICAL VEHICLE REQUIREMENTS	21
11.2 STATIONING BRAKES	21
11.3 SIDEBARS	21
11.4 ACTIVATION OF TWIST SYSTEM	21
11.5 LOWERING THE STRETCHER	21
11.6 STRETCHER LIFTING	22
11.7 LOADING AND UNLOADING THE STRETCHER FROM THE MEDICAL VEHICLE	22
11.8 TRENDELENBURG PLATFORM	23
11.9 TRENDELENBURG/FOWLER PLATFORM	23
11.10BACKREST ADJUSTMENT	23
12. CLEANING AND MAINTENANCE	24

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1. MODELS

The standard following models can undergo change, revision and implementation without any notice.

- CINCO MAS 65
- CINCO MAS 75

The models above, are available with different colors of patient board (yellow, black, red), with Trendelenburg or Trendelenburg and Fowler and, depending on model, can include accessories such as the mattress.

2. INTENDED USE

Ambulance stretchers are the primary means of transporting sick and/or injured people in transit safely and comfortably in a lying position. It is not foreseen that the patient be able to intervene on the device, nor is it intended for prolonged storage or use as a hospital bed. The stretcher must be used with the dedicated fixing systems. It is not foreseen that the patient be able to intervene on the device.

TARGET PATIENTS

There are no particular indications related to the patient group.

The product configuration is able to accommodate any subject as long as he/she is within the maximum capacity of the device. If paediatric subjects must be transported, it will be in the role of the rescuer to determine whether the strap systems are suitable for immobilization or if it will be necessary to use another device.

PATIENT SELECTION CRITERIA

The intended patients are those for whom ambulance transport is necessary.

CONTRAINDICATIONS AND UNWANTED SIDE EFFECTS

No particular contraindications or side effects are known with relation to use of the device, as long as it is used in accordance with the user manual.

USERS AND INSTALLERS

The intended users are persons trained in first aid procedures and in use of medical equipment in an emergency medical service (EMS) environment.

Possible users also include emergency vehicle outfitters who may use the product before commissioning or during maintenance of the vehicle on which the stretcher is used.

USER TRAINING

Note: despite all efforts, laboratory tests, trials, and instructions for use, standards do not always reproduce practice, so the results obtained under actual conditions of product use in the natural environment may sometimes differ significantly.

The best instructions are the continuous practice of use under the supervision of competent and trained personnel.

- Regardless of your level of experience with similar devices in the past, you should carefully read and understand the contents of this manual before installing, operating, or servicing this product. In case of any questions, please contact Spencer Italia S.r.l. for the necessary clarifications.
- The product must be used only by personnel trained in the use of this product and not on other similar products.
- The suitability of the users for use of this product can be attested by the training registration, in which trained persons, trainers, date and place are specified. **This documentation must be kept for at least 10 years after the end of the product's life and must be made available to the competent authorities and/or the Manufacturer when requested. In the absence of such documentation, the relevant bodies will apply any foreseen sanctions.**
- Do not allow untrained persons to help while using the product, as they may cause injury to themselves or others.
- The product must be put into use only by personnel trained in the use of this product and not on other similar products.

Note: Spencer Italia S.r.l. is always available for training courses.

Document user training using the form in Annex A of this manual or other appropriate forms.

INSTALLER TRAINING

Device installation must be carried out by qualified personnel who have been trained and certified by Spencer Italia S.r.l.

The installer must strictly adhere to these instructions as well as to the state of the art for vehicle installations.

3. REFERENCE STANDARD

As a distributor or end user of products manufactured and/or marketed by Spencer Italia S.r.l. you are strictly required to have basic knowledge of any legal requirements applying to the devices contained in this supply that are in power in the goods final destination Country (including laws and norms regarding technical specifications and/or safety requirements) and therefore you are strictly required to have the necessary knowledge to guarantee all aspects regarding the total conformity of the products to the regulations in the relevant territory.

REFERENCE	TITLE OF DOCUMENT
EU Regulation 2017/745	EU Regulation on Medical Devices
UNI EN ISO 1865-1	Patient handling equipment used in road ambulances - Part 1: Specification for general stretcher systems and patient handling equipment
UNI EN 1789	Medical vehicles and their equipment - Road ambulances

4. INTRODUCTION

4.1 USE OF THE MANUAL

This manual is intended to provide the health care operator with the all the necessary information for its safe and appropriate use as well as adequate maintenance of the device.

Note: The manual is an integral part of the device. It must be kept for the duration of the device and must accompany the device in case of change of ownership or destination. If the operating instructions received relate to products not received, you must immediately contact the manufacturer before use.

The Spencer product manuals can be downloaded from the website or can be requested by <https://www.spencer.it/en/emergency-medical-solutions/user-manuals> or by contacting the manufacturer. Exceptions are items whose essentiality for reasonable and predictable use is such as to make it unnecessary to prepare instructions in addition to the following warnings and directions on the label.

Regardless of the level of experience gained in the past with similar devices, it is recommended that you carefully read this manual before installing, operating or using the product or any maintenance.

4.2 LABELLING AND TRACKING CONTROL OF THE DEVICE

Each device has got an identification label positioned on the device itself and/or on its box, which includes identification data about the manufacturer, the product, the CE mark, the serial number (SN) or lot number (LOT). **It must never be removed or covered.**

In case of damage or loss, request a duplicate from the manufacturer. Failure to do so will interrupt the validity of the guarantee as the device can no longer be traced.

The Regulation 2017/745/UE requires manufacturers and distributors of medical devices to keep track of the device location. If the device was in a different location to the address where it was sent or to where it had been sold, donated, lost, stolen, exported, or destroyed, permanently removed from use, or if the device had not been delivered directly from Spencer Italia S.r.l., inform the customer service (see § 4.4).

IT

EN

DE

EL

ES

HU

4.3 SYMBOLS

Symbol	Meaning	Symbol	Meaning
	Product complies with the requirements of Regulation 2017/745/EU		General or specific warnings
	Front leg unlocking		Rear leg unlocking / unlocking of the stretcher from the fastener
	Medical device		See instructions for use
	Manufacturer		Lot number
	Date of manufacture		Product code
	Unique Device Identifier		Serial Number

	Production identification Alphanumeric code that identifies the production units of the device, composed of: (01)0805771123 company prefix 000 progressive G51 6 control number (11)200626 date of production (YYMMDD) (10) 1234567890 SN
(01)8057711230006 (11) 200626 (10) 1234567890	

4.4 WARRANTY AND SUPPORT

Spencer Italia S.r.l. guarantees that products are without defects for a period of **one year from the date of purchase**.

For any information regarding the correct interpretation of the instruction manual, the use, maintenance, installation and restore of the product, contact the Spencer customer care service ph. +39.0521.541111, fax +39.0521.541222, e-mail service@spencer.it, or visit <https://www.spencer.it/en/medical-emergency-company> to find the nearest service centre for assistance.

In order to facilitate the assistance service, please always indicate the lot number (LOT) or serial number (SN) shown on the label applied on the box or on the device.

Conditions for warranty and assistance can be viewed on <https://www.spencer.it/en/support>

Note: Record and store; these instructions, lot (LOT) or serial number (SN) if any, date and place of purchase, date of first use, date servicing, user names and comments.

In order to manage the maintenance and service procedures of your equipment, Spencer provides you with the SPENCER SERVICE portal <https://www.spencer.it/en/services/after-sales> through which you can view the data of your products or their placing on the market, monitor and update the schedules of periodic revisions and read and manage the extraordinary maintenance operations.

5. WARNINGS

The warnings, notes and other important safety information are indicated in this section and are clearly visible throughout the entire manual.

At least every 6 months, is important to verify if updated User Manuals are available or there is any change involving your product. This information is freely available on the website <https://www.spencer.it/en/emergency-medical-solutions/user-manuals>.

User training

Note: Laboratory testing, post production tests, instruction manuals cannot always consider every possible scenario for use. This means that in some cases the performance of the product could be notable different from results to date obtained. **The best instructions are continuous use under the supervision of trained and competent staff.**

- Regardless of the level of experience gained previously with similar devices, it is recommended that you carefully read this manual before installing, operating or using the product or carrying out any maintenance. If in doubt, contact Spencer Italia S.r.l. to obtain the necessary clarifications.
- The product must be used by trained personnel only, having attended specific training for this device and not for similar products.
- The suitability of the user to use the product may be attested by the records of training, where the names of those trained, of the trainers, dates and place are indicated. **This register which certifies the eligibility of the operators to use the Spencer device must be kept for a period of 10 years after the disposal of the device itself. This register will be made available to the competent Authorities and/or manufacturer if requested. In the absence of such documentation, sanctions will be applied.**
- Do not allow any untrained person to help during the use of the device, because they could cause damage to the patient or to themselves.
- The product should be operated only by personnel trained in the use of this product and not in others similar.

Note: Spencer Italia S.r.l. is always at your disposal to organise product training.

Installers training

The installation of the device must be carried out only by qualified and trained staff authorized by Spencer Italia S.r.l. Timing and procedures for conducting training courses shall be agreed between the customer and our Commercial offices.

Product functionality

Use of the device in anyway other than described in this manual is forbidden.

- Before each use of device the perfect operating state of the device must be checked as specified in the instruction manual. If any damage or abnormalities which could in any way influence the correct functioning and the safety of the device are detected, the device must be immediately removed from service and the manufacturer must be contacted.
- If any failure or incorrect functioning of the device is detected, it must be immediately substituted with a similar item so that the rescue procedures are guaranteed without any interruption.
- The appliance must not in any way be tampered with (modification, adjustment, addition, replacement). In such cases all responsibility will be denied for any malfunctions or injuries caused by the appliance itself; moreover CE certification and product warranty will be considered void.
- Those who modify or have modified, prepare or have prepared medical appliances in such a way that they no longer serve the purpose for which they were intended, or no longer supply the intended service, must satisfy the valid conditions for the introduction onto the market.
- During use, position and adjust the device taking care not to cause any obstruction to rescuers and or any other rescue equipment.
- Ensure that all the necessary precautions are taken in order to avoid the hazards that can arise as the result of contact with blood or body fluids, when applicable.
- For devices intended for the transport of patients, always respect the maximum load capacity of the device, as indicated in this user's manual. Maximum load capacity means the total weight distributed according to the human anatomy. In determining the total loading weight of the product, the operator must consider the weight of the patient, the equipment and the accessories. Moreover, the operator must consider that the overall dimensions of the patient do not reduce the functionality of the device.
- For devices intended for the transport of patients, make sure, before lifting, that operators have appropriate physical conditions, as indicated in the manual.
- The maximum weight supported by each rescuer must comply with the requirements prescribed by the law of the land, in the field of Health and Safety at Work.**
- The warranty seals, where present, must not be removed;** in such case, the manufacturer will no longer recognize the product warranty and will accept no responsibility in case of incorrect operation or damage caused by the product itself.

- Avoid contact with sharp objects
- The installation of the device must be carried out only by qualified and trained staff authorized by Spencer Italia S.r.l. Timing and procedures for conducting training courses shall be agreed between the customer and our Commercial offices.
- Operating temperature: from -10°C to + 50°C

Storage

- The device should not be exposed to or come into contact with any source of combustion or inflammable agents. Store in a cool, dry, dark place and do not expose to direct sun.
- Do not store the device underneath any heavy objects which could cause structural damage.
- Store and transport the device in its original packaging. Failure to do so makes the warranty void.
- Storage temperature: from -20°C to +60°C

Maintenance/cleaning

Spencer Italia S.r.l. disclaims any liability for any damage, direct or indirect, which is a result of improper use of the product and replacement parts and / or otherwise of any repairs made by an entity other than the authorized Spencer service centres; this will also invalidate the warranty.

- The operator must always wear adequate personal protection such as gloves and mask etc. during all checking, maintenance and cleaning procedures.
- Establish a maintenance program and periodic testing, identifying an employee responsible for overseeing. The person to whom the ordinary maintenance of the device is entrusted must ensure the basic requirements foreseen by the manufacturer in the user's manual.
- **The frequency of inspection is determined by factors such as legal requirements, the type of use, frequency of use, environmental conditions during use and storage.**
- Spencer Italia S.r.l. disclaims any liability for any damage, direct or indirect, that is the result of incorrect repairs or use of products made by Spencer Italia S.r.l. Repairs must necessarily be carried out by an authorized Spencer Italia service centre, which in using original spare parts will provide a quality repair service in strict accordance with the technical specifications given by the manufacturer. Spencer Italia S.r.l. disclaims any liability for any damage, direct or indirect, which is a result of improper use of spare parts and/or otherwise of any repairs made by an entity other than the Spencer service centres authorized to repair or make substitutions on this product and parts and/or otherwise of any repairs made by an entity other than the Spencer service centres authorized to do so; the warranty will also be invalidated.
- Use only original components, spare parts and or accessories, approved by Spencer Italia S.r.l., in order to carry out any operation without causing any alteration or modification to the device.
- For any operations that are not carried out directly by the manufacturer but by an authorised centre, we have to underline that a report regarding all operations carried out must be requested. This will permit both Spencer Italia S.r.l. and the end user to keep a log book regarding the operations carried out on the device.
- All maintenance and revision must be recorded and documented with the corresponding report for technical assistance; documentation shall be maintained for at least 10 years from the end of life and must be made available to the competent authorities and/or the manufacturer if requested.
- The cleaning schedule for reusable products must be performed in accordance with the directions provided by the manufacturer in the user manual, in order to avoid the risk of cross-infection due to the presence of secretions and/or residuals.
- The device and all its components, after washing, should be allowed to dry completely before storing.
- If required, lubrication must be carried out after cleaning and complete drying.
- The fastener intended to be used with CINCO MAS stretcher, has an integrated system indicating the need for maintenance both for stretcher and fastener, related to the loading and unloading cycles. The required maintenance intervals shall be strictly respected.

Regulatory requirements

As a distributor or end user of products manufactured and/or marketed by Spencer Italia S.r.l., you are strictly required to have a basic knowledge of any legal requirements applying to the devices contained in this supply that are in power in the final destination Country (including laws and norms regarding technical specifications and/or safety requirements) of the goods and therefore you are also strictly required to have the necessary knowledge to guarantee all aspects regarding the total conformity of the products to the regulations in the relevant territory.

- Promptly notify Spencer Italia S.r.l. (already during the first product enquiry) when requesting in details regarding any revisions to be made by manufacturer in order to guarantee the conformity of the products to the territory's legal specifications (including those resulting from rules and/or norms of other kind).
- Act, with all due care and diligence, and contribute to ensure conformity to general safety requirements of all devices marketed in the territory, by providing final users with all necessary information for carrying out periodical checks on their devices, exactly as specified in the relevant user's manual.
- **Actively contribute to safety checks on product sold**, by communicating any relevant risk analysis information both to the manufacturer and to any competent authorities so that the necessary action can be promptly taken.
- The distributor or final user is aware that in the event of any failure to conform to the above mentioned requirements you will be deemed fully responsible for all damages that might occur. Therefore Spencer Italia S.r.l. expressly disclaims any responsibility and/or liability for your non-compliance with the present regulatory provisions.

General warnings for medical devices

The user must carefully read not only these general warnings, but also those listed below.

- It is not foreseen that the use of the device is prolonged beyond the time necessary for the first responders to the complete their operation and the subsequent stages of transport to the nearest rescue point.
- When the device is being used, the assistance of qualified staff must be guaranteed and at least one operator must be present.
- Follow the procedures and protocols approved by the internal organization.
- The activities of disinfection and sterilization should be carried out in accordance with the parameters given in the validated cycle as specified in the technical standards.
- With reference to Regulation 2017/745/UE, we remind both public and private operators, that in the exercise of their activity detect an accident involving a medical product are required to notify the Ministry of Health, under the terms and in the manner established by the relative ministerial decrees and also to the manufacturer. Health care providers whether public or private are required to communicate to the manufacturer, any other inconvenience that may allow for the adoption of measures to ensure the protection and health of patients and users.

6. SPECIFIC WARNINGS

For the use of CINCO MAS, the user must have read, understood and follow carefully all the instructions described in the user manual.

- Follow the procedures approved by the Emergency Medical Service for the immobilization and transport of patients.
- Follow the procedures approved by the Emergency Medical Service for the positioning and transport of patients.
- Do not use if the device or its parts are pierced, torn, frayed or excessively worn out.
- Make sure, before lifting, that the operators have a firm grip on the device.
- Avoid pulling the device on rough surfaces.
- Do not lift by crane or other mechanical lifters.
- Do not use drying machines.
- The device is intended for the transport of patients and **cannot be used as a parking device**.
- Do not use with devices other than those approved by the manufacturer.
- Have practice with an empty stretcher/ chair, in order to make sure you become familiar with the manoeuvres.
- For loading techniques of the patient, for particularly heavy patients, for working on uneven ground or in special and unusual circumstances, the presence of more operators is recommended (not only 2 as expected under standard conditions).
- Before placing the patient on the stretcher, make sure that the patient is adequately immobilized. Failure to immobilize correctly can cause serious damage.

IT

EN

DE

EL

ES

HU

- Make sure that the sheet does not interfere with any moving mechanism and control of the stretcher.
- Do not move the device if the weight is not evenly distributed.
- Always use belts anchored to the stretcher frame to ensure patient safety.
- Use only peripheral main frame for moving the stretcher and not the sidebars, tabletops or other parts not suitable for such purpose.
- Avoid extreme force during the loading of the stretcher on the ambulance. Too much force may have negative effect on the coupling system.
- Avoid extreme force during the loading of the stretcher on the ambulance. Too much force may have negative effect and can adversely affect the functionality of the stretcher.
- **Hold the stretcher firmly if the patient is on it.**
- **The stationing brakes are aids for the operator, they can not in any way substitute his supervision.**
- Pay a lot of attention to possible obstacles (water, ice, debris, etc.) on the route of the stretcher/ chair, because they could cause loss of balance of the operator and compromise the proper functioning of the device. If the path free cannot be made free of obstacles, choose an alternative path.
- Water, ice or debris on the stairs can cause loss of balance for the operator and jeopardize the proper functioning of the device. In order to avoid injury, clear the way of obstacles or choose an alternative route.
- For obstacles higher than 10mm the device must be raised, making sure to hold it by the main perimetral frame and not from the sidebars/ platforms or other parts not intended for such purpose.
- Once the wheels of the loading carriage have been positioned on the floor of the medical vehicle, the wheels of the front leg must have a distance from the ground of at least 5/6 cm, that allow safe opening and locking of the front legs. After every use check the height of the loading floor of the ambulance; if altered the stretcher settings must be verified by the manufacturer or by an authorized service centre. Otherwise we assume no responsibility on the improper functioning or any damage caused by the device itself.
- If the vehicle is equipped with pneumatic or hydraulic suspensions, the loading height adjustment must be carried out taking into account the worsening condition of use and/or the operation one provided by the vehicle builder.
- Problems during use and/or safety risks associated with such system, cannot be attributed to the manufacturer.
- **Improper installation of the loading platform can result in enervation and consequent damage of welding of the front legs.**
- **Improper installation of the loading platform may cause undesired operation of the device and cause harm to the patient and to the operator.**
- Do not alter or modify the stretcher arbitrarily to make it fit into the emergency vehicle: the modification may cause unforeseeable functioning and damage to the patient or operators. In any case the warranty will be void, relieving the manufacturer from any liability.
- CINCO MAS can be compliant with the standard EN 1789 if used with the Spencer fastener Sensor Lock or other Spencer fasteners intended for this stretcher. Therefore, is not permitted the use of fastening systems not approved by the manufacturer. Fastening systems that have not been approved may alter the structural and functional characteristics of the device.
- Do not operate on the adjustable heights system before the weight of the stretcher including patient and accessories has been carefully evaluated. Operators must be able to support the full load during the transition between the available heights. Wrong assessments could result in the sudden fall of the stretcher with consequent risks for the patient and for operators.
- Do not place magnetic objects between the stretcher and the fastener as they could interfere with its locking and unlocking system.
- Do not place limbs and / or objects between the legs and the frame, near the leg movement pistons and generally between moving parts as this could result in crushing injuries.
- Do not hit the piston of the front leg as a shock may cause unwanted release of the leg.
- If the twist system has been activated, before loading the stretcher on the ambulance, make sure that the wheels have been locked again as this procedure could be risky as well as difficult with the front wheels unlocked.
- Handling with four swivel wheels can be very difficult, especially in case of inclined or uneven ground. Carefully evaluate the conditions of use before unlocking the pivoting of the front wheels.

6.1 PHYSICAL REQUIREMENTS OF OPERATORS

CINCO MAS is a device intended for professional use only. Each operator must be trained to transport patients safely and efficiently. Do not allow untrained persons to help operators during use of the product, as this may cause injury to themselves or to other people.

The operators that use the device must have the physical ability to use it and good muscle coordination, as well as presenting strong back, arms and legs to raise and support and be able to grasp firmly the device with both hands.

Operators must be able to provide the necessary assistance to the patient.

Users should be able to lift and handle safely the weight of the stretcher, the patient and any other equipment used with the device

During patient loading procedures for extremely heavy patients, operations on rough terrains and in particular situations, more operators may be needed (not only 2 as expected in normal conditions)

The abilities of all operators must be considered before determining their role in the employment of the stretcher.

7. RESIDUAL RISK

The residual risks listed below have been identified exclusively in reference to the intended use of the device.

- The use by untrained personnel may result in injury to the patient, rescuer or others.
- Inadequate disinfection procedures may involve cross-infection risks.
- Partial opening of the legs, may result in a fall to the ground of the device. Make sure the legs are properly locked before making any movement of the stretcher and leg's pistons are fully extended and stable
- Improper locking of the stretcher on the fastener, or an improper positioning, can cause dangerous movement, especially if a violent deceleration of the vehicle happens. This may cause serious injury to patient and operators. Always verify the proper insertion of the locking system.
- Incorrect closure of the front telescopic handles of Cross Chair, may cause obstruction in the movement of the stretcher. Verify that the handles are fully closed during transport.
- Failure to comply with warnings for operators, may cause danger of crushing cause by handling mechanisms.
- Accidental activation of the adjustable heights unlocking command, could result in sudden fall of the stretcher with consequent injuries for patient and/or operators. Make sure the release handle is never accidentally activated
- Before activating the adjustable height operation lever, operators must prepare to support the full load of the stretcher with patient and accessories. Operation of this command without the application of sufficient force, would result in sudden fall of the stretcher with consequent injuries for patient and/or operators
- **Failure to read and understand the instruction manual of this device, can have consequences for the patient and operators.**

IT

EN

DE


EL

ES

HU

8. TECHNICAL DATA AND COMPONENTS

Note: Spencer Italia S.r.l. reserves the right to make changes to specifications without prior notice.

N°	DESCRIPTION	MATERIAL	N°	DESCRIPTION	MATERIAL
1	Sidebars	Al/Steel/Nylon	12	Leg support pad in closed position	PE/Nylon
2	Patient board (feet end)	PE	13	Trendelenburg selection lever	Steel
3	Main perimetral frame	Steel	14	Sidebar unlocking lever	Nylon/ Steel
4	Rear legs unlocking lever and lever to unlock the stretcher from the fastener	Al	15	Handle to operate the Twist System	Silicone/Al
5	Front legs unlocking lever	Al	16	Frame to unlock the adjustable height mechanism	Al/Rubber
6	Fifth wheel	Al/Nylon	17	Loading carriage wheels	PU
7	Rear piston	Steel /Al/Nylon	18	Patient backrest	PE
8	Rear pivoting wheels with brake	PU/Nylon/Steel	19	Trendelenburg/Fowler selection lever	Steel
9	Front wheels with twist system	PU/Nylon/Steel	20	Backrest adjustment piston	Steel /Nylon
10	Front leg	Steel	21	Front leg sliding wheels	Rubber/PP
11	Front piston	Steel/Al/Nylon			

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Specifications	Dimensions
Length (mm)	1970
Width (mm)	570
Loading carriage wheel height (cm)	65/75 ± 1
Loading carriage and leg's wheels (mm)	Ø 200
Loading capacity (kg)	250
Weight without standard equipment (kg)	45

1 Maximum load capacity means the total weight distributed according to the human anatomy. In determining the total load of the weight on the product, the operator must consider the weight of the patient, the equipment and the accessories.

9. INSTALLATION AND START-UP

Before the first use verify that:

- The packaging is intact and has protected the device during transport
- Check that are present all the components included in the accompanying list.
- Functionality of the device
- The sanitary vehicle is equipped with a Spencer fastening system dedicated to the stretcher.
- The patient board of the stretcher is flat.
- The loading floor of the stretcher wide and long enough to accommodate the stretcher and its accessories without becoming an obstacle.
- The wheels of the front leg, during loading and unloading from the vehicle, shall have a distance from the ground of at least 6 cm, allowing a safe opening and locking of the front legs (see the image at point 11.6).
- The fasteners must maintain the device fixed to the body of the vehicle.
- Do not modify for any reason the stretcher in the structural, leverage and traction parts. This type of intervention could cause damage to patient and/or rescuers.


⚠ If conditions mentioned above are not satisfied, the device cannot be considered secure, compromising its safe use and making it a possible cause of injury for patient, operators and damage to the device.

To facilitate the loading of the stretcher on the ambulance, it is recommended to eliminate sharp edges on the edge of the loading floor. The stretcher must be secured with Spencer fastening systems to prevent any kind of movement during transport on the sanitary vehicle. **Practice with a stretcher without patient before the regular start-up. For instruction about use after start-up, follow the operations described on paragraph 12.**

If the above conditions are met, the device may be considered ready for use; otherwise you must immediately remove the device from service and contact the manufacturer.

Do not alter or modify in any way the appliance; any such interference could cause malfunctions and injury to the patient and/or rescuer. Moreover the manufacturer will no longer recognize the product warranty and will accept no responsibility.

10. FUNCTIONAL CHARACTERISTICS

Element	Description	Function
A	Sidebar unlocking levers	By actuating the two levers at the same time, the sidebar is lowered.
B	Leg's release levers	This lever is identified by the symbol  . By actuating it while loading the stretcher, the rear legs are unlocked. When operated when the stretcher is coupled to the Sensor lock fastener, the stretcher is released from the fastener. It must also be used to bring the rear legs in intermediate height.
C	Front leg unlocking lever	This lever is identified by the symbol  . By actuating it, the front legs are unlocked in order to load the stretcher inside the vehicle.

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D	Rear wheels brake bracket	Is the steering wheel assembly. They are self-positioning, have a rotation of 360° and are equipped with stationing brake.
E	Front wheels bracket	They integrate the twist system which, if activated using the lever F, makes pivoting the front wheels
F	Twist activation lever	Operating it, the twist system is activated
G	Handle for intermediate height	By activating it, the leg at the head side is unlocked, thus allowing the stretcher to be moved to one of the intermediate heights
H	Backrest piston	Allows the adjustment of the backrest angle

11. INSTRUCTIONS FOR USE

Before transferring, lifting or transporting the patient, primary medical evaluations have to be performed. Once the diagnosis has been assured, it is preferable (if possible) to suggest to the patient that he actively collaborates during transfer onto the stretcher, and to make sure the patient is fully aware of all risks. Before loading the patient, place the device near the patient; always maintain the safety belts fastened during the manoeuvres.

11.1 MEDICAL VEHICLE REQUIREMENTS

Cinco Mas stretcher, is designed to enter and leave the medical vehicle compartment. The vehicle must have:

- A levelled loading floor
- Loading floor wide and long enough to permit a safe positioning of the stretcher. The stretcher must not be an obstacle for the operators.

The wheels of the front legs, when loading/unloading and the loading carriage is still on the loading floor of the vehicle, must be at a minimum distance of 6cm from the ground, to allow safe opening of the front legs.

If the above measure is not met, the device cannot be considered secure with consequent risks of injury for patient and operators and damage to the device.

11.2 STATIONING BRAKES

To engage the parking brakes, just press with a foot the strips on the brackets of the rear wheels. To disengage, simply press on the opposite side of the brake, they will return automatically in original position.


11.3 SIDEBARS

The stretcher is equipped with sidebars, necessary to contain the patient on the stretcher.

NEVER MOVE THE PATIENT AND NEVER LEAVE HIM ON THE STRETCHER WITHOUT FIRST LIFTING THE SIDES. Failure to observe this warning could result in serious damage to the patient.

To lower the sidebars, simultaneously pull towards yourself the levers located at the bottom and identified by the word PULL. The sidebar will open automatically. To close the sidebar, lift it back to its original position and check that it is correctly locked by repeatedly pulling and pushing against them. When closing, make sure that nothing interferes with the locking systems. For example, the sheets may prevent them from closing properly.

11.4 ACTIVATION OF TWIST SYSTEM

This system allows to unlock the pivoting of the front wheels. This rotation is released by lifting one of the yellow levers placed at the sides of the stretcher.


When the rotation is no longer needed, lower the lever and push the stretcher firmly forward, the wheels will be locked again.

To use the twist system, at least one additional operator must be present to control the head side of the stretcher as this is much more difficult to maneuver with all four wheels unlocked.

Make sure you have carefully assessed the ground conditions before operating the twist system, the stretcher could move unexpectedly and uncontrollably **lowering of the stretcher.**


11.5 LOWERING THE STRETCHER

To facilitate the transfer of patient, it's advisable to bring the stretcher in intermediate height or lowered position (only if the patient is not on the device).

To load an immobilized or hospitalized patient on the stretcher at intermediate height, proceed as follows:

- Operate the rear leg release lever, hold it and lift the stretcher slightly to unlock the rear leg movement. Lower the stretcher about 10cm, release the lever and accompany down the stretcher until it reaches the first intermediate position. Ensure that the stretcher has reached a stable position.

IT

EN


DE

EL

ES

HU


- To lower the front part of the stretcher, it is necessary to operate the shaped handle located just below the front part of the frame (command G par. 10). Lift the shaped handle holding firmly the stretcher. **It is necessary to support the weight of the stretcher, the patient and any equipment applied to the stretcher itself.** After starting the descent movement, release the shaped handle while holding the frame of the stretcher, until it reaches the first available intermediate position. Make sure that the stretcher is stable.
- To position the front part at another intermediate height, operate the handle as previously done. This time, however, it is necessary to lift the stretcher slightly at the front to unlock the mechanism. Always holding the grip firmly on the stretcher, supporting its weight, accompany it to the next intermediate position.

Procedure to load/unload not hospitalized patient

If the patient to be loaded has motion difficulties but has the ability to position himself on the stretcher without the need to be lifted, proceed as follows:

- 1 - Lower the back of the stretcher (handle side) by 2 or a maximum of 3 positions.
- 2 - Activate the brakes of the rear wheels to avoid unwanted movements of the stretcher
- 3 - Lower the front part of the stretcher (patient head side) by one position only and adjust the backrest as required.
- 4 - Assist the patient during the ascent on the stretcher, making sure that the first leg to be lifted, rests on the footrest of the patient surface. When correctly positioned, help him to lift the other leg.
- 5 - Immobilize the patient with the belts and lift the sidebars
- 6 - Once the patient is lying on the stretcher and properly immobilized, lift the stretcher carefully and return it back to the horizontal position.

Lowering the stretcher to the ground

This procedure allows the stretcher to reach the minimum distance from the ground.

-  **It is forbidden to make this manoeuvre when patient is on the device.**
- Raise the back of the stretcher until the loading carriage wheels have reached the ground.
- Keeping the loading carriage wheels on the ground, and ready to support the weight of the device, operate both levers and push up stretcher to unlock front legs. Once the legs are unlocked, accompany it to the ground.

Note: When the stretcher is fully lowered, the brakes doesn't perform their function. Make sure that the stretcher is held in place by at least one operator.


IT

11.6 STRETCHER LIFTING

To bring back the stretcher to the standard height, starting from any of the previous configurations, the operators must be coordinated while lifting the front and rear of the assembly, in order to ensure the proper alignment of the patient's board. Carry out these operations only if the immobilization of the patient has been carried out and if the sidebars are lifted.

- In the rear part, grab the part of the frame near the footrest of the patient board.
- In the front part, similarly grasp the frame above the loading trolley.
- Lift the assembly until the proper locking of the legs has been obtained.

-  **For the lifting procedures, use only the main frame of the stretcher.**
- **For the lifting, don't grab the area where the patient's lie or other parts as they are not designed for this purpose.**

EN

DE

EL


ES

11.7 LOADING AND UNLOADING THE STRETCHER FROM THE MEDICAL VEHICLE

 The loading and unloading phases of the stretcher are among the most delicate during the use of the stretchers. It is necessary to follow with extreme attention all the warnings given in this manual, practising in simulated conditions before placing the product into service. The patient shall always be properly immobilized. **Always make sure that the conditions of use always allow the correct opening of the front legs before carrying out any loading and/or unloading operation.**

To load the stretcher on the vehicle proceed as follows

- The operator which is loading the stretcher, must be aligned with the stretcher and its fastener so that both devices are in front of him and aligned as required by their intended use.
- Advance the wheels of the loading carriage inside the vehicle until the front legs of the stretcher have reached the bodywork of the vehicle.
- If a slide/ramp is used, the wheel of the loading trolley must rest on the flat part of the ambulance floor before the front legs are in contact with the vehicle bodywork.
- Ensure that the wheels of the front legs are at least 5cm from the ground. Failure to comply with this specification could cause serious damage and/or injuries when unloading the stretcher from the vehicle.

HU

- Operate the right lever (identified by the symbol ▲ and the red colour), to unlock the movement of the front legs, continuing to push the stretcher inside the vehicle until the rear legs have also reached the bumper of the vehicle.
- Only after making sure that the rear legs are in contact with the bumper of the vehicle, operate the release lever of the rear legs (identified by the symbol ▼ and the green colour also present on the legs plastics), continuing to push towards the inside of the vehicle compartment.

⚠ In this procedure, a considerable part of the weight of the assembly, will weight on the operator, which must be strong enough to support and accompany the device in all its movements.

- Lock the stretcher on the medical vehicle using the Spencer fixing system installed on the vehicle.



⚠ **Warning: DURING THE NORMAL MOVEMENTS OF THE STRETCHER, NEVER OPERATE THE HANDLES WHICH ARE DEDICATED TO THE UNLOCKING OF THE LEGS BECAUSE THE STRETCHER COULD FALL TO THE GROUND. THESE LEVERS MUST ONLY BE OPERATED TO LOAD THE STRETCHER ON THE AMBULANCE OR TO LOWER THE STRETCHER AS DESCRIBED ABOVE.**

To unload the stretcher from the medical vehicle, proceed as follow:

- If the stretcher is used with the Sensor Lock fastener, release it by operating the rear leg release lever located on the stretcher and identified by the symbol ▼. If the stretcher is used with the S-MAX fastener, the release must be done manually by operating the dedicated lever placed on the fastener.
- ▲ Pull the stretcher outside of the vehicle, using the handle point in the rear of the stretcher near the footrest. Support the weight of the device until the locking mechanism of the rear legs has been inserted. **As for the loading procedures, the operator must be able to support the weight of the device.**
- Pull the stretcher out from the vehicle until the front legs are completely opened. **The loading carriage can't leave the loading floor until the locking of the front legs hasn't been verified.**
- Terminate by extracting the stretcher from the vehicle.

11.8 TRENDELENBURG PLATFORM

If your stretcher has the possibility to place the platform in the Trendelenburg position, grab the patient platform at the footrest, then pull the selection lever shown in the figure outwards. There are 3 different possible inclinations in addition to the horizontal one.

11.9 TRENDELENBURG/FOWLER PLATFORM


The Trendelenburg/Fowler movement, if present, allows the lifting of the patient platform in different positions with consequent raising of the patient's lower limbs.

Change the platform inclination from horizontal to Fowler position

If the stretcher allows the platform to be positioned in the Fowler position, lift the intermediate part in correspondence with the patient's knees, so that the pins on the sides of the patient's platform are positioned in the slot of the metal rod (no. 19 par. 8), closest to the sidebars. Make sure that the pins are well inserted and solid, thus ensuring the stability of the platform in this position.

Changing the platform inclination from a horizontal to Trendelenburg position

With one hand support and raise the platform by grasping the footrest, select the Trendelenburg position by placing the screw of the platform, in the slot of the metal rod closest to the control handle. Make sure that the screw fits into the correct seat and is attached to it.

11.10 BACKREST ADJUSTMENT

Always inform the patient when an adjustment is to be made.


Changing the backrest inclination from a horizontal to a vertical position

Lift the backrest to the first position, where the backrest will automatically lock. Follow the same procedure to reach the next position, always checking the correct insertion of the locking system.

Changing the backrest inclination from a vertical to a horizontal position

Support the backrest structure with one hand (to avoid sudden movement) and release weight by moving the backrest upwards.

At the same time, pull the pull of the backrest piston upwards until the safety mechanism is released, then accompany it with the other hand, lower the backrest towards the horizontal position to the desired height, then return the knob to the rest position. To reach the next positions, lift and accompany the movement of the backrest. Always unload the weight of the backrest before operating the knob. An incorrect procedure can irreversibly damage the backrest lifting piston

12. CLEANING AND MAINTENANCE

12.1 CLEANING

Failure to carry out the correct cleaning routine could increase the risk of cross infection, due to presence of body fluids and/or residuals.
The operator must always wear adequate personal protection such as gloves and mask, etc. during all checking and cleaning procedures.

The exposed metal parts are usually treated and/or painted in order to increase their resistance to external agents. Clean the exposed parts with water and delicate soap; **never use solvents or stain removers.**

Rinse thoroughly with warm water making sure that you have removed all traces of detergent, which could degrade or compromise the integrity and durability of the device. **The use of high pressure water should be avoided** because water penetrates the joints and removes the lubricant, creating the risk of corrosion of components. Allow the device to dry thoroughly before storing. Drying after washing or after use in wet environments must be natural and not forced; do not use flames or other sources of direct heat.

If you need to carry out **disinfection**, use products that do not have corrosive or solvent action on the materials constituting the device. Be sure to take every precaution to ensure that there is no risk of cross-infection or contamination for patients and / or operators.

12.2 PRECAUTIONARY MAINTENANCE

Establish a maintenance program and periodic testing routine, identifying an employee responsible for this. The person to whom the ordinary maintenance of the device is entrusted must ensure the basic requirements foreseen by the manufacturer in following paragraphs are inspected.

All maintenance and periodic servicing activities must be registered and kept together with the servicing reports. These documents have to be kept for a period of 10 years after the disposal of the device itself. This register will be made available to the competent authorities and/or manufacturer if requested.

In order to manage the maintenance and service procedures of your equipment, Spencer provides you with the SPENCER SERVICE portal <https://www.spencer.it/en/services/after-sales> through which you can view the data of your products or their placing on the market, monitor and update the schedules of periodic revisions and read and manage the extraordinary maintenance operations.

Routine maintenance of the device must be carried out by operators in possession of specific qualifications, trained and experienced in the use and maintenance of the device.

The operator must always wear adequate personal protection such as gloves and mask, etc. during all checking and cleaning procedures.

Checks to be carried out before and after each use and at deadline indicated above, are as follows:

- General functionality of the device
- Cleanliness of the device (remember that the failure of cleaning may cause the risk of cross infections)
- Correct fixation of all nuts, bolts and screws
- Absence of cuts, holes, tears on the structure, including the straps
- None of the tubes or metal sheets present bends or cracks
- Welded areas are intact, without any cracks or breaks
- Moving parts, wheels, levers, handles are intact and work properly
- Lubrication of moving parts
- State of use of wheels and breaking system
- The wheels are correctly fixed, they are stable and turn properly
- The wheels are free from debris.
- The device opens and locks properly
- The device opens and closes properly
- Functioning of springs
- The stretcher easily enters the ambulance
- The emergency vehicle is equipped with a Spencer fastening system intended for the stretcher
- The coupling between the fastener and the stretcher is suitable to guarantee the safety of the fixing.

The inspection frequency is determined by factors such as legal requirements, the type of use, frequency of use, environmental conditions during use and storage.

Please note that you must do the cleaning as described in this manual and verify functionality before and after each use. Spencer Italia S.r.l. declines any responsibility for the improper functioning or damages caused to the patient or user by the use of devices not subject to routine maintenance and will void the warranty and the compliance to the Regulation 2017/745/UE.

Use only accessories/original spare parts approved by Spencer Italia S.r.l., otherwise we will accept no responsibility for the incorrect functioning and/or damage caused by the use of any device which has not been repaired, or certified on expiry date by the manufacturer or by one of the manufacturer's authorised service centres. Warranty will be considered void in compliance with the Regulation 2017/745/UE.

12.3 PERIODIC MAINTENANCE

The device must be serviced by the manufacturer or by an authorised centre, every year.

If the correct revision is not carried out, the device **MUST BE PUT OUT OF SERVICE**, because the CE branding will no longer be considered valid and consequently it is no longer compliant with the safety standards declared by the manufacturer at time of purchase.

Spencer Italia S.r.l. will take no responsibility for the incorrect functioning or any damage caused by a device that has not undergone regular revision.

Mattress and belts shall be replaced every two years.

For any operations that are not carried out directly by the manufacturer but by an authorised centre, we have to underline that a report regarding all operations carried out must be requested. This will permit both Spencer Italia S.r.l. and the end user to keep a log book regarding the operations carried out on the device.

12.4 SPECIAL SERVICING

Only the manufacturer or centres with written authorisation are authorised to complete any special servicing operations.

For any operations that are not carried out directly by the manufacturer but by an authorised centre, we have to underline that a report regarding all operations carried out must be requested. This will permit both Spencer Italia S.r.l. and the end user to keep a log book regarding the operations carried out on the device.

The end user is authorised to replace only the spare parts indicated in the paragraph 15.

12.5 LIFE SPAN

The device, if used as indicated in the following instruction manual, has a life span of 5 years starting from the purchase date. The life span can be extended following the annual revision.

General revisions must be carried out by the manufacturer or by a centre authorised by the manufacture. **If such annual revisions are not carried out, the device MUST BE DISPOSED ACCORDING TO THE PROCEDURES SPECIFIED IN PARAGRAPH 16 AND THIS EVENT MUST BE NOTIFIED TO THE MANUFACTURER.**

The life span can be extended, by unquestionable judgment of Spencer or authorized service center, if the safety requirements of the device are still guaranteed.

Spencer Italia S.r.l. will accept no responsibility for the incorrect functioning and/or damage caused by the use of any device which has not been serviced by the manufacturer or authorized center, or of any device for which the life span is expired.

13. TROUBLESHOOTING

PROBLEM	CAUSE	SOLUTION
The leg's release levers don't work properly or are difficult to operate	The moving mechanisms are damaged	Put the stretcher immediately out of service and contact the service centre
	Means of connection between components have been lost.	Put the stretcher immediately out of service and contact the service centre
The fixation of the stretcher on the fastener is not adequate	Wear or damage of the components that constitute the locking mechanisms.	Put the stretcher immediately out of service and contact the service centre
	The fifth wheel has not been properly inserted in the fastener	Position the stretcher properly and check that the fifth wheel enters the dedicated part on the fastener.
Damage to the structure	Improper use	Put the stretcher immediately out of service and contact the service centre
It's not possible to place the stretcher in intermediate position.	The moving mechanisms are damaged	Put the stretcher immediately out of service and contact the service centre
	Something interferes with the moving systems	Look for objects that could interfere with the moving systems
	Levers have not been properly operated	Carefully follow the instructions to place the stretcher in intermediate position
During unloading procedures, front legs don't reach the locking position.	The moving mechanisms are damaged	Put the stretcher immediately out of service and contact the service centre
	The height of the loading floor is not suited to the device, the security height described in the instruction manual is not respected.	Adjust the loading floor height to meet the requirements described in this manual. If the loading floor does not allow adjustments, put the device immediately out of service and contact the service centre.
When the stretcher release lever is operated, the stretcher does not move from the fastener and remains attached.	The magnet controlling the release, located inside the fifth wheel holder, has moved or been lost	Unlock manually the stretcher and finish the unloading operation. At the end of service, check the positioning of the magnet, then restore its original position. If the problem persists, immediately put the device out of service and contact the service centre.
	The fastener Sensor Lock is not powered, is locked or is faulty.	Check the power supply of the fastener. If the problem persists, immediately put the device out of service and contact the service centre.
	The stretcher is used with a fastener that does not provide automatic release.	If necessary, request a different fastener

14. ACCESSORIES

14.1 STANDARD EQUIPMENT

Depending on the model, Cinco Mas stretcher can be supplied with belts with automatic reel and mattress.

14.2 ACCESSORIES

ST42706C	10G FASTENER SENSOR LOCK	CB09025C	TANKER YELLOW OXYGEN TANK HOLDER
IF01047C	TRACK 4-30 – TELESCOPIC IV HOLDER	ST00491A	STX 90 TELESCOPIC HEADREST FOR STRETCHER
IF01049B	TRACK 5 TELESCOPIC IV HOLDER	CR90013C	LOADING CARRIAGE EXTENDER FOR CINCO MAS
EN90003C	END-T 10G INSTRUMENTS TABLE	CR90014B	LOADING CARRIAGE HEIGHT REDUCER 4cm FOR CINCO MAS
CB09028C	TANKER SUPPORT FOR SAFEBAR		

15. SPARE PARTS

ST00497B	DNA STRAP belt with rewinding system	ST42021A	Black Wheel ø200mm w/bearing
ST00498B	DNA STRAP thoracic belt w/re-winding system	ST70606A	Black Wheel ø200 mm without bearing

16. DEMOLITION

When the device is no longer suitable for use, if they haven't been contaminated by any particular agents, they can be disposed of as normal solid waste, otherwise follow the current regulations for demolition.

Warning

The information contained in this document could be modified without any warning and is not to be intended as a commitment on behalf of Spencer Italia S.r.l. Spencer products are exported to many countries and the same identical regulations are not always valid. For this reason there could be differences between the description here described and the product actually delivered. Spencer continually strives to reach the perfection of all items sold. We therefore hope you will understand if we reserve the right, at any time, to modify the shape, equipment, lay-out or technical aspects that are herein described.

© Copyright Spencer Italia S.r.l.

All rights reserved. No part of this document can be photocopied, reproduced or translated into another language without the written approval of Spencer Italia S.r.l.

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If the correct revision is not carried out, the device **MUST BE PUT OUT OF SERVICE**, because the CE branding will no longer be considered valid and consequently it is no longer compliant with the safety standards declared by the manufacturer at time of purchase. Spencer Italia S.r.l. will take no responsibility for the incorrect functioning or any damage caused by a device that has not undergone regular revision. Mattress and belts shall be replaced every two years. For any operations that are not carried out directly by the manufacturer but by an authorised centre, we have to underline that a report regarding all operations carried out must be requested. This will permit both Spencer Italia S.r.l. and the end user to keep a log book regarding the operations carried out on the device.

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For any operations that are not carried out directly by the manufacturer but by an authorised centre, we have to underline that a report regarding all operations carried out must be requested. This will permit both Spencer Italia S.r.l. and the end user to keep a log book regarding the operations carried out on the device. The end user is authorised to replace only the spare parts indicated in the paragraph 15.

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