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
EN

1.	MODELS	16
2.	INTENDED USE	16
3.	REFERENCE STANDARDS	16
4.	INTRODUCTION	16
4.1	Using the manual	16
4.2	Device labelling and traceability	16
4.3	Symbols	17
4.4	Warranty and service	17
5.	WARNINGS	17
6.	SPECIFIC WARNINGS	18
6.1	Physical requirements for operators	19
7.	RESIDUAL RISK	19
8.	TECHNICAL DATA AND COMPONENTS	19
9.	COMMISSIONING	20
10.	OPERATING CHARACTERISTICS	20
11.	PROPER USE	21
11.1	Emergency vehicle requirements	21
11.2	Parking brakes	21
11.3	Sidebars	21
11.4	Twist system activation	21
11.5	Lowering the stretcher	21
11.6	Lifting the stretcher	22
11.7	Loading and unloading the stretcher from the medical vehicle	23
11.8	Trendelenburg position	23
11.9	Fowler position	24
11.10	Backrest adjustment	24
11.11	Lowering the loading trolley	24
11.12	S-Wings application	24
11.13	Restraint straps	24
12.	CLEANING AND MAINTENANCE	25
12.1	Cleaning	25
12.2	Routine maintenance	25
12.3	Periodic overhaul	26
12.4	Special maintenance	26
12.5	Life span	26
13.	TROUBLESHOOTING TABLE	26
14.	ACCESSORIES	27
15.	SPARE PARTS	27
16.	DISPOSAL	27

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PT

1. MODELS

The following basic models may be subject to implementation or change without notice.

- AGILE – ROLL-IN STRETCHER

The above-listed models may be marketed in different variants and, depending on the model, may include accessories. The models are available with a height of 60/68/75.

2. INTENDED USE

Ambulance stretchers are the primary means of transporting sick and/or injured people in transit safely and comfortably in a lying position. It is not foreseen that the patient be able to intervene on the device, nor is it intended for prolonged storage or use as a hospital bed. The stretcher must be used with the dedicated fixing systems. It is not foreseen that the patient be able to intervene on the device.

TARGET PATIENTS

There are no particular indications related to the patient group.

The product configuration is able to accommodate any subject as long as he/she is within the maximum capacity of the device. If paediatric subjects must be transported, it will be in the role of the rescuer to determine whether the strap systems are suitable for immobilization or if it will be necessary to use another device.

PATIENT SELECTION CRITERIA

The intended patients are those for whom ambulance transport is necessary.

CONTRAINDICATIONS AND UNWANTED SIDE EFFECTS

No particular contraindications or side effects are known with relation to use of the device, as long as it is used in accordance with the user manual.

USERS AND INSTALLERS

The intended users are persons trained in first aid procedures and in use of medical equipment in an emergency medical service (EMS) environment.

Possible users also include emergency vehicle outfitters who may use the product before commissioning or during maintenance of the vehicle on which the stretcher is used.

User training

Note: Despite all efforts, laboratory tests, trials, and instructions for use, standards do not always reproduce practice, so the results obtained under actual conditions of product use in the natural environment may sometimes differ significantly.

The best instructions are the continuous practice of use under the supervision of competent and trained personnel.

- Regardless of your level of experience with similar devices in the past, you should carefully read and understand the contents of this manual before installing, operating, or servicing this product. In case of any questions, please contact Spencer Italia S.r.l. for the necessary clarifications.
- The product must be used only by personnel trained in the use of this product and not on other similar products.
- The suitability of the users for use of this product can be attested by the training registration, in which trained persons, trainers, date and place are specified. **This documentation must be kept for at least 10 years after the end of the product's life and must be made available to the competent authorities and/or the Manufacturer when requested. In the absence of such documentation, the relevant bodies will apply any foreseen sanctions.**
- Do not allow untrained persons to help while using the product, as they may cause injury to themselves or others.
- The product must be put into use only by personnel trained in the use of this product and not on other similar products.

Note: Spencer Italia S.r.l. is always available for training courses. Document user training using the form in Annex A of this manual or other appropriate forms.

Installer training

The device must be installed by qualified personnel who are trained and authorised to use and install the device.

The installer must strictly adhere to these instructions as well as to the state of the art for vehicle installations.

3. REFERENCE STANDARDS

As Distributor or End-User of the products manufactured and/or marketed by Spencer Italia S.r.l., users are strictly required to be familiar with the legal provisions in force in the country of destination of the goods, applicable to the devices to be supplied (including regulations relating to technical specifications and/or safety requirements) and, therefore, to understand the requirements necessary to ensure compliance of the products themselves with all legal requirements of the territory.

REFERENCE	DOCUMENT TITLE
REGULATION (EU) 2017/745	REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 5 April 2017 on medical devices
UNI EN ISO 1865-1	Patient handling equipment used in road ambulances - Part 1: General stretcher systems and patient handling equipment
UNI EN ISO 1865-3	Patient handling equipment used in road ambulances - Part 3: Heavy duty stretchers
UNI EN 1789	Medical vehicles and their equipment - Road ambulances

4. INTRODUCTION

4.1 USING THE MANUAL

The purpose of this manual is to provide healthcare professionals with the information necessary for safe and appropriate use and maintenance of the device. **For correct reading of this manual, please refer to the pictures at the beginning of the manual.**

Note: The Manual is an integral part of the device and therefore it must be kept for the entire life of the device and must accompany it in any changes of use or ownership. If any instructions for use for products other than the one received are present, please contact the Manufacturer immediately before use.


Spencer product User Manuals can be downloaded from the site www.spencer.it/en/emergency-medical-solutions/user-manuals or by contacting the Manufacturer. Exceptions are those items whose essentiality and reasonable and predictable use are such that it is not necessary to draw up instructions, in addition to the following warnings and indications on the label.

Regardless of your level of experience with similar devices in the past, it is advisable to carefully read and understand the contents of this manual before installing, operating, or servicing this product.

4.2 DEVICE LABELLING AND TRACEABILITY

Each device is provided with a label, placed on the device itself and/or on the packaging, which contains the Manufacturer's identification data, product, CE marking or serial number (SN). **This must never be removed or covered.**

The following image shows the positioning and main information contained on the label.

In the event of damage or removal, request a duplicate from the Manufacturer, or else the warranty will be void as the device can no longer be traced.

Regulation 2017/745/EU requires manufacturers and distributors of medical devices to keep track of their location. If the device is in a location other than the address to which it was shipped or sold, or if it was donated, lost, stolen, exported or destroyed, permanently removed from use, or if the device was not delivered directly by Spencer Italia S.r.l., inform Customer Service (see § 4.4).

4.3 SYMBOLS

Symbol	Meaning	Symbol	Meaning
CE	Device in compliance with EU Regulation 2017/745	⚠	Danger – Indicates a hazardous situation that may result in a situation directly related to serious injury or death.
MD	Medical device	📖	See the user manual.
🏭	Manufacturer	REF	Product code
📅	Date of manufacture	SN	Serial Number
🚫	Hazard of crushing fingers	👤	Required reading of Instructions for use

Unique Device Identifier

Production identification

Alphanumeric code that identifies the production units of the device, composed of:

(01)0805771123XXXX	Company prefix (XXXX=GTIN)
(11) YYMMDD	Date of manufacture
(21) 2100443030	(21) 2100443030 ; 21=>SN

4.4 WARRANTY AND SERVICE

Spencer Italia S.r.l. guarantees that products are free from defects for a period of **one year from the date of purchase**.

For information regarding correct interpretation of the instructions for use, maintenance, installation or return, please contact Spencer Customer Service tel. +39 0521 541111, fax +39 0521 541222, e-mail service@spencer.it.

To facilitate service, always indicate the lot number (LOT) or serial number (SN) on the label attached to the package or device itself.

Warranty and service conditions are available at www.spencer.it/en/support/terms-and-conditions

Note: Record and keep with these instructions: lot (LOT) or serial number (SN), if present, place and date of purchase, date of first use, date of checks, user name and comments.

5. WARNINGS

Warnings, notes, and other important safety information are provided in this section and are clearly visible throughout the manual.

At least every 6 months, it is important to check for updated instructions and any changes involving your product. This information is freely available on the website www.spencer.it/en on the specific product page.

Product features

Use of the product for any purpose other than that described in the User Manual is prohibited.

- Before each use, always check the conditions of the product, as specified in the User Manual. In the event of faults/damage that could compromise its functionality/safety, immediately remove it from service and contact the Manufacturer.
- If the product is found to be malfunctioning, immediately use a similar device to ensure continuity of ongoing operations.
- The product must not be tampered with (modification, tweaking, additions, repair). Should these operations be performed, we decline any responsibility for incorrect operation or any damage caused by the product itself; moreover, the CE certification (when required by law) and the product warranty shall be null and void.
- Anyone who modifies or causes products manufactured by Spencer Italia S.r.l. to be modified or re-processed in such a way that they no longer serve the intended purpose or no longer provide the intended performance shall comply with the conditions valid for the first placing on the market.
- When using the devices, position and adjust them in such a way that they do not hinder operator works or the use of any other equipment.
- Be sure to take every precaution to avoid hazards from contact with blood or body secretions, if applicable.
- Always comply with the maximum capacity indicated in the User Manual. Maximum load capacity means the total weight distributed according to human anatomy. When determining the total weight load on the product, the operator should consider the weight of the patient, equipment and accessories. Moreover, the operator should assess whether the overall size of the patient reduces the functionality of the product.
- Make sure that operators are in good physical condition before lifting, as listed in the User Manual.
- **The maximum weight, which weighs on each operator, must comply with local health and safety requirements.**
- Avoid contact with sharp objects.
- Device installation must be carried out by qualified personnel trained and enabled by Spencer Italia S.r.l. The times and methods of such shall be agreed upon between the customer and our Sales Offices.
- Operating temperature: from -10°C to + 50°C.

Storage

- The product must not be exposed or come into contact with thermal sources of combustion or flammable agents, but must instead be stored in a dry, cool place, away from light and sun.

IT
EN
DE
FR
ES
PT

- Do not store the product under other more or less heavy materials that may damage the structure of the product.
- Store and transport the product with its original packaging, otherwise the warranty shall be invalidated.
- Storage temperature: -20°C to +60°C.

Regulatory requirements:

- As Distributor or End-User of the products manufactured and/or marketed by Spencer Italia S.r.l., users are strictly required to be familiar with the legal provisions in force in the country of destination of the goods, applicable to the devices to be supplied (including regulations relating to technical specifications and/or safety requirements) and, therefore, to understand the requirements necessary to ensure compliance of the products themselves with all legal requirements of the territory.
- Promptly and in detail notify Spencer Italia S.r.l. (already in the quotation request phase) about possible fulfillments by the Manufacturer necessary for the compliance of products with specific legal requirements of the territory (including those deriving from regulations and/or regulatory provisions of another nature).
 - Act with due care and diligence to help ensure compliance with the general safety requirements of the devices placed on the market, providing end-users with all the information necessary to carry out periodic revisions on the supplied devices, exactly as indicated in the User Manual.
 - **Participate in safety checks on products** placed on the market, transmitting information regarding product risks to the Manufacturer as well as to the Competent Authorities for their respective actions.
 - Without prejudice to the above, the Distributor or End-User shall assume wider liability related to non-compliance with non-fulfilment of the above-mentioned obligations, with consequent obligation to indemnify and/or hold Spencer Italia S.r.l. harmless from any possible injurious effect.
 - With reference to EU Regulation 2017/745, please note that public or private operators who, when exercising their activity, detect an incident involving a medical product are required to notify the Ministry of Health, within the terms and in the manner established by one or more ministerial decrees, and notify the Manufacturer. Public or private health care professionals are required to notify the Manufacturer of any other incident that may allow the adoption of measures to ensure the protection and health of patients and users.

General warnings for medical devices

The user must carefully read the following in addition to the general warnings.

- The application of the device should not last longer than the time required for first aid operations and subsequent transport to the nearest rescue point.
- Qualified personnel and at least two operators must be present during use of the device.
- Follow the internal procedures and protocols approved by your organisation.
- Disinfection operations must be carried out in accordance with the validated cycle parameters, as stated in the specific technical standards.

6. SPECIFIC WARNINGS

To use the product, you must also have read, understood and carefully follow all the instructions in the user manual.

- Follow approved Emergency Medical Service procedures for patient immobilization and transportation.
- Follow approved Emergency Medical Service procedures for patient positioning and transportation.
- Do not use if the device or parts of it are punctured, torn, frayed, or excessively worn.
- Before any type of handling, make sure that operators have a secure grip on the device.
- Do not drag the device on uneven surfaces.
- Do not lift with a crane or other mechanical lift.
- Do not use drying machines.
- The device is a transport device and cannot be used as a stationing device.
- Do not use with devices other than those expressly approved by the manufacturer.
- Practice with a device without a patient to make sure you are familiar with the manoeuvres.
- For patient loading techniques for particularly heavy patients, for operations on steep terrain or in special and unusual circumstances, the presence of multiple operators is recommended (not just 2 as expected in standard conditions).
- Make sure that the patient is adequately immobilized before loading him/her on the stretcher. A failure to immobilize can cause serious damage to the patient.
- Make sure that the sheet does not interfere with any stretcher handling or control mechanisms.
- Do not move the device if the weight is not properly distributed.
- Always use straps that are anchored to the stretcher frame to ensure patient safety.
- Use only the perimeter frame or handles to move the stretcher and not the sidebars, beds or other points not intended for this purpose.
- Avoid excessive force when loading the stretcher into the ambulance: unnecessary force can cause damage and can adversely affect the fastening system.
- Avoid excessive force when loading the stretcher into the medical vehicle: unnecessary force can cause damage and can adversely affect the functionality of the stretcher itself.
- **Hold the device firmly if the patient is lying on it.**
- **Parking brakes are aids for the operator but they do not in any way replace his or her supervision.**
- Pay the utmost attention to any obstacles (water, ice, debris, etc.) present on the route, as they may cause the operator to lose balance and compromise proper functioning of the device. If you cannot clear the way, choose an alternative route.
- Condensation, water, ice and accumulations of dust can affect correct operation of the device, making it unpredictable and causing a sudden change in the weight the operators have to support.
- For height differences greater than 10mm, the device must be lifted, taking care to grasp it from the structure and not from the sidebars/beds or other points not intended for this purpose.
- Once the wheels of the loading trolley have been placed on the ambulance floor, the wheels of the front leg must have a ground clearance of between at least 5 and 10 cm, which allows the front leg to be safely opened and locked. Check the height of the ambulance bed after each use; if it is altered, the trolley must be adjusted immediately by the manufacturer or by a specialised technician authorised by the manufacturer. If not, the manufacturer cannot be held responsible for correct functioning of the device or any damage caused by it.
- If the vehicle is equipped with pneumatic or hydraulic suspensions, the loading height must be adjusted taking into account the worst-case conditions and/or the working conditions foreseen by the outfitter.
- Any problems regarding the use and/or safety risks associated with this system are not the responsibility of the manufacturer.
- **Inadequate installation of the loading floor can lead to buckling and subsequent damage to the welds of the front legs.**
- **Improper installation of the loading floor can lead to abnormal operation of the device and cause harm to the patient and user.**
- Do not alter or modify the stretcher to adapt it to the emergency vehicle: doing so could result in unpredictable operation and damage to the patient or rescuers and shall void the manufacturer's warranty and release the Manufacturer from all liability.
- The product can only comply with regulation EN 1789 if used with the Spencer 10g fixing system compatible with this stretcher. The use of fasteners not approved by the Manufacturer is therefore prohibited. Non-approved fixing systems may alter the structural and functional characteristics of the device.
- Do not place limbs and/or objects between the legs and the frame, in the vicinity of the leg lifting pistons or in general between moving parts as this may result in crushing injuries.
- Direct contact with the patient's skin is not foreseen. In the event of direct contact of the patient's skin with the surfaces of the stretcher, use a protective sheet that complies with biocompatibility requirements in order to protect both the patient and the surfaces of the device.
- When loading the stretcher onto a vehicle or lowering it to the ground, hands must return to a position on the tubular frame, holding it firmly, after operating the controls. Incorrect positioning of hands could result in them being crushed between the frame and the wheel bracket.
- If the device is used in conjunction with restraint systems such as spine boards and/or vacuum mattresses, make sure that the patient is secured to the stretcher and the restraint system so that the patient is safe while the vehicle is moving. In case of doubt about the procedures to be used, refer to the operating protocols of your relative 118 service.

6.1 PHYSICAL REQUIREMENTS FOR OPERATORS

The device is intended for professional use only. All operators must be trained to transport patients safely and efficiently. Do not allow untrained persons to help while using the product, as they may cause injury to themselves or others.

Operators using the device should be physically able to use the device and have good muscle coordination, as well as strong back, arms, and legs for lifting and holding and be able to grasp the device firmly with both hands.

Operators must be able to provide the necessary patient care.

Users must be able to safely lift and handle the weight of the stretcher and patient together and any other equipment used with the device.

For patient loading techniques for particularly heavy patients, for operations on steep terrain or in special and unusual circumstances, the presence of multiple operators is recommended (not just 2 as expected in standard conditions).

The ability of each operator must be assessed before the definition of rescuer roles in use of the stretcher.


7. RESIDUAL RISK

The residual risks listed below have been identified solely with reference to the intended use of the device.

- Use by untrained personnel may result in injury to the rescuer or third parties.
- Inadequate disinfection procedures may lead to risk of cross-infection.
- Partial opening of the legs may result in the device falling to the ground. Make sure that the legs are properly locked before any handling and that the pistons are fully extended and stable.
- If the stretcher is not locked on the fixing system or if it is not correctly positioned, this can result in dangerous movements, especially in the event of strong medical vehicle deceleration, which can cause damage to the patient and the operators. Always make sure that the locking system is properly engaged.
- Failure to observe the warnings for operators can lead to risks of crushing caused by the moving mechanisms.
- Accidental operation of the variable height system may result in the stretcher falling causing injury to the patient and/or operators. Make sure that the unlocking handle is not accidentally activated.
- The operators should prepare to fully load the stretcher with the patient and accessories before operating the variable height unlocking control. Operating this control without applying sufficient support force will result in the stretcher falling unexpectedly and causing injury to the patient or operators.
- Incorrect positioning of hands when loading the stretcher into the ambulance or when lowering it to the ground could result in them being crushed between the frame and the wheel bracket. Hands should only hold the tubular frame after operating the controls.
- Failure to read and understand the product instructions may result in injury to the patient and operators.

8. TECHNICAL DATA AND COMPONENTS

Note: Spencer Italia S.r.l. reserves the right to make changes to specifications without notice.

No.	DESCRIPTION	MATERIAL
1	Backrest release control	Nylon
2	Folding loading trolley	Stainless steel
3	Folding trolley unlocking control	Steel/PU
4	Loading trolley wheels	PU
5	Intermediate height unlocking lever	Nylon
6	Twist system unlocking lever	Nylon
7	Front leg piston	Steel
8	Front leg	Stainless steel
9	Twist bracket with brake	Nylon
10	Front wheel	Steel/PU
11	Rear wheel	PU

No.	DESCRIPTION	MATERIAL
12	Self-aligning rear bracket with brake	Steel
13	Rear leg	Steel
14	Rear leg locking release lever	Nylon/Steel
15	Fifth wheel	Nylon/Steel
16	Front leg locking release lever	Nylon/Steel
17	Variable height unlocking handle	Nylon/Steel
18	Trendelenburg/Fowler lever	Nylon/Steel
19	Patient bed footrest	Nylon/HDPE
20	Sidebar unlocking lever	TPU
21	Trendelenburg-Fowler positioning support strap	TPU
22	Sidebar	Aluminum/PA

IT

EN

DE

FR

ES

PT

FEATURE	DIMENSION
Length [mm]	1975
Width [mm]	575
Weight (2) [kg]	56.5
Load capacity (1) [kg]	260
Wheels [mm]	200
Loading trolley underwheel height [cm]	60 / 68 / 75, depending on the model
Frame material	Stainless steel
Lunghezza sponde [mm]	680
Altezza sponde [mm]	190

(1) Maximum load capacity means the total weight distributed according to human anatomy. When determining the total weight load on the product, the operator should consider the weight of the patient, equipment and accessories.

(2) The weight may vary by no more than ± 1 kg depending on the model.

9. COMMISSIONING

For first use, check that:

- Packaging is intact and has protected the device during transportation
- Check that all parts included in the packing list are present
- General functionality of the device
- The medical vehicle is equipped with a Spencer fixing system dedicated to the stretcher
- The stretcher floor is properly levelled
- The stretcher floor is wide and long enough to accommodate the stretcher and its accessories without becoming an obstacle
- The wheel of the front leg is between at least 5 and 10 cm from the ground during vehicle loading and unloading procedures in order to allow for the safe opening and locking of the front leg - see image in point 11.6.
-  A ground clearance of less than 5 cm for the front leg wheel during loading may prevent the full extension of the stretcher's front leg, potentially causing the stretcher to fall to the ground when unloaded from the vehicle. Clearances greater than 10 cm may make loading significantly more difficult for operators and cause severe jolts to the stretcher during unloading, potentially resulting in damage to the stretcher and harm to the patient.
- The fastenings must keep the stretcher attached to the structure of the vehicle.
- Do not modify the stretcher structural, leverage or traction parts in any way as this could cause damage to the patient and/or the rescuers.

 **Failure to take the above measures will preclude safe use of the device, resulting in risk of damage to the patient, operators and the device itself.**

To facilitate insertion of the stretcher into the ambulance, eliminate sharp edges on the edge of the ambulance loading floor. The stretcher should be secured using Spencer hooks in such a way as to prevent any movement during transport in the ambulance, even in difficult driving conditions. Practice with a stretcher without a patient before putting the device into service. **For subsequent use, perform the operations specified in paragraph 12.**


If the above conditions are met, the device may be considered ready for use; otherwise, you must immediately remove the device from service and contact the Manufacturer.

Do not alter or modify the device arbitrarily, as doing so could result in unpredictable operation and damage to the patient or rescuers and will void the warranty and release the Manufacturer from all liability.

 **MEASURING OF THE LOADING FLOOR MUST BE CARRIED OUT AFTER HAVING SET THE AMBULANCE ON AN EVEN SURFACE WITH NO BUMPS AND/OR UNEVENNESS AND AFTER ACCOMMODATING 2 (TWO) OPERATORS AT THE REAR OF THE SANITARY UNIT, WITH A SIMULATED LOAD OF UP TO 250 KG.**

If the vehicle is equipped with pneumatic or hydraulic suspensions, the loading height must be adjusted taking into account the worst-case conditions and/or the working conditions foreseen by the outfitter. Any problems regarding the use and/or safety risks associated with this system are not the responsibility of the manufacturer.

10. OPERATING CHARACTERISTICS

ELEMENT	DESCRIPTION	FUNCTION
A	Backrest adjustment unlocking control	Adjusts the backrest inclination chosen from one of the pre-set positions.
B	Folding trolley unlocking control	When engaged, unlocks the loading trolley, allowing it to be lowered to facilitate stretcher movement in confined spaces.

C	Variable height unlocking lever	When engaged, unlocks front leg movement, allowing the stretcher to be positioned at intermediate heights.
D	Twist system unlocking lever	When engaged, it pivots the front wheels. This control is present on both sides of the stretcher.
E	Brake	Located on the front and rear wheel brackets; brakes the wheels when activated
F	Rear leg unlocking handle (green)	When engaged, it unlocks the movement of the rear leg, allowing it to be closed.
G	Front leg unlocking (red) and fixing system release handle	When engaged, it unlocks the movement of the front leg, allowing it to be folded. When used with a fixing system, the lever also unlocks the stretcher from the fixing.
H	Variable height unlocking handle	When engaged, unlocks front leg movement, allowing the stretcher to be positioned at intermediate heights.
I	Trendelenburg/Fowler selection lever	Allows the foot side patient bed to be positioned, resulting in the Trendelenburg and Fowler position.
L	Sidebar unlocking levers	When engaged, unlocks the sidebar and opens it. It is not necessary to operate any controls to close the sidebar again, simply return it to its original position.

11. PROPER USE

Primary medical evaluations must be carried out before the patient is moved, lifted or transported. Once the diagnosis has been made, it is preferable to advise the patient to actively contribute to the transition from bed to stretcher/chair, while at the same time informing them of the risks they may be running. Before loading the patient, bring the device as close as possible to the patient.

11.1 EMERGENCY VEHICLE REQUIREMENTS

The stretcher is designed to enter and exit the medical compartment of an ambulance. The vehicle requirements should be:

- Level stretcher floor
- Stretcher floor that is wide and long enough to accommodate the stretcher without becoming an obstacle


The wheels of the front legs must be at a safe distance of between at least 5 and 10 cm from the ground during loading/unloading or when the loading trolley is still on the vehicle, to allow the front leg to open safely and safe use of the stretcher.

Failure to take the above measure will preclude safe use of the device, resulting in risk of damage to the patient, operators and the device itself.

11.2 PARKING BRAKES

To engage the parking brakes, simply press on the strips on the wheel brackets with your foot. To disengage the parking brakes, simply press down on the opposite side of the parking brake pedal and they will return to their original position with a slight click. Never leave the patient unattended, even if the parking brakes have been engaged.

11.3 SIDEBARS


The stretcher is equipped with sidebars, which are necessary for restraining the patient on the stretcher.

NEVER MOVE THE PATIENT AND NEVER LEAVE THE PATIENT ON THE STRETCHER WITHOUT FIRST RAISING THE SIDEBARS. Failure to heed this warning could cause serious damage.

To lower the sidebars, pull the levers located on the lower part of the machine, identified by the word **PULL** (letter G par.10), towards you simultaneously. The sidebar will open automatically.

To close the sidebar, lift it back to its original position and check that it is correctly hooked by making a few repeated pulling movements. Make sure that nothing interferes with the locking systems when closing. Sheets, for example, could prevent proper closure.

11.4 TWIST SYSTEM ACTIVATION

This system allows the front wheels to pivot. This rotation is unlocked by turning the control on the sides of the front leg.


When rotation is no longer necessary, reposition the control in the indicated locked position and push the stretcher firmly forwards; the wheels will now be locked again.

When using the twist system, at least one additional operator must be present to steer the head end of the stretcher, as the stretcher is much more difficult to manoeuvre with all four wheels unlocked.

Be sure to carefully assess the ground conditions before operating the twist system, as the stretcher may move in an unexpected and uncontrolled manner.

11.5 LOWERING THE STRETCHER

To facilitate patient transfer operations, it is advisable to bring the stretcher to one of its intermediate heights. There are 4 heights possible in addition to fully raised and fully lowered.

IT

EN

DE

FR



ES

PT

 It is the user's responsibility to evaluate which is the best position for transferring the patient in relation to the conditions of use.

Procedure for loading/unloading immobilized or hospitalized patients onto the stretcher at intermediate height:

- Engage the brakes on the rear wheels, taking care to position the wheels parallel to the stretcher. Engage the rear leg unlocking lever, hold it in place and lift the stretcher slightly to unlock the rear leg moving mechanism. Lower the stretcher about 5 cm, release the unlocking lever and support the stretcher until the first intermediate height position is reached. Make sure that the stretcher has reached a stable position. .
- To lower the front part of the stretcher, engage the front leg release lever (Letter G par.10) or pull the variable height unlocking handle towards you (Letter H par.10). The two controls perform the same function for this procedure. At the same time, if the patient is not present, apply a slight downward push to the front of the frame to open the front leg. If the patient is present, pushing is not necessary but it is necessary to support the weight of the stretcher by firmly grasping the front of the frame. **The weight of the stretcher, patient and any equipment attached to the stretcher must be supported.** After starting the lowering motion, release the handle while keeping your grip on the frame until the stretcher reaches the desired intermediate height position. Make sure that the position reached is stable.


The next heights can be reached by following the same procedure. Before engaging the unlocking controls, relieve some of the weight supported by the stretcher in order to ensure correct release of the movement mechanisms.

The device is equipped with an assistance system called "assisted lowering" during descent which reduces the load on the operator. This system is an aid for the operator and in no way replaces the support of the operator him/herself, who must therefore always be able and ready to support the load during all phases of handling, both during lowering and lifting as well as in case of unforeseen events.

Warning: the load on the operator increases gradually but not proportionally with each lowering step.

CAUTION: Do not place your hands on the sides of the loading trolley but only at the front to avoid hands being crushed between the loading trolley and the front leg.

When loading a non-hospitalized patient:


- 1 - Assist the patient as he/she climbs onto the stretcher, making sure that the first leg to be lifted rests on the footrest of the patient bed. When correctly positioned, help him/her raise the other leg.
- 2 - Immobilize the patient with the straps and raise the sidebars.
- 3 - Once the patient is lying on the stretcher and correctly immobilized, carefully lift the stretcher to a horizontal position.

Caution: The patient must always be seated, laid down or positioned by the operators with extreme care. The patient must therefore always be accompanied until completely laid down, thus avoiding impact on the patient's body which could jeopardise its integrity, as well as avoiding possible damage or aggravation of the patient's injuries.

11.6 LIFTING THE STRETCHER

To bring the stretcher back to the standard height from the lowered position or any of the intermediate positions, the operators must coordinate their movements to raise both the front and back of the stretcher simultaneously, ensuring the proper alignment of the patient bed. Do not carry out these operations until you have checked that the patient has been properly restrained with the strap and that the sidebars are raised.

- When lifting, firmly grip the stretcher frame with both hands.
- Warning: The leg unlocking levers must not be engaged during this procedure.
- Raise the assembly until the locking mechanisms are properly engaged.
- **Do not use the boards or other areas that are not designed for this purpose.**

11.7 LOADING AND UNLOADING THE STRETCHER FROM THE MEDICAL VEHICLE


⚠ Stretcher loading and unloading are among the most delicate steps during the use of stretchers. It is necessary to follow all the warnings in this manual very carefully, practising under simulated conditions before putting the product into service. The patient must always be properly immobilized. Always ensure that the conditions of use allow the front legs to open properly before any loading and/or unloading operations.

TO LOAD THE STRETCHER ONTO THE VEHICLE:

- The operator performing the loading should align with the stretcher and the fixing system so that both devices are in front of him/her and in the intended alignment for fixing.
- Move the wheels of the loading trolley forward into the vehicle until the front legs of the stretcher come into contact with the bodywork of the vehicle.
- In the event of a slide/ramp, the wheel of the loading trolley should be resting on the horizontal part of the bed before the front legs come into contact with the vehicle.
- Ensure that the wheels of the front legs are between at least 5 and 10cm from the ground. Failure to do so may result in serious damage and/or injury when unloading the stretcher from the vehicle.
- Engage the front leg unlocking lever while continuing to push the stretcher inside the vehicle until the rear legs are brought to a stop.
- Only after making sure that the rear legs are in contact with the bumper of the vehicle, engage the lever to release the rear legs, continuing to push towards the inside of the vehicle compartment.

⚠ During these operations, a part of the weight of the assembly will rest on the operator, who must therefore be able to support and accompany all device movements.


- Lock the stretcher on the medical vehicle using the Spencer stretcher hook installed on the vehicle.

⚠ Warning: DURING NORMAL STRETCHER MOVEMENTS, NEVER ENGAGE THE LEG LOCKING LEVERS, AS DOING SO COULD CAUSE THE STRETCHER TO FALL TO THE GROUND. THESE LEVERS MUST ONLY BE ENGAGED TO LOAD THE STRETCHER ON THE AMBULANCE OR TO BRING THE STRETCHER ONTO THE GROUND AS DESCRIBED ABOVE.

To unload the stretcher from the medical vehicle:

- If the stretcher is used with the S-Max Pro fixing system, it is released by pulling the lever on the hook as indicated by arrow 2 in the image. Immediately after engaging the levers, pull the stretcher slightly towards you and after moving the stretcher about 2 cm, place your hands back on the tubular frame.
- Pull the stretcher towards the outside of the vehicle, gripping the frame by its rear part. Support the weight of the stretcher until you feel the rear leg locking system engage. **As with the loading phase, the operator must be able to support the weight of the device.**
- Do not grab the footrest or any other area not intended for handling as this may cause damage to the operator, the patient or the device.
- **⚠** Pull the stretcher outwards until the front legs are fully open. **Do not lower the loading trolley from the loading surface until you have checked that the front leg locking system is properly engaged.**
- Finish extracting the stretcher from the vehicle.

IT

EN

DE


FR

ES

PT


11.8 TRENDLENBURG POSITION

To position the foot end of the patient bed in the Trendelenburg position, lift the patient bed by gripping it near the footrest so that the bolt slides into the designated groove until it reaches the upper housing. To return the bed to a horizontal position, slightly lift both levers on the sides of the bed, allowing the bolt to exit the housing and slide back to its initial position.

11.9 FOWLER POSITION

To position the patient bed in the Fowler position, grasp one of the two handles located on the sides of the bed and lift it. At the same time, press down on the patient bed near the footrest, ensuring the bolt slides into the groove and locks into the lower housing of the lever. To return the bed to a horizontal position, slightly lift both levers on the sides of the bed, allowing the bolt to exit the housing and slide back to its initial position.

11.10 BACKREST ADJUSTMENT


Pull the knob of the backrest piston to raise the backrest until it automatically locks into the first position. Follow the same procedure to reach the next positions, always checking the correct engagement of the locking system.

Changing the backrest inclination from a vertical to horizontal position:

Support the backrest with one hand (to prevent unexpected movements) and take away the weight from the backrest by moving it upward.

At the same time, pull the knob on the backrest piston upward until the safety mechanism unlocks. Then, supporting it with the other hand, lower the backrest toward the horizontal position until you reach the desired height and then release the knob to the rest position.

To reach subsequent positions, lift and support the backrest. Always lift the weight of the backrest before pulling the knob. Incorrect procedure could permanently damage the backrest piston.


11.11 LOWERING THE LOADING TROLLEY

To facilitate housing of the stretcher in particularly narrow spaces (e.g. lifts), the maximum length of the stretcher can be reduced following the steps below:

- Raise the backrest of the stretcher to a vertical position.
- Hold the loading trolley steady with one hand and with the other unlock the lever by pressing downwards as shown in the image to the side.
- Slowly lower the loading trolley with one hand, bringing it to rest on the front legs.
- To bring the trolley back into a horizontal position, guide it upwards and, engaging the lever again, align the trolley with the dedicated stops. Then release the lever and check that the loading trolley is locked in the desired position.

-  **CAUTION:** Never place your hands near the joint points to avoid the risk of crushing.

11.12 S-WINGS APPLICATION

S-Wings bariatric supports can be mounted on the stretcher. These are an accessory that increases the surface area of the stretcher with the aim of offering greater support to obese patients. To apply the device, refer to the manual supplied with this accessory.

11.13 RESTRAINT STRAPS

To apply the two-piece restraint straps, identify a portion of the frame at chest height and one at the patient's leg level.


The restraint straps should be applied in a position that allows for adequate immobilization of the patient.

Secure to the frame by means of a slipknot.

The restraint straps supplied with the stretcher must be fastened immediately after positioning the patient. To do so, insert the male buckle into the female buckle.

Strap positioning

To apply the thoracic belt (A), position the upper attachment points on the backrest and the lower attachment points on the patient bed, as shown in the figure. For the femoral straps (B) and ankle straps (C), position the attachment points on the lower part of the patient bed, as shown in the figure.

The male buckle must be passed through the appropriate slots in the metal ends of the chest straps.

Tension the straps after attaching them. When the patient is to be moved from the stretcher at the end of use, release the straps by pressing the release button marked "PRESS". Then open the stretcher sidebar on the side towards which you want to lower the patient from the stretcher.

To apply them:

- 1 – Remove the strap from the attachment point on the re-winding system.
- 2 – Insert the DNA system, keeping the belt taut, into the special cavity on the body of the bed. The male and female straps must be on the left and right of the platform, respectively.
- 3 – Put the belts back into the hole on the DNA re-winding system and complete the installation.
- 4 – Proceed in the same way for thorax systems.

After positioning the patient, fasten the straps and check that they are correctly applied.

 **DNA straps are not equipped with a pre-tensioner. Therefore, they must be manually adjusted after fastening, tensioning the straps and bringing the re-winding system to the end of the stroke.** Incorrect adjustment makes the device unsafe and ineffective.

12. CLEANING AND MAINTENANCE

Spencer Italia S.r.l. declines all responsibility for any direct or indirect damage which is the consequence of improper use of the product and spare parts and/or in any case of any repair carried out by a person other than the Manufacturer, who uses internal and external technicians authorised to do so; moreover, doing so will invalidate the warranty.

- The operator must wear suitable personal protective equipment, such as gloves, goggles, etc. during all checking, maintenance and cleaning operations.
- Establish a maintenance schedule, periodic inspections and extend the average life span, if foreseen by the Manufacturer in the User Manual, identifying a reference person who meets the basic requirements set forth in the User Manual.
- The frequency of inspections is determined by factors such as legal requirements, type of use, frequency of use, and environmental conditions during use and storage.
- Repairs of products manufactured by Spencer Italia S.r.l. must be carried out by the Manufacturer, who shall make use of specialised internal or external technicians

IT

EN

DE

FR

ES

PT

who, using original spare parts, shall provide quality repair service in strict compliance with the technical specifications indicated by the Manufacturer. Spencer Italia S.r.l. declines any responsibility for any direct or indirect damage resulting from improper use of spare parts and/or any repair work carried out by unauthorised parties.

- Use only original or Spencer Italia S.r.l. approved components/spare parts and/or accessories in order to carry out all operations without causing alterations or modifications to the product.
- All maintenance and overhaul activities must be recorded and documented with the relevant technical operation reports. This documentation must be kept for at least 10 years after the end of the product's life and must be made available to the competent authorities and/or the Manufacturer when requested.
- Cleaning, provided for reusable products, must be carried out in accordance with the Manufacturer's instructions in the User Manual in order to avoid the risk of cross-infection due to the presence of body fluids and/or residues.
- The product and all its components must be washed and left to dry completely before storage.
- If the product requires lubrication, this must be done after cleaning and complete drying.
- Document maintenance using the form in Annex B of this manual or other appropriate forms.

12.1 CLEANING

Failure to carry out the correct cleaning operations could increase the risk of cross-infection due to presence of body fluids and/or residues.
The operator must wear suitable personal protective equipment, such as gloves, goggles, etc. during all checking and cleaning operations.

Metal parts exposed to external agents undergo surface treatments and/or coating in order to obtain better resistance. Clean the exposed parts with water and delicate soap. **Never use solvents or stain removers. Do not use cleaning agents containing sodium hypochlorite as corrosion of components may occur.**

Rinse thoroughly with lukewarm water, making sure you have removed all traces of soap, which may deteriorate or compromise conditions and durability. **Avoid using high pressure water**, as it penetrates the joints and removes lubricant, creating the risk of corrosion on components. Let dry completely before storing. Drying after washing or after use in a wet environment must be natural and not forced. Do not use flames or other direct heat sources.

If disinfecting, use products that do not have a solvent or corrosive action on materials constituting the device. Be sure to take all precautions to ensure that there is no risk of cross-infection or contamination of patients and operators..

12.2 ROUTINE MAINTENANCE

Establish a maintenance program and periodic checks, identifying a designated reference person. The person entrusted with device maintenance must ensure the basic requirements envisaged by the manufacturer in the following paragraphs. All routine and special maintenance activities and all general overhauls must be recorded and documented with the relevant technical intervention reports. This documentation must be kept for at least 10 years after the end of the device's life and must be made available to the competent authorities and/or the Manufacturer when requested.

Routine device maintenance must be carried out by operators with specific qualifications, training and instruction in the use and maintenance of the device.

The operator must wear suitable personal protective equipment, such as gloves, goggles, etc. during all checking, maintenance and cleaning operations.

Checks to be carried out before and after each use and at deadline indicated above, are as follows:

- General functionality of the device
- Cleanliness of the device (remember that the failure of cleaning may cause the risk of cross infections)
- Correct tightening of all nuts, bolts and screws
- Absence of cuts, holes, tears on the structure, including the straps
- None of the tubes or metal sheets present bends or cracks
- Welded areas are intact, without any cracks or breaks
- Moving parts, wheels, levers, handles are intact and work properly
- Lubrication of moving parts
- Conditions of use of wheels and brake system
- The wheels are correctly fixed, are stable and turn properly
- The wheels are free of debris
- The device opens and locks properly
- The device opens and closes properly
- Functioning of springs
- The stretcher easily enters the ambulance
- The medical vehicle is equipped with a Spencer fixing system intended for the stretcher
- The coupling between the fixing system and the stretcher is suitable to guarantee the safety of the fixing .

The frequency of inspections is determined by factors such as legal requirements, type of use, frequency of use, and environmental conditions during use and storage.

Please note that cleaning must be carried out as described in this manual and that functionality must be verified before and after each use. Spencer Italia S.r.l. declines any responsibility for improper functioning or damages caused to the patient or user by the use of devices not subject to routine maintenance and will void the warranty and compliance to Medical Device Directive 2017/745/EU.

Use only original or Spencer Italia S.r.l. approved components/replacement parts and/or accessories to carry out any operation without causing alterations or modifications to the device. Otherwise, we decline all responsibility regarding incorrect operation or any damage caused by the device to the patient or the operator, invalidating the warranty and invalidating compliance with Regulation 2017/745/EU Medical Devices.

12.3 PERIODIC OVERHAUL

The device must be overhauled every year by the Manufacturer, who uses internal and external technicians specialised and authorized by the Manufacturer itself.

In the absence of an above-mentioned overhaul, the device must be DECOMMISSIONED, as compliance with Regulation 2017/745/EU will cease and, despite the CE marking, the device will no longer meet the safety requirements guaranteed by the Manufacturer at the time of supply.

Spencer Italia srl declines any responsibility for incorrect operation or damage caused by use of devices that have not been regularly serviced.

The mattress and straps shall be replaced every two years.

Only overhaul activities carried out by specialised technicians authorised by the Manufacturer are considered valid by Spencer Italia S.r.l.

12.4 SPECIAL MAINTENANCE

Special maintenance can only be carried out by the Manufacturer, who uses internal and external technicians specialised and authorised by the Manufacturer itself.

Only maintenance activities carried out by specialised technicians authorised by the Manufacturer are considered valid by Spencer Italia S.r.l..

The end-user can replace only the spare parts indicated in § 15.

12.5 LIFE SPAN

The device, if used as described in the following instructions, has a life span of 5 years from the date of purchase, which can be extended following annual overhauls.

Overhauls must be carried out by the Manufacturer, who uses internal and external technicians specialised and authorised by the Manufacturer. **In the absence of such annual overhauls, the device must be DISPOSED OF ACCORDING TO THE INFORMATION IN PARAGRAPH 16 AND NOTIFICATION MUST BE GIVEN TO THE MANUFACTURER.** The life span may be extended at the sole discretion of the Manufacturer or authorised centre if the safety requirements of the device are still met.

Spencer Italia S.r.l. will accept no responsibility for incorrect operation or damage caused by the use of devices that have not been overhauled by the Manufacturer or authorized centre, or that have exceeded the maximum allowable life span.

13. TROUBLESHOOTING TABLE

PROBLEM	CAUSE	REMEDY
The leg unlocking levers don't work properly or are difficult to operate.	The moving mechanisms are damaged.	Immediately remove the device from service and contact customer service.
	Means of connection between components have been lost.	
The stretcher is not correctly attached to the fixing system.	Wear or damage of the components that constitute the locking mechanisms.	Position the stretcher properly and check that the fifth wheel enters the dedicated part on the fixing system.
	The fifth wheel has not been properly inserted in the fixing system.	
Damage to the structure	Improper use	Immediately remove the device from service and contact customer service.
	The moving mechanisms are damaged.	
It is not possible to set the stretcher in intermediate position.	Something is interfering with the moving systems.	Make sure that nothing is interfering with the mechanisms.
	Levers have not been properly operated.	
The front legs do not lock when unloading from the medical vehicle.	The moving mechanisms are damaged.	Immediately remove the device from service and contact customer service.
	The height of the loading floor is not suited to the device, the safety height is not being respected.	
During positioning at the intermediate heights, no aid is provided by the system used for this purpose		Adjust the loading floor height to meet the requirements described in this manual. If the loading floor does not allow adjustments, immediately remove the device from service and contact the service centre.
	Damaged or incorrectly adjusted system	
		Take the device out of service and contact customer service

14. ACCESSORIES

ST42720	S-MAX PRO FIXING SYSTEM	CB09012	OXYPACK S - 2-LITRE CYLINDER HOLDER FOR STRETCHER IN PVC
IF01047	TRACK 4-30 TELESCOPIC IV HOLDER	ST00480	FIXO-KID - PAEDIATRIC RESTRAINT SYSTEM
EN90003C	END-T 10G INSTRUMENTS TABLE	CR90024	EXTENSION KIT FOR LOWERING SERIES STARK/AGILE
ST00491A	STX 90 TELESCOPIC HEADREST FOR STRETCHERS	CR90016	EXTENSION FOR MOVEMENT CROSS/STAR-K/AGILE

15. SPARE PARTS

ST00498	DNA STRAP - YELLOW CHEST STRAP WITH METAL FITTINGS	ST70000	QMX 777 ANATOMICAL HEAT-SEALED MATTRESS IN BLACK PVC
ST00497	DNA STRAP - YELLOW 2-PIECE STRAP WITH METAL FITTINGS	RIQM007	QMX 776-777 MATTRESS FASTENING STRAP

16. DISPOSAL

When devices and their accessories are no longer suitable for use, they can be disposed of as normal municipal solid waste if they have not been contaminated by special agents. Otherwise, follow the regulations in force regarding disposal.

Notice

The information contained in this document is subject to change without notice and is to be intended as a commitment by Spencer Italia S.r.l. subject to change. The images are included as examples and may vary slightly from the actual device.

© Copyright Spencer Italia S.r.l.

All rights reserved. No part of the document may be photocopied, reproduced or translated into another language without prior written consent from Spencer Italia S.r.l.

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