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
Front anchorage

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

Rear anchorage

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IT

Instructions for use
FIELD STRETCHER ANCHORAGES

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EN

1. MODELS	8
2. INTENDED USE	8
2.1 INTENDED USE	8
2.2 TARGET PATIENTS	8
2.3 PATIENT SELECTION CRITERIA	8
2.4 CONTRAINDICATIONS AND UNWANTED SIDE EFFECTS	8
2.5 USERS AND INSTALLERS	8
3. REFERENCE STANDARDS	8
4. INTRODUCTION	8
4.1 DEVICE LABELLING AND TRACEABILITY	8
4.2 SYMBOLS	8
4.3 WARRANTY AND SERVICE	9
5. WARNINGS/DANGERS	9
6. SPECIFIC WARNINGS	9
7. RESIDUAL RISK	9
8. TECHNICAL DATA AND COMPONENTS	9
9. PROPER USE	10
10. CLEANING AND MAINTENANCE	10
10.1 CLEANING	10
10.2 ROUTINE MAINTENANCE AND PERIODIC OVERHAUL	11
10.3 LIFE SPAN	11
11. TROUBLESHOOTING TABLE	11
12. ACCESSORIES	11
13. SPARE PARTS	11
14. DISPOSAL	11

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FR

1. MODELS

The following basic models may be subject to implementation or change without notice.

- FRONT FIELD STRETCHER ANCHORAGE
- REAR FIELD STRETCHER ANCHORAGE
- REAR FIELD STRETCHER ANCHORAGE 506

2. INTENDED USE

2.1 INTENDED USE

Fastenings for field stretchers are parking devices necessary to prevent unwanted movement of the stretcher. They should be installed in designated locations within the rescue area. Unsupervised parking is not permitted.

2.2 TARGET PATIENTS

The target patients are those for whom use of the field stretcher is foreseen.

2.3 PATIENT SELECTION CRITERIA

The selection criteria are those foreseen for use of the field stretcher.

2.4 CONTRAINDICATIONS AND UNWANTED SIDE EFFECTS

No particular contraindications or side effects are known with relation to use of the device, as long as it is used in accordance with the user manual.

2.5 USERS AND INSTALLERS

The intended users are those foreseen for use of field stretchers. These devices are not intended for lay people.

Operators must be able to provide the necessary patient care.

3. REFERENCE STANDARDS


REFERENCE	DOCUMENT TITLE
EU Regulation 2017/745	EU Regulation on Medical Devices

4. INTRODUCTION

4.1 DEVICE LABELLING AND TRACEABILITY

Each device is provided with a label, placed on the device itself and/or on the packaging, which contains the Manufacturer's identification data, product, CE marking, serial number (SN) or lot number (LOT). This must never be removed or covered.

EU Regulation 2017/745 requires manufacturers and distributors of medical devices to keep track of their location.

4.2 SYMBOLS

Symbol	Meaning	Symbol	Meaning
CE	Device in compliance with EU Regulation 2017/745	⚠	Danger – Indicates a hazardous situation that may result in a situation directly related to serious injury or death.
MD	Medical device	📖	See the user manual
🏭	Manufacturer	LOT	Lot number

	Date of manufacture		Serial number
	Caution: Federal law restricts the sale of this device by or on the order of a licensed professional (US market only)		Product code

4.3 WARRANTY AND SERVICE

Spencer Italia S.r.l. guarantees that products are free from defects for a period of **one year from the date of purchase**.
 Spencer Customer Service tel. +39 0521 541154, fax +39 0521 541222, e-mail service@spencer.it.
 Warranty and service conditions are available at www.spencer.it/en/support/terms-and-conditions

5. WARNINGS/DANGERS

PRODUCT FEATURES

Use of the product for any purpose other than that described in the User Manual is prohibited.

-  The product must not be tampered with or modified without the manufacturer's authorisation.
- Avoid contact with sharp or abrasive objects.
- Operating temperature: from -5°C to +50°C.
- Storage temperature: -10°C to +60°C.

GENERAL WARNINGS FOR MEDICAL DEVICES

- Do not use if the device or parts of it are punctured, damaged, or excessively worn.
- Do not alter or modify the device arbitrarily, as doing so could result in unpredictable operation and damage to the patient or rescuers and shall void the manufacturer's warranty and release the manufacturer from all liability.
- Participate in safety checks on products placed on the market, transmitting information regarding product risks to the Manufacturer as well as to the Competent Authorities for their respective actions.

6. SPECIFIC WARNINGS



To use the product, you must also have read, understood and carefully follow all the instructions in the Instructions for use.

-  Perform simulations of use with a stretcher and a patient simulating load and accessories before putting the device into service.
- Accumulations of dust can affect correct operation of the device.

7. RESIDUAL RISK

No residual risks, or rather risks that could arise despite compliance with all warnings in these instructions for use, have been identified.

8. TECHNICAL DATA AND COMPONENTS

FRONT ANCHORAGE		Component ID	Description
		1	Coupling screws Made of steel, they attach the right and left arm to the main anchorage.
		2/3	Guide arm Made of steel; present on the left (2) and right (3) side of the anchorage. They are equipped with holes for floor installation.
		4	Anchorage Made of steel and attached to the two guide arms, this is the locking element in the front direction for the stretcher.
		4	Unlocking lever Made of cast iron; unlocks sliding of the sliding carriage.
REAR ANCHORAGES		Component ID	Description
		1	Locking lever Made of cast iron; when placed vertically, this is the locking element in the rear direction for the stretcher. It is supplied in two different heights.
		2/3	Locking lever for 506 Made of cast iron; performs the same function as the one described above but is specific to the 506 stretcher.
		4	Sliding carriage Made of cast iron; adjusts the position of the locking lever
		4	Unlocking lever Made of cast iron; unlocks sliding of the sliding carriage.
		5	Base Made of cast iron; equipped with 4 holes used for fixing to the floor.

	FRONT ANCHORAGE	REAR ANCHORAGE
Dimensions	455 x 430 x 67 ± 5mm	131 x 115 x 161 (high lever) / 114 (low lever) / 74 (lever for 506)
Material	Steel	Cast iron


IT
EN
FR

9. PROPER USE

Find a level area where you can place the device to ensure its stability.
Set the stretcher in the desired position.
Place the front fastener against the front anchorage crossbar of the stretcher.
Mark the position of the holes using the front fastener as a template. Then drill holes in the floor and install the front stop by screwing it into the ground.


The front fastener should be attached using the following hardware:

ID	SPECIFICATIONS
A	M8 Hex head screw
B	ø8mm washer
C	ø8x32mm washer
D	M8 self-locking nut
E	M6 Hex head screw
F	ø6mm washer
G	ø6x24mm washer
H	M6 self-locking nut

NOTE: The length of the screws must be assessed by the installer in relation to the surface on which the product is to be screwed.

Raise the locking lever of the rear fastener to the vertical position and push it to the fully closed position.


Position the rear fastener against the rear crossbar of the stretcher, making sure that the stretcher is up against the previously installed front fastener.
After checking that the stretcher is firmly anchored, mark the position of the holes using the fastener itself as a template.
Then drill holes in the floor and install the stretcher by screwing the fastener into the ground.

The front fastener must be secured using M8 TCEI screws, ø8mm wide-band washers and M8 low self-locking nuts.

NOTE: The length of the screws must be assessed by the installer in relation to the surface on which the product is to be screwed.

To unlock the rear fastener, press and hold the unlocking lever, then pull the locking lever outwards until it comes out of the guide in the base, allowing the locking lever to be brought to a horizontal position.

This will allow the stretcher to be pulled out of the anchor position.

10. CLEANING AND MAINTENANCE

10.1 CLEANING

Failure to carry out the correct cleaning operations could increase the risk of cross-infection due to presence of body fluids and/or residues.
The operator must wear suitable personal protective equipment, such as gloves, goggles, etc. during all checking and cleaning operations.

Clean the exposed parts with water and delicate soap. **Never use solvents or stain removers.**
Do not use cleaning agents containing chlorine / sodium hypochlorite as corrosion of components may occur.
Rinse thoroughly with lukewarm water, making sure you have removed all traces of soap, which may deteriorate or compromise conditions and durability. **Avoid using high**

pressure water, as it penetrates the joints and removes lubricant, creating the risk of corrosion on components. Let dry completely before storing. Drying after washing or after use in a wet environment must be natural and not forced. Do not use flames or other direct heat sources.
If **disinfecting**, use products that do not have a solvent or corrosive action on materials constituting the device. Be sure to take all precautions to ensure that there is no risk of cross-infection or contamination of patients and operators.

10.2 ROUTINE MAINTENANCE AND PERIODIC OVERHAUL

Routine maintenance is the responsibility of the user, who must carry out checks before and after each use and in any case at least every 3 months:

- General functionality of the device
 - Cleanliness of the device (remember that the failure of cleaning may cause the risk of cross infections)
 - Fulfilment of the requirements of the instructions for use in section Warnings and Specific Warnings.
 - Fulfilment of the requirements of the manual in section on Proper use.
- No periodic overhaul is foreseen for the device.

10.3 LIFE SPAN

The device, if used as described in the following instructions, has a life span of 5 years from the date of purchase, which can be extended following annual overhauls.

Spencer Italia S.r.l. will accept no responsibility for incorrect operation or damage caused by the use of devices that have not been overhauled by the Manufacturer or authorized centre, or that have exceeded the maximum allowable life span.

11. TROUBLESHOOTING TABLE

PROBLEM	CAUSE	REMEDY
The unlocking lever does not activate.	The unlocking lever has not been activated.	Press the unlocking lever all the way down.
	Possible spring breakage.	Check the conditions of the spring. If it is damaged, remove the device from service and contact the manufacturer.
	The unlocking lever is broken.	Remove the device from service and contact the manufacturer.

If a problem or fault is detected that does not correspond to the above, please contact Spencer Italia srl customer care service.

12. ACCESSORIES

There are no accessories for these devices.

13. SPARE PARTS

There are no spare parts available for these devices.

14. DISPOSAL

When devices and their accessories are no longer suitable for use, they can be disposed of as normal municipal solid waste if they have not been contaminated by special agents. Otherwise, follow the regulations in force regarding disposal.

Warning

The information contained in this document is subject to change without notice and is to be intended as a commitment by Spencer Italia S.r.l. subject to change.
The images are included as examples and may vary slightly from the actual device.

© Copyright Spencer Italia S.r.l.

All rights reserved. No part of the document may be photocopied, reproduced or translated into another language without prior written consent from Spencer Italia S.r.l.

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