**Istruzioni d'uso**  
**CARRERA INCUBATOR**  
**Barella autocaricante**

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## CONTENTS

---

<b>1.</b>	<b>MODELS .....</b>	<b>19</b>
<b>2.</b>	<b>INTENDED USE.....</b>	<b>19</b>
<b>3.</b>	<b>REFERENCE STANDARDS.....</b>	<b>19</b>
<b>4.</b>	<b>INTRODUCTION .....</b>	<b>20</b>
4.1	USING THE MANUAL .....	20
4.2	DEVICE LABELLING AND TRACEABILITY .....	20
4.3	SYMBOLS .....	21
4.4	WARRANTY AND SERVICE.....	21
<b>5.</b>	<b>WARNINGS .....</b>	<b>21</b>
<b>6.</b>	<b>SPECIFIC WARNINGS.....</b>	<b>22</b>
6.1	PHYSICAL REQUIREMENTS FOR OPERATORS .....	23
<b>7.</b>	<b>RESIDUAL RISK.....</b>	<b>23</b>
<b>8.</b>	<b>TECHNICAL DATA AND COMPONENTS .....</b>	<b>24</b>
<b>9.</b>	<b>COMMISSIONING .....</b>	<b>25</b>
<b>10.</b>	<b>OPERATING CHARACTERISTICS .....</b>	<b>26</b>
<b>11.</b>	<b>PROPER USE .....</b>	<b>26</b>
11.1	EMERGENCY VEHICLE REQUIREMENTS .....	26
11.2	INCUBATOR INSTALLATION .....	26
11.3	PARKING BRAKES.....	27
11.4	LOWERING THE STRETCHER .....	27
11.5	LIFTING THE STRETCHER.....	28
11.6	LOADING AND UNLOADING THE STRETCHER FROM THE MEDICAL VEHICLE .....	28
<b>12.</b>	<b>CLEANING AND MAINTENANCE .....</b>	<b>29</b>
12.1	CLEANING .....	29
12.2	ROUTINE MAINTENANCE .....	30
12.3	PERIODIC OVERHAUL.....	30
12.4	SPECIAL MAINTENANCE .....	30
12.5	LIFE SPAN.....	30
<b>13.</b>	<b>TROUBLESHOOTING TABLE .....</b>	<b>31</b>
<b>14.</b>	<b>ACCESSORIES .....</b>	<b>31</b>
<b>15.</b>	<b>SPARE PARTS .....</b>	<b>31</b>
<b>16.</b>	<b>DISPOSAL .....</b>	<b>31</b>

## 1. MODELS

The following basic models may be subject to implementation or change without notice.

- LOW CARRERA INCUBATOR – ROLL-IN STRETCHER
- HIGH CARRERA INCUBATOR – ROLL-IN STRETCHER
- CARRERA INCUBATOR 63 – ROLL-IN STRETCHER
- CARRERA INCUBATOR 68 – ROLL-IN STRETCHER

## 2. INTENDED USE

Incubator stretchers are trolleys for transporting incubators in ambulances. Certified 10g when used with dedicated fastening systems in relation to the coupling between the fastening system and the trolley. Compliance with the requirements of EN 1789 for fastening between the incubator and the trolley is the responsibility of the installer.

### TARGET PATIENTS

The intended patients are those for whom transport with an incubator is necessary.

### PATIENT SELECTION CRITERIA

The intended patients are those for whom transport with an incubator is necessary.

### CONTRAINDICATIONS AND UNWANTED SIDE EFFECTS

No particular contraindications or side effects are known with relation to use of the device, as long as it is used in accordance with the user manual.

### USERS AND INSTALLERS

The intended users are persons trained in first aid procedures and in use of medical equipment in an emergency medical service (EMS) environment. Possible users also include emergency vehicle outfitters who may use the product before commissioning or during maintenance of the vehicle on which the stretcher is used.

### USER TRAINING

*Note: despite all efforts, laboratory tests, trials, and instructions for use, standards do not always reproduce practice, so the results obtained under actual conditions of product use in the natural environment may sometimes differ significantly.*

The best instructions are the continuous practice of use under the supervision of competent and trained personnel.

- Regardless of your level of experience with similar devices in the past, you should carefully read and understand the contents of this manual before installing, operating, or servicing this product. In case of any questions, please contact Spencer Italia S.r.l. for the necessary clarifications.
- The product must be used only by personnel trained in the use of this product and not on other similar products.
- The suitability of the users for use of this product can be attested by the training registration, in which trained persons, trainers, date and place are specified. This documentation must be kept for at least 10 years after the end of the product's life and must be made available to the competent authorities and/or the Manufacturer when requested. In the absence of such documentation, the relevant bodies will apply any foreseen sanctions.
- Do not allow untrained persons to help while using the product, as they may cause injury to themselves or others.
- The product must be put into use only by personnel trained in the use of this product and not on other similar products.

*Note: Spencer Italia S.r.l. is always available for training courses.*

### INSTALLER TRAINING

Device installation must be carried out by qualified personnel who have been trained and certified by Spencer Italia S.r.l.

The installer must strictly adhere to these instructions as well as to the state of the art for vehicle installations.

## 3. REFERENCE STANDARDS

As Distributor or End-User of the products manufactured and/or marketed by Spencer Italia S.r.l., users are strictly required to be familiar with the legal provisions in force in the country of destination of the goods, applicable to the devices to be supplied (including regulations relating to technical specifications and/or safety requirements) and, therefore, to understand the requirements necessary to ensure compliance of the products themselves with all legal requirements of the territory.

REFERENCE	DOCUMENT TITLE
EU REGULATION 2017/745	EU REGULATION 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL OF 5 April 2017 on medical devices
UNI EN ISO 1865-1	Patient handling equipment used in road ambulances - Part 1: General stretcher systems and patient handling equipment
UNI EN 1789	Medical vehicles and their equipment - Road ambulances

## 4. INTRODUCTION

### 4.1 USING THE MANUAL

The purpose of this manual is to provide healthcare professionals with the information necessary for safe and appropriate use and maintenance of the device.


*Note: The Manual is an integral part of the device and therefore it must be kept for the entire life of the device and must accompany it in any changes of use or ownership. If any instructions for use for products other than the one received are present, please contact the Manufacturer immediately before use.*

The instructions for use of Spencer products can be viewed and downloaded from [User Manuals Archives | Spencer](#) on the Spencer website or you can contact the Manufacturer. Exceptions are those items whose essentiality and reasonable and predictable use are such that it is not necessary to draw up instructions, in addition to the following warnings and indications on the label.

Regardless of your level of experience with similar devices in the past, it is advisable to carefully read and understand the contents of this manual before installing, operating, or servicing this product.

### 4.2 DEVICE LABELLING AND TRACEABILITY

Each device is provided with a label, placed on the device itself and/or on the packaging, which contains the Manufacturer's identification data, product, CE marking, serial number (SN) or lot number (LOT). **This must never be removed or covered.**

In the event of damage or removal, request a duplicate from the Manufacturer, or else the warranty will be void as the device can no longer be traced.

Regulation 2017/745/EU requires manufacturers and distributors of medical devices to keep track of their location. If the device is in a location other than the address to which it was shipped or sold, or if it was donated, lost, stolen, exported or destroyed, permanently removed from use, or if the device was not delivered directly by Spencer Italia S.r.l., inform Customer Service (see § 4.4).

### 4.3 SYMBOLS

SYMBOL	MEANING						
	Device in compliance with EU Regulation 2017/745						
	Medical device						
	Manufacturer						
	Date of manufacture						
	See the user manual						
	Serial number						
	Product code						
	Danger – Indicates a hazardous situation that may result in a situation directly related to serious injury or death.						
	Unique Device Identifier Production identification Alphanumeric code that identifies the production units of the device, composed of:						
	<table border="1"> <tbody> <tr> <td>(01)0805771123XXXX</td> <td>Company prefix (XXXX=GTIN)</td> </tr> <tr> <td>(11)YYMMDD</td> <td>Date of manufacture</td> </tr> <tr> <td>(NN) 1234567890</td> <td>(NN)1234567890 NN=10 =&gt;LOT/ NN=21=&gt;SN</td> </tr> </tbody> </table>	(01)0805771123XXXX	Company prefix (XXXX=GTIN)	(11)YYMMDD	Date of manufacture	(NN) 1234567890	(NN)1234567890 NN=10 =>LOT/ NN=21=>SN
(01)0805771123XXXX	Company prefix (XXXX=GTIN)						
(11)YYMMDD	Date of manufacture						
(NN) 1234567890	(NN)1234567890 NN=10 =>LOT/ NN=21=>SN						

### 4.4 WARRANTY AND SERVICE

Spencer Italia S.r.l. guarantees that products are free from defects for a period of **one year from the date of purchase**.

For information regarding correct interpretation of the instructions for use, maintenance, installation or return, please contact Spencer Customer Service tel. +39 0521 541111, fax +39 0521 541222, e-mail [service@spencer.it](mailto:service@spencer.it).

To facilitate service, always indicate the lot number (LOT) or serial number (SN) on the label attached to the package or device itself.

Warranty and service conditions are available at [General Terms and Conditions | Spencer](#).

To ensure product traceability and protect the maintenance and assistance procedures for your devices, Spencer has made the SPENCER SERVICE portal available to you at [Report an Issue | After Sales Support | Spencer](#).

*Note: Record and keep with these instructions: lot (LOT) or serial number (SN), if present, place and date of purchase, date of first use, date of checks, user name and comments.*

### 5. WARNINGS

Warnings, notes, and other important safety information are provided in this section and are clearly visible throughout the manual.

At least every 6 months, it is important to check for updated instructions and any changes involving your product. This information is freely available on the website [User Manuals Archives | Spencer](#)

#### PRODUCT FEATURES

Use of the product for any purpose other than that described in the instructions for use is prohibited.

- Before each use, always check the conditions of the product, as specified in the instructions for use. In the event of faults/damage that could compromise its functionality/safety, immediately remove it from service and contact the Manufacturer.
- If the product is found to be malfunctioning, immediately use a similar device to ensure continuity of ongoing operations.
- The product must not be tampered with (modification, tweaking, additions, repair). Should these operations be performed, we decline any responsibility for incorrect operation or any damage caused by the product itself; moreover, the CE certification (when required by law) and the product warranty shall be null and void.
- Anyone who modifies or causes to be modified, or intervenes on products manufactured by Spencer Italia S.r.l. in such a way as to alter their function or compromise their intended performance, must comply with the conditions laid down for the first placing on the market.
- When using the devices, position and adjust them in such a way that they do not hinder operator works or the use of any other equipment.
- Be sure to take every precaution to avoid hazards from contact with blood or body secretions, if applicable.
- Always comply with the maximum capacity indicated in the instructions for use. Maximum load capacity means the total weight distributed according to human anatomy. When determining the total weight load on the product, the operator should consider

the weight of the patient, equipment and accessories. Moreover, the operator should assess whether the overall size of the patient reduces the functionality of the product.

- Make sure that operators are in good physical condition before lifting, as listed in the instructions for use.
- The maximum weight, which weighs on each operator, must comply with local health and safety requirements.
- Avoid contact with sharp objects.
- Device installation must be carried out by qualified personnel trained and enabled by Spencer Italia S.r.l. The times and methods of such shall be agreed upon between the customer and our Sales Offices.
- Operating temperature: from -10°C to + 50°C.

#### STORAGE

- The product must not be exposed or come into contact with thermal sources of combustion or flammable agents, but must instead be stored in a dry, cool place, away from light and sun.
- Do not store the product under other more or less heavy materials that may damage the structure of the product.
- Store and transport the product with its original packaging, otherwise the warranty shall be invalidated.
- Storage temperature: -20°C to +60°C.

#### REGULATORY REQUIREMENTS

As Distributor or End-User of the products manufactured and/or marketed by Spencer Italia S.r.l., users are strictly required to be familiar with the legal provisions in force in the country of destination of the goods, applicable to the devices to be supplied (including regulations relating to technical specifications and/or safety requirements) and, therefore, to understand the requirements necessary to ensure compliance of the products themselves with all legal requirements of the territory.

- Promptly and in detail notify Spencer Italia S.r.l. (already in the quotation request phase) about possible fulfilments by the Manufacturer necessary for the compliance of products with specific legal requirements of the territory (including those deriving from regulations and/or regulatory provisions of another nature).
- Act with due care and diligence to help ensure compliance with the general safety requirements of the devices placed on the market, providing end-users with all the information necessary to carry out periodic revisions on the supplied devices, exactly as indicated in the instructions for use.
- Participate in safety checks on products placed on the market, transmitting information regarding product risks to the Manufacturer as well as to the Competent Authorities for their respective actions.
- Without prejudice to the above, the Distributor or End-User shall assume wider liability related to non-compliance with non-fulfilment of the above-mentioned obligations, with consequent obligation to indemnify and/or hold Spencer Italia S.r.l. harmless from any possible injurious effect.
- With reference to EU Regulation 2017/745, please note that public or private operators who, when exercising their activity, detect an incident involving a medical product are required to notify the Ministry of Health, within the terms and in the manner established by one or more ministerial decrees, and notify the Manufacturer. Public or private health care professionals are required to notify the Manufacturer of any other incident that may allow the adoption of measures to ensure the protection and health of patients and users.

#### GENERAL WARNINGS FOR MEDICAL DEVICES

- The application of the device should not last longer than the time required for first aid operations and subsequent transport to the nearest rescue point.
- Qualified personnel and at least two operators must be present during use of the device.
- Follow the internal procedures and protocols approved by your organisation.
- Disinfection operations must be carried out in accordance with the validated cycle parameters, as stated in the specific technical standards. Inadequate disinfection procedures may lead to risk of cross-infection.

#### 6. SPECIFIC WARNINGS

- Follow approved Emergency Medical Service procedures for patient positioning and transportation.
- Use by untrained personnel may result in injury to the rescuer or third parties.
- Do not use if the device or parts of it are punctured, torn, frayed, or excessively worn.
- Before any type of handling, make sure that operators have a secure grip on the device.
- Do not drag the device on uneven surfaces.
- Do not lift with a crane or other mechanical lift.
- Do not use drying machines.
- The device is a transport device and cannot be used as a stationing device.
- Do not use with devices other than those expressly approved by the manufacturer.
- Practice with a device without a patient to make sure you are familiar with the manoeuvres.
- For loading techniques, for operations on steep terrain or in special and unusual circumstances, the presence of multiple operators is recommended (not just two as expected in standard conditions).

- Make sure that nothing interferes with any stretcher handling or control mechanisms.
- Do not move the device if the weight is not properly distributed.
- Use only the perimeter frame to move the stretcher and not other points not intended for this purpose.
- Avoid excessive force when loading the stretcher into the ambulance: unnecessary force can cause damage and can adversely affect the fastening system and the functioning of the stretcher itself.
- Hold the device firmly if there is a patient in the incubator.
- Parking brakes are aids for the operator but they do not in any way replace his or her supervision.
- Pay the utmost attention to any obstacles (water, ice, debris, etc.) present on the route, as they may cause the operator to lose balance and compromise proper functioning of the device. If you cannot clear the way, choose an alternative route.
- Condensation, water, ice and accumulations of dust can affect correct operation of the device, making it unpredictable and causing a sudden change in the weight the operators have to support.
- For height differences greater than 10mm, the device must be lifted, taking care to grasp it from the structure and not from the sidebars/beds or other points not intended for this purpose.
- Once the wheels of the loading trolley have been placed on the ambulance floor, the wheels of the front leg must be between 5 cm and 10 cm from the ground, so that the front leg can be opened and locked securely. Check the height of the ambulance bed after each use; if it is altered, the trolley must be adjusted immediately by the manufacturer or by a specialised technician authorised by the manufacturer. If not, the manufacturer cannot be held responsible for correct functioning of the device or any damage caused by it.
- Partial opening of the legs may result in the device falling. Make sure that the legs are properly locked before any handling and that the pistons are fully extended and stable.
- Accidental operation of the variable height system, if present on the model, may result in the stretcher falling, causing injury to the patient and/or operators. Make sure that the unlocking handle is not accidentally activated.
- The operators should prepare to fully load the stretcher with the patient and accessories before operating the variable height unlocking control, if present on the model. Operating this control without applying sufficient support force will result in the stretcher falling unexpectedly and causing injury to the patient or operators.
- If the vehicle is equipped with pneumatic or hydraulic suspensions, the loading height must be adjusted taking into account the worst-case conditions and/or the working conditions foreseen by the outfitter.
- Any problems regarding the use and/or safety risks associated with this system are not the responsibility of the manufacturer.
- Inadequate installation of the loading floor can lead to buckling and subsequent damage to the welds of the front legs.
- Improper installation of the loading floor can lead to abnormal operation of the device and cause harm to the patient and user.
- Do not alter or modify the stretcher to adapt it to the emergency vehicle: doing so could result in unpredictable operation and damage to the patient or rescuers and shall void the manufacturer's warranty and release the Manufacturer from all liability.
- The product can only comply with regulation EN 1789 if used with the dedicated fixing system. The use of fasteners not approved by the Manufacturer is therefore prohibited. Non-approved fixing systems may alter the structural and functional characteristics of the device. Compliance with EN 1789 is limited to the interaction between the trolley and its fixing system. The assessment and verification of conformity of the entire incubator system, i.e. the trolley complete with incubator and its accessories, is the responsibility of the installer of such systems.
- If the stretcher is not locked on the fixing system or if it is not correctly positioned, this can result in dangerous movements, especially in the event of strong medical vehicle deceleration, which can cause damage to the patient and the operators. Always make sure that the locking system is properly engaged.
- Do not place limbs and/or objects between the legs and the frame, in the vicinity of the leg lifting pistons or in general between moving parts as this may result in crushing injuries.

### 6.1 PHYSICAL REQUIREMENTS FOR OPERATORS

The device is intended for professional use only. All operators must be trained to transport patients safely and efficiently. Do not allow untrained persons to help while using the product, as they may cause injury to themselves or others.

Operators using the device should be physically able to use the device and have good muscle coordination, as well as strong back, arms, and legs for lifting and holding and be able to grasp the device firmly with both hands.

Operators must be able to provide the necessary patient care. Users must be able to safely lift and handle the weight of the stretcher and patient together and any other equipment used with the device. The ability of each operator must be assessed before the definition of rescuer roles in use of the stretcher.

## 7. RESIDUAL RISK

No residual risks have been identified, i.e. risks that could arise despite compliance with all the warnings contained in these instructions for use.


## 8. TECHNICAL DATA AND COMPONENTS

Note: Spencer Italia S.r.l. reserves the right to make changes to specifications without notice.

Definition of height under wheel depending on model:

Model	Height under wheel [cm]
Low Carrera incubator	50
High Carrera incubator	60
Carrera incubator 63	63
Carrera incubator 68	68

Components:

No.	Description	Material
1	Perimeter frame	Steel
2	Rear leg unlocking handle (green)	Nylon
3	Front leg unlocking handle (red)	Nylon
4	Rear wheel brake	Iron
5	Rear wheel with bracket	PU
6	Rear piston	Steel
7	Rear leg	Steel
8	Front wheel	PU
9	Front piston	Steel
10	Loading trolley wheels	PU
11	Front crossbar for incubator fixing	Steel
12	Rear crossbar for incubator fixing	Steel

Feature	Dimension
Length [mm]	1970
Width [mm]	570
Wheel diameter [mm]	Ø 200 ± 5%
Load capacity [kg]	150
Weight <sup>without standard accessories</sup> [kg]	32

<sup>2</sup>The weight may vary by no more than ±2 kg depending on the model.

## 9. COMMISSIONING


For first use, check that:

- Packaging is intact and has protected the device during transportation
- Check that all parts included in the packing list are present
- General functionality of the device
- The medical vehicle is equipped with a Spencer fixing system dedicated to the stretcher
- The stretcher floor is properly levelled
- The stretcher floor is wide and long enough to accommodate the stretcher and its accessories without becoming an obstacle
- During loading and unloading from the vehicle, the front leg wheel must be between 5 cm and 10 cm from the ground to allow the front leg to open and lock safely - see image in point 11.6.
- The fastenings must keep the stretcher attached to the structure of the vehicle
- Do not modify the stretcher structural, leverage or traction parts in any way as this could cause damage to the patient and/or the rescuers.
- Failure to take the above measures will preclude safe use of the device, resulting in risk of damage to the patient, operators and the device itself.
- To facilitate insertion of the stretcher into the ambulance, eliminate sharp edges on the edge of the ambulance loading floor. The stretcher should be secured using Spencer hooks in such a way as to prevent any movement during transport in the ambulance, even in difficult driving conditions. Practice with a stretcher without a patient before putting the device into service.
- For subsequent use, perform the operations specified in paragraph 12.
- If the above conditions are met, the device may be considered ready for use; otherwise, you must immediately remove the device from service and contact the Manufacturer.

Do not alter or modify the device arbitrarily, as doing so could result in unpredictable operation and damage to the patient or rescuers and will void the warranty and release the Manufacturer from all liability.

## 10. OPERATING CHARACTERISTICS

Functions:

Element	Description	Function
A	Rear leg unlocking handle (green)	When engaged, it unlocks the movement of the front leg, allowing it to be closed.
B	Front leg unlocking handle (red)	When engaged, it unlocks the movement of the front leg, allowing it to be folded.
C	Rear wheel brake	When engaged, it brakes the rear wheels.
D	Crossbars for incubator fixing	Crossbars for fixing the incubator on the roll-in stretcher.

## 11. PROPER USE

Primary medical evaluations must be carried out before the patient is moved, lifted or transported. Once the diagnosis has been made, it is preferable to advise the patient to actively contribute to the transition from bed to stretcher/chair, while at the same time informing them of the risks they may be running. Before loading the patient, bring the device as close as possible to the patient.

### 11.1 EMERGENCY VEHICLE REQUIREMENTS

The stretcher is designed to enter and exit the medical compartment of an ambulance. The vehicle requirements should be:

- Level stretcher floor
- Stretcher floor that is wide and long enough to accommodate the stretcher without becoming an obstacle

The wheels of the front legs must be at a safe distance of at between 5 cm and 10 cm included from the ground during loading/unloading or when the loading trolley is still on the vehicle, to allow the front leg to open safely.

Failure to take the above measure will preclude safe use of the device, resulting in risk of damage to the patient, operators and the device itself.

### 11.2 INCUBATOR INSTALLATION


Spencer incubator trolleys are equipped with dedicated supports for installation of the main incubators on the market.

These supports are shown in red in the image below and the distance between them can be adjusted by loosening the screws located in the area shown by the arrows.


Once the distance between the supports has been changed, make sure that the screws have been tightened correctly.

To install the incubator or dedicated interfaces, drill the plates according to incubator manufacturer instructions.

The incubator installer is solely responsible for its correct fastening.  
 The supports can be installed at a distance of between 35 cm and 120 cm, included, from each other.  
 The supports have a usable drilling surface of 40x3 cm (red region).

### 11.3 PARKING BRAKES

To engage the parking brakes, simply press on the strips on the rear wheel brackets with your foot. To disengage the parking brakes, simply press down on the opposite side of the parking brake pedal and they will return to their original position with a slight click.


Never leave the patient unattended, even if the parking brakes have been engaged.

### 11.4 LOWERING THE STRETCHER


Lowering the stretcher to the ground  
 This procedure allows the stretcher to reach the minimum distance from the ground. This manoeuvre cannot be carried out when an incubator is installed on the device.

1. Raise the stretcher from the side of the feet until the loading trolley wheels rest on the ground.
2. With the loading trolley wheels resting on the ground and ready to bear the weight of the stretcher, engage both leg unlocking levers, push towards the stretcher to unlock the front leg moving system and then the stretcher towards the ground. The stretcher is now in the lowest position.

*Note: The parking brakes will not perform their function when the stretcher is fully lowered. Make sure that the stretcher is held in place by at least one operator.*

## 11.5 LIFTING THE STRETCHER

To bring back the stretcher back to the standard height, starting from any of the previous configurations, the operators must coordinate while simultaneously lifting the front and rear of the stretcher, ensuring the proper alignment of the patient bed. Only perform these operations after checking that the incubator has been correctly positioned and secured.

1. At the rear, grasp the part of the frame near the leg movement handles.
  2. At the front, similarly grab the frame above the loading trolley.
  3. Raise the assembly until the locking mechanisms are properly engaged.
- Only use the frame for lifting.
  - Do not use other areas that are not designed for this purpose.


Loading the incubator:

- 1) Make sure that the incubator is ready for transport, with all internal devices correctly secured and support systems functioning. Bring the incubator close to the stretcher and place it carefully.
- 2) Complete the positioning of the incubator and secure it using the appropriate anchoring systems.
- 3) Once the incubator has been correctly positioned and secured, carefully lift the stretcher, returning it to a horizontal and stable position to ensure safety during transport.

## 11.6 LOADING AND UNLOADING THE STRETCHER FROM THE MEDICAL VEHICLE

Stretcher loading and unloading are among the most delicate steps during the use of stretchers. It is necessary to follow all the warnings in this manual very carefully, practising under simulated conditions before putting the product into service. The patient must always be properly immobilized.

Always ensure that the conditions of use allow the front legs to open properly before any loading and/or unloading operations.

To load the stretcher onto the vehicle:

The operator performing the loading should align with the stretcher and the fixing system so that both devices are in front of him/her and in the intended alignment for fixing.

1. Move the wheels of the loading trolley forward into the vehicle until the front legs of the stretcher come into contact with the bodywork of the vehicle.
2. In the event of a slide/ramp, the wheel of the loading trolley should be in a horizontal position before the front legs come into contact with the vehicle bodywork.
3. Ensure that the wheels of the front legs are between at least 5 and 10cm from the ground. Failure to do so may result in serious damage and/or injury when unloading the stretcher from the vehicle.
4. Engage the right red lever (letter B par.10) to unlock movement of the front legs while continuing to push the stretcher inside the vehicle until the rear legs are brought to a stop.
5. Only after making sure that the rear legs are in contact with the bumper of the vehicle, engage the green lever (letter A par.10) to release the rear legs, continuing to push towards the inside of the vehicle compartment.

During these operations, a part of the weight of the assembly will rest on the operator, who must therefore be able to support and accompany all device movements.

- Lock the stretcher on the medical vehicle using the Spencer fixing system installed on the vehicle.
- **WARNING: DURING NORMAL STRETCHER MOVEMENTS, NEVER ENGAGE THE LEG LOCKING LEVERS, AS DOING SO COULD CAUSE THE STRETCHER TO FALL TO THE GROUND. THESE LEVERS MUST ONLY BE ENGAGED TO LOAD THE STRETCHER ON THE AMBULANCE OR TO REACH THE INTERMEDIATE HEIGHT AS DESCRIBED ABOVE.**

To unload the stretcher from the medical vehicle:

1. Free the stretcher from its S-MAX fastening system by operating the dedicated lever on the fastening system.
  2. Pull the stretcher towards the outside of the vehicle, gripping the rear part of the frame near the leg movement handles. Support the weight of the stretcher until you feel the rear leg locking system engage. As with the loading phase, the operator must be able to support the weight of the device.
- Do not grab any other area not intended for handling as this may cause damage to the operator, the patient or the device.
  - Pull the stretcher outwards until the front legs are fully open. Do not lower the loading trolley from the loading surface until you have checked that the front leg locking system is properly engaged.
  - Finish extracting the stretcher from the vehicle.

## **12. CLEANING AND MAINTENANCE**

Spencer Italia S.r.l. declines all responsibility for any direct or indirect damage which is the consequence of improper use of the product and spare parts and/or in any case of any repair carried out by a person other than the Manufacturer, who uses internal and external technicians authorised; to do so; moreover, doing so will invalidate the warranty.

- The operator must wear suitable personal protective equipment, such as gloves, goggles, etc. during all checking, maintenance and cleaning operations.
- Establish a maintenance schedule, periodic inspections and extend the average life span, if foreseen by the Manufacturer in the instructions for use, identifying a reference person who meets the basic requirements set forth in the instructions for use.
- The frequency of inspections is determined by factors such as legal requirements, type of use, frequency of use, and environmental conditions during use and storage.
- Repairs of products must be carried out by Spencer Italia S.r.l., who shall make use of specialised internal or external technicians who, using original spare parts, shall provide quality repair service in strict compliance with the technical specifications indicated by the Manufacturer. Spencer Italia S.r.l. declines any responsibility for any direct or indirect damage resulting from improper use of spare parts and/or any repair work carried out by unauthorised parties.
- Use only original or Spencer Italia S.r.l. approved components/spare parts and/or accessories in order to carry out all operations without causing alterations or modifications to the product. Otherwise, we decline all responsibility regarding incorrect operation or any damage caused by the device to the patient or the operator, invalidating the warranty and invalidating compliance with Regulation 2017/745/EU Medical Devices.
- All maintenance and overhaul activities must be recorded and documented with the relevant technical operation reports. This documentation must be kept for at least ten years after the end of the product's life and must be made available to the competent authorities and/or the Manufacturer when requested.
- Cleaning, provided for reusable products, must be carried out in accordance with the Manufacturer's instructions in the instructions for use in order to avoid the risk of cross-infection due to the presence of body fluids and/or residues.
- The product and all its components must be washed and left to dry completely before storage.
- If the product requires lubrication, this must be done after cleaning and complete drying.

### **12.1 CLEANING**

Metal parts exposed to external agents undergo surface treatments and/or coating in order to obtain better resistance. Clean the exposed parts with water and delicate soap. Never use solvents or stain removers.

Do not use cleaning agents containing sodium hypochlorite as corrosion of components may occur.

Rinse thoroughly with lukewarm water, making sure you have removed all traces of soap, which may deteriorate or compromise conditions and durability. Avoid using high pressure water, as it penetrates the joints and removes lubricant, creating the risk of corrosion on components. Let dry completely before storing. Drying after washing or after use in a wet environment must be natural and not forced. Do not use flames or other direct heat sources.

If disinfecting, use products that do not have a solvent or corrosive action on materials constituting the device. Be sure to take all precautions to ensure that there is no risk of cross-infection or contamination of patients and operators.

## 12.2 ROUTINE MAINTENANCE

Establish a maintenance program and periodic checks, identifying a designated reference person with specific qualifications, training and instruction in the use and maintenance of the device. The person entrusted with device maintenance must ensure the basic requirements envisaged by the manufacturer in the following paragraphs.

All routine and special maintenance activities and all general overhauls must be recorded and documented with the relevant technical intervention reports. This documentation must be kept for at least ten years after the end of the device's life and must be made available to the competent authorities and/or the Manufacturer when requested.

To ensure product traceability and protect the maintenance and assistance procedures for your devices, Spencer has made the SPENCER SERVICE portal available to you ([Request Assistance | Authorized Customer Service | Spencer](#)).

The operator must wear suitable personal protective equipment, such as gloves, goggles, etc. during all checking, maintenance and cleaning operations.

Checks to be carried out before and after each use and at deadline indicated above, are as follows:

- General functionality of the device
- Conditions of cleanliness of the device.
- Correct tightening of all nuts, bolts and screws.
- Absence of cuts, holes, tears or abrasions on the entire structure
- None of the tubes or metal sheets present bends or cracks.
- Welded areas are intact, without any cracks or breaks.
- Moving parts, wheels, levers, handles are intact and work properly
- Lubrication of moving parts.
- Conditions of use of wheels and brake system.
- The wheels are correctly fixed, are stable and turn properly
- The wheels are free of debris.
- The locking systems on the device work properly once it has been opened.
- The device opens and closes properly.
- The stretcher easily enters the ambulance.
- The medical vehicle is equipped with a Spencer fixing system intended for the stretcher
- The coupling between the fixing system and the stretcher is suitable to guarantee the safety of the fixing

The frequency of inspections is determined by factors such as legal requirements, type of use, frequency of use, and environmental conditions during use and storage.

Please note that cleaning must be carried out as described in this manual and that functionality must be verified before and after each use. Spencer Italia S.r.l. declines any responsibility for improper functioning or damages caused to the patient or user by the use of devices not subject to routine maintenance and will void the warranty and compliance to Medical Device Directive Regulation 2017/745/EU.

## 12.3 PERIODIC OVERHAUL

The device must be overhauled every year by the Manufacturer, who uses internal and external technicians specialised and authorized by the Manufacturer itself.

In the absence of an above-mentioned overhaul, the device must be DECOMMISSIONED, as compliance with Regulation 2017/745/EU will cease and, despite the CE marking, the device will no longer meet the safety requirements guaranteed by the Manufacturer at the time of supply. Spencer Italia srl declines any responsibility for incorrect operation or damage caused by use of devices that have not been regularly serviced. The front leg spring must be replaced every year.

Failure to replace it within the prescribed period could compromise the safety of the stretcher.

Only overhaul activities carried out by specialised technicians authorised by the Manufacturer are considered valid by Spencer Italia S.r.l.

## 12.4 SPECIAL MAINTENANCE

Special maintenance can only be carried out by the Manufacturer, who uses internal and external technicians specialised and authorised by the Manufacturer itself. Only maintenance activities carried out by specialised technicians authorised by the Manufacturer are considered valid by Spencer Italia S.r.l.. The end-user can replace only the spare parts indicated in § 15.

## 12.5 LIFE SPAN

The device, if used as described in the following instructions, has a life span of five years from the date of purchase, which can be extended following annual overhauls. Overhauls must be carried out by the Manufacturer, who uses internal and external technicians specialised and authorised by the Manufacturer. In the absence of such annual overhauls, the device must be DISPOSED OF ACCORDING TO THE INFORMATION IN PARAGRAPH 16 AND NOTIFICATION MUST BE GIVEN TO THE MANUFACTURER.

The life span may be extended at the sole discretion of the Manufacturer or authorised centre if the safety requirements of the device are still met.

Spencer Italia S.r.l. will accept no responsibility for incorrect operation or damage caused by the use of devices that have not been overhauled by the Manufacturer or authorized centre, or that have exceeded the maximum allowable life span.

### 13. TROUBLESHOOTING TABLE

PROBLEM	CAUSE	REMEDY
The leg unlocking levers don't work properly or are difficult to operate.	The moving mechanisms are damaged.	Immediately remove the device from service and contact customer service.
	Means of connection between components have been lost.	Immediately remove the device from service and contact customer service.
The stretcher is not correctly attached to the fixing system.	Wear or damage of the components that constitute the locking mechanisms.	Immediately remove the device from service and contact customer service.
	The fifth wheel has not been properly inserted in the fixing system.	Position the stretcher properly and check that the fifth wheel enters the dedicated part on the fixing system.
Damage to the structure	Improper use	Immediately remove the device from service and contact customer service.
The front legs do not lock when unloading from the medical vehicle.	The moving mechanisms are damaged.	Immediately remove the device from service and contact customer service.
	The height of the loading floor is not suited to the device, the safety height is not being respected.	Adjust the loading floor height to meet the requirements described in this manual. If the loading floor does not allow adjustments, immediately remove the device from service and contact the service centre.
When the stretcher release lever is operated, the stretcher does not move from the fixing system and remains attached.	The magnet controlling the release, located inside the fifth wheel holder, has moved or been lost.	Manually unlock the stretcher and finish the unloading operation. At the end of service, check the positioning of the magnet, then restore its original position. If the problem persists, immediately remove the device from service and contact the service centre.
	The Sensor Lock fixing system is not powered, is locked or is faulty.	Check the fixing system power supply. If the problem persists, immediately remove the device from service and contact the service centre.
	The stretcher is used with a fastener that does not provide automatic release.	If necessary, request a different fixing system.

### 14. ACCESSORIES

ST42707	S-MAX 10G FASTENER FOR STRETCHER
EN90003	END-T TRAY FOR STRETCHER 10G PATENTED

### 15. SPARE PARTS

There are no spare parts for this device.

### 16. DISPOSAL

When devices and their accessories are no longer suitable for use, they can be disposed of as normal municipal solid waste if they have not been contaminated by special agents. Otherwise, follow the regulations in force regarding disposal.

**Notice**

The information contained in this document is subject to change without notice and is to be intended as a commitment by Spencer Italia S.r.l. subject to change. The images are included as examples and may vary slightly from the actual device.

**© Copyright Spencer Italia S.r.l.**

All rights reserved. No part of the document may be photocopied, reproduced or translated into another language without prior written consent from Spencer Italia S.r.l.

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