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

© Copyright Spencer Italia s.r.l.

All rights reserved. No part of the document may be photocopied, reproduced or translated into another language without prior written consent from Spencer Italia s.r.l.

## Instructions for use ENERGY Termoregulated bag for infusion heating

The image is purely indicative and may not fully reflect the characteristics of the product

Spencer Italia s.r.l. – Via Provinciale n° 12 – 43038 Sala Baganza (PR) – Italy

First issue:  
Rev. 0.1  
Code

04/06/2025  
03/07/2025  
CCI5374

<b>INDICE</b>		
<b>1.</b>	<b>MODELS</b>	<b>18</b>
<b>2.</b>	<b>INTENDED USE</b>	<b>18</b>
2.1	INTENDED USE AND CLINICAL BENEFITS	18
2.2	TARGET PATIENTS	18
2.3	PATIENT SELECTION CRITERIA	18
2.4	CONTRAINDICATIONS AND SIDE EFFECTS	18
2.5	USERS AND INSTALLERS	18
2.6	USER TRAINING	18
2.7	INSTALLER TRAINING	18
<b>3.</b>	<b>REFERENCE STANDARDS</b>	<b>19</b>
<b>4.</b>	<b>INTRODUCTION</b>	<b>19</b>
4.1	USING THE MANUAL	19
4.2	DEVICE LABELLING AND TRACEABILITY	19
4.3	SYMBOLS	19
4.4	WARRANTY AND SERVICE	20
<b>5.</b>	<b>WARNINGS</b>	<b>20</b>
<b>6.</b>	<b>SPECIFIC WARNINGS</b>	<b>21</b>
<b>7.</b>	<b>RESIDUAL RISKS</b>	<b>22</b>
<b>8.</b>	<b>TECHNICAL DATA AND COMPONENTS</b>	<b>23</b>
<b>9.</b>	<b>COMMISSIONING AND INSTALLATION</b>	<b>24</b>
9.1	INSTALLATION	25
9.2	START-UP	27
<b>10.</b>	<b>OPERATING CHARACTERISTICS</b>	<b>28</b>
10.1	CONNECTION OF THE DEVICE	28
10.2	POSITIONING OF THE BAGS	28
10.3	USE	28
10.4	ALARMS AND FAILURE SIGNALS	29
<b>11.</b>	<b>CLEANING AND MAINTENANCE</b>	<b>29</b>
11.1	CLEANING	29
11.2	ROUTINE MAINTENANCE	29
11.3	PERIODIC OVERHAUL	30
11.4	SPECIAL MAINTENANCE	30
11.5	LIFE SPAN	30
<b>12.</b>	<b>TROUBLESHOOTING TABLE</b>	<b>30</b>
<b>13.</b>	<b>ACCESSORIES</b>	<b>30</b>
<b>14.</b>	<b>SPARE PARTS</b>	<b>31</b>
<b>15.</b>	<b>DISPOSAL</b>	<b>31</b>

## 1. MODELS

This document applies to the following references:

IF03030 ENERGY – TERMOREGULATED BAG FOR INFUSION HEATING

The above mentioned variants can be subject to implementation or change without notice

## 2. INTENDED USE

### 2.1 INTENDED USE AND CLINICAL BENEFITS

Thermoregulated bag for use in ambulances, for heating saline solutions during intravenous administration to the patient with the aim of reducing the effects linked to the infusion of a liquid at a temperature lower than body temperature.

The device heats the liquid to a temperature between 37 and 41 degrees.

The time required for heating depends on initial environmental conditions. Once heated, the device can be taken outside the ambulance.

The device is not intended to heat blood or plasma.

### 2.2 TARGET PATIENTS

There are no particular indications related to the patient group

### 2.3 PATIENT SELECTION CRITERIA

The device does not have a patient selection criteria. The selection depends on the need to carry out infusion with warmed solution which could be necessary for any patient.

Any selection is related only to the infusion set (tube and needle) to be used. The infusion set is not part of the device .

The device has the following indications:

- Infusion of warmed saline solution in case of hypothermia
- Infusion of warmed saline solution when the operator/physician evaluates that the infusion of colder liquids could result in patient discomfort or body temperature reduction

### 2.4 CONTRAINDICATIONS AND SIDE EFFECTS

There are no known contraindications.

### 2.5 USERS AND INSTALLERS

Intended users are doctors, paramedics or nurses trained for the use of the device and in infusion procedures.

Intended installers are typically ambulance builders.

Operators using the device should be physically able to use the device. Operators' ability must be assessed before the definition of roles in use of the device.

Operators must be able to provide the necessary patient care.

These devices are not intended for lay people.

Energy is a device intended for professional use only. Do not allow untrained persons to help while using the product as they may cause injury to themselves or others.

Despite all efforts, laboratory tests, trials, and instructions for use, standards do not always reproduce practice, so the results obtained under actual conditions of product use in the natural environment may sometimes differ significantly.

The best instructions are the continuous practice of use under the supervision of competent and trained personnel.

Operators using the device should be physically able to use the device. Operators' ability must be assessed before the definition of roles in use of the device.

Operators must be able to provide the necessary patient care.

### 2.6 USER TRAINING

Regardless of your level of experience with similar devices in the past, the contents of this manual must be carefully read and understood prior to installing, operating, or servicing this product. In case of any questions, please contact Spencer Italia S.r.l. for the necessary clarifications.

- The product must be used only by personnel trained specifically in the use of this product and not on other similar products.
- The suitability of the users for use of this product can be attested by the training registration, in which trained persons, trainers, date and place are specified. **This documentation must be kept for at least 10 years after the end of the product's life and must be made available to the competent authorities and/or the Manufacturer when requested. In the absence of such documentation, the relevant bodies will apply any foreseen sanctions.**
- Do not allow untrained people to help while using the product, as they may cause injury to themselves or others.
- The product must be put into use only by personnel trained in the use of this product and not on other similar products.
- Training must be registered for all personnel involved in the use of the device in all phases of its life cycle. This training must be repeated periodically. The training must be carried out by personnel with proven knowledge in the use and maintenance of the device.
- Note: Spencer Italia S.r.l. is always available for training courses.

### 2.7 INSTALLER TRAINING

The installer of the device must be able to ensure that all equipment, systems, containers and connections comply with all applicable safety requirements and standards. This implies knowledge of all applicable regulatory requirements and standards.

The installer should:

- Ensure proper installation of devices inside the ambulance
- Ensure that during installation phases the device has been not damaged and does not come in contact with substances that could compromise safety of patient or operators.

- Ensure the suitability of connections according to the performances provided for the device
- Perform all post-installation test procedures as specified in this manual

### 3. REFERENCE STANDARDS

As Distributor or End-User of the products manufactured and/or marketed by Spencer Italia S.r.l., users are strictly required to be familiar with the legal provisions in force in the country of destination of the goods, applicable to the devices to be supplied (including regulations relating to technical specifications and/or safety requirements) and, therefore, to understand the requirements necessary to ensure compliance of the products themselves with all legal requirements of the territory.

REFERENCE	DOCUMENT TITLE
EU Regulation 2017/745	EU Regulation on Medical Devices
IEC 60601-1+AMD1:2012 + AMD2:2020	Medical electrical equipment - Part 1: General requirements for basic safety and essential performance
IEC 60601-1+2:2014+AMD1:2020	Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard: Electromagnetic disturbances - Requirements and tests
IEC 60601-1-8:2006/A1:2012/A2:2020	Medical electrical equipment - Part 1-8: General requirements for basic safety and essential performance - Collateral Standard: General requirements, tests and guidance for alarm systems in medical electrical equipment and medical electrical systems
IEC 60601-1-12:2014/Amd 1:2020	Medical Electrical Equipment - Part 1-12: General requirements for basic safety and essential performance - Collateral Standard: Requirements for medical electrical equipment and medical electrical systems used in the emergency medical services environment
UNI EN 1789:2024	Medical vehicles and their equipment - Road ambulances

### 4. INTRODUCTION

#### 4.1 USING THE MANUAL

- The purpose of this manual is to provide healthcare professionals with the information necessary for safe and appropriate use and maintenance of the device.
- Note: The Manual is an integral part of the device and therefore it must be kept for the entire life of the device and must accompany it in any changes of use or ownership. If any instructions for use for products other than the one received are present, please contact the Manufacturer immediately before use.

Spencer products User Manuals can be downloaded from the site <http://www.spencer.it> or by contacting the Manufacturer. Exceptions are those items whose essentiality, reasonable and predictable use are such that it is not necessary to draw up instructions, in addition to the following warnings and indications on the label.

Regardless of your level of experience with similar devices in the past, it is advisable to carefully read and understand the contents of this manual before installing, operating, or servicing this product.

#### 4.2 DEVICE LABELLING AND TRACEABILITY

Each device is provided with a label, placed on the device itself and/or on the packaging, which contains the Manufacturer's identification data, product, CE marking, serial number (SN) or lot number (LOT). **This must never be removed or covered.**

EU Regulation 2017/745 requires manufacturers and distributors of medical devices to keep track of their location. If the device is in a location other than the address to which it was shipped or sold, or if it was donated, lost, stolen, exported or destroyed, permanently removed from use, or if the device was not delivered directly from Spencer Italia S.r.l., inform Customer Service (see § 4.4).

#### 4.3 SYMBOLS

Symbol	Meaning
	Device in compliance with EU Regulation 2017/745
	Medical device
	Manufacturer
	Date of manufacture
	Mandatory symbol "read instruction for use".
	Insulation class – Class II
	Applied part – BF applied part

	Direct current
	Heater active
	Status indicator – device powered
	Active alarm
	Failure/Alarm
<b>IP33</b>	Protection of enclosures for electrical devices First digit: protection against ingress of particulate greater than 2.5 mm diameter Second digit: Protected against Spraying water water with 60° angle
	Warning for the correct disposal of the product according to the European Directive 2012/19/UE
<b>LOT</b>	Lot number
<b>SN</b>	Serial number
<b>REF</b>	Product code

Production identification

Alphanumeric code that identifies the production units of the device, composed of:

<b>UDI</b>		(01)0805771123xxxx	company prefix (xxxx=GTIN)
		(11) 200626	date of production (YYMMDD)
		(NN) 1234567890	(NN)1234567890 NN=10 =>LOT/ NN=21=>SN

#### 4.4 WARRANTY AND SERVICE

Spencer Italia S.r.l. guarantees that products are free from defects for a period of **one year from the date of purchase**.

For information regarding correct interpretation of the instructions for use, maintenance, installation or return, please contact Spencer Customer Service tel. +39 0521 541154, fax +39 0521 541222, e-mail [service@spencer.it](mailto:service@spencer.it).

To facilitate service, always indicate the lot number (LOT) or serial number (SN) on the label attached to the package or device itself.

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To ensure product traceability and safeguard the maintenance and service procedures of your devices, Spencer has made the SPENCER SERVICE portal available to you (Request Assistance | Authorized Service Centers | Spencer).

*Note:*

*Record and keep with these instructions: lot (LOT) or serial number (SN), if present, place and date of purchase, date of first use, date of checks, user name and comments.*

## 5. WARNINGS

Warnings, danger, notes, and other important safety information are provided in this section and are clearly visible throughout the manual.

### Product features

- Use of the product for any purpose other than that described in the User Manual is prohibited.
- Before each use, the condition of the product should always be checked, as specified in the User Manual. In the event of faults/damage that could compromise its functionality/safety, remove immediately from service and contact the Manufacturer.
- If the product is found to be malfunctioning, use immediately a similar device to ensure continuity of ongoing operations. Non-compliant devices must be taken out of service.
- The product must not be tampered with or modified without the manufacturer's authorisation (modification, tweaking, additions, repair, use of non-approved accessories), as they may constitute imminent danger of injury to persons and material damage. Should these

operations be performed, we decline any responsibility for incorrect operation or any damage caused by the product itself; moreover, the CE marking and the product warranty shall be null and void.

- Be sure to take every precaution to avoid hazards from contact with blood or body secretions, if applicable.
- Avoid contact with sharp or abrasive objects.
- For the use of the device, it is recommended the adoption of an internal procedure or operational management which defines:
  - a) Control of documents and records
  - b) Training methods and related plans
  - c) Responsibilities relating to preventive and extraordinary maintenance
  - d) Definition of communication methods for personnel involved in the use and maintenance of the device.

#### Storage

- The product must not be exposed or come into contact with thermal sources of combustion or flammable agents and must instead be stored in a dry, cool place, away from light and sun.
- Do not store the product under other items no matter how heavy. They could damage the product.
- Store and transport the product with its original packaging, otherwise the warranty will be invalidated.

#### Regulatory requirements:

As Distributor or End-User of the products manufactured and/or marketed by Spencer Italia s.r.l., users are strictly required to be familiar with the legal provisions in force in the country of destination of the goods, applicable to the devices to be supplied (including regulations relating to technical specifications and/or safety requirements) and, therefore, to understand the requirements necessary to ensure compliance of the products themselves with all legal requirements of the territory.

- Promptly and in detail notify Spencer Italia s.r.l. (from the quotation request phase) about possible fulfilments by the Manufacturer necessary for the compliance of products with specific legal requirements of the territory (including those deriving from regulations and/or regulatory provisions of another nature).
- Act with due care and diligence to help ensure compliance with the general safety requirements of the devices placed on the market, providing end-users with all the information necessary to carry out periodic revisions on the supplied devices, exactly as indicated in the User Manual.
- Participate in safety checks on products placed on the market, transmitting information regarding product risks to the Manufacturer as well as to the Competent Authorities for their respective actions.
- Without prejudice to the above, the Distributor or End-User shall assume wider liability related to non-compliance with non-fulfilment of the above-mentioned obligations, with consequent obligation to indemnify and/or hold Spencer Italia s.r.l. free from any possible injurious effect.
- With reference to EU Regulation 2017/745, please note that public or private operators who, when exercising their activity, detect an incident involving a medical product are required to notify the Competent authority of the Member State in which the used and/or patient is established, within the terms and in the manner established by one or more ministerial decrees, and notify the Manufacturer. Public or private health care professionals are required to notify the Manufacturer of any other incident that may allow the adoption of measures to ensure the protection and health of patients and users.

#### General warnings for medical devices

The user must read carefully the following in addition to the general warnings.

- Qualified personnel must be present during use of the device.
- Do not use if the device or parts of it are punctured, torn, frayed, or excessively worn.
- Follow the internal procedures and protocols approved by your organisation.
- Do not alter or modify the device arbitrarily, as doing so could result in unpredictable operation and damage to the patient or rescuers and shall void the manufacturer's warranty and release the manufacturer from all liability.

## 6. SPECIFIC WARNINGS

To use the device, you have must also read, understood and carefully follow all the instructions in the user manual.

- The device is intended to be used on medical vehicles and not for home therapy.
- Do not use if the device is damaged or is in a poor state of cleaning.
- To be used only by trained staff.
- Do not wash or clean the device with pressurized water or air.
- Do not use drying machines.
- Condensation, water, ice and dust accumulation can affect the proper operation of the device making it dangerous for patient and operators.
- Regularly check the fixation of the device (where present); if altered or damaged is necessary to restore their safety conditions before the use of the device. Otherwise we assume no responsibility on the improper functioning or any damage caused by the device itself.
- In the event of a malfunction of the device, infusion therapy should be immediately restored by means of alternative devices able to carry out successfully and effectively such procedures
- Before each use, check the integrity of the device as specified in the user's manual. In case of malfunction or damage that may compromise the functioning and safety of the device, patient or operator, the device must be taken out of service and the manufacturer or a service center must be contacted.
- During use, the assistance of qualified personnel must be ensured.
- Do not leave the patient without the assistance of at least one physician or practitioner with clinical experience of the ongoing therapy when the device is in use.
- The device should not be exposed to or come into contact with any source of combustion or inflammable agents.
- It must be stored in a dry, cool place, protected from light and it must not be exposed to sunlight.
- The device shall be used in a ventilated environment.
- The use of the device in environmental and supply conditions other than those specified in this manual, will undermine the safety of operations and of the device itself.
- Do not lubricate any part of the device. It's not required by any kind of maintenance.


- Do not use the device if the safety conditions described in this manual are not met.
- Do not use the device if it has not been subject to scheduled maintenance or maintenance is required after normal use.
- Follow the maintenance program as specified by the manufacturer.
- Respect the maintenance and revision deadlines and the replacing of components described in this manual or provided by the manufacturer.
- The installation must ensure easy access to the device, visibility of all the information displayed on it, as well as audibility of the alarms
- The infusion tubes must be as linear as possible, avoiding bends that could reduce the flow rates required for the therapy
- Use only accessories approved by the manufacturer
- Do not use the device for warming blood, plasma, or medications. Only soft-bag saline solution containers equipped with a hanging hole may be warmed. The use of rigid plastic or glass containers may affect heating performance
- Portable and mobile RF equipment can affect the operation of the device
- The installation and placing of the device must consider what is described in tables of paragraph 9 in order to ensure that the device maintains its basic safety and performance.
- Use of cables or power supplies other than those approved by the manufacturer, can have adverse effect on the electromagnetic performance of the device.
- The use of RF equipment, including antennas, can adversely affect the device. The device or its supply cables should be kept away from other devices according to parameters listed in paragraph 9.
- The use of accessories other than those approved by the manufacturer, may result in increased electromagnetic emissions or reduction of the immunity level of the device.
- The device should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the device should be observed to verify normal operation in the configuration being used.
- Failure to follow warnings related to electromagnetic compatibility, can compromise the essential performance or basic safety of the device, effecting the proper operation of components of the software and cause unexpected changes or behavior or false indications.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) from any part of the device, including cables. Failure to observe this could degrade the performance of the equipment.
- In order to maintain the basic safety and performance of the device, it is essential to check before and after each use, the integrity of the electrical connections and to verify that the environment in which the device is installed has not been changed for example by installation of additional devices.
- Avoid contact with abrasive surfaces and/or edges that could compromise the integrity of the bag. Damage to the bag may cause heat leakage, resulting in reduced device efficiency.
- The user must not have impairments that prevent proper reading and interpretation of information and proper operation of controls.
- The installation must be performed ensuring appropriate distances between devices that could cause each other electromagnetic interference as specified in paragraph 9
- Ensure that the ambulance electrical installation is compliant with the requirements described in this user manual and as required in the standard EN 1789. Connection to a power supply not compliant with such requirements will affect the proper operation of the device.
- The use in presence of flammable and/or anesthetic gas may cause fire risks
- For patient safety, although not directly related to device safety, it is recommended to strictly follow hygienic, clinical, and safety practices during the application of the infusion needle and the connected systems.

## **7. RESIDUAL RISKS**

Ignoring the activation of alarm signals may pose a risk of haemolysis for the patient.


**8. TECHNICAL DATA AND COMPONENTS**

Note: Spencer Italia s.r.l. reserves the right to make changes to specifications without notice.

Id	Description
1	Transport handle – Made of polypropylene, it allows for carrying and attaching the device.
2	Cigarette lighter plug – Made of ABS, it is intended for connection to the power source represented by the cigarette lighter socket installed in the vehicle
3	Pass-through holes – By opening the closure flaps, they allow the infusion tubes connected to the bags housed inside the bag, to pass through.
4	Control unit – A nylon casing that houses the control electronics and from which the wiring for the heating element inside the bag and the power supply originates.
5	Hooks – Internal components of the bag designed for attaching the bag to be heated.

**Components of accessory fastener**

Id	Description
7	Buckle Strap – The polypropylene strap with buckle allows the device to be held during transport in an ambulance.
8	Fixing bracket – Made of steel, it is the support anchored in the vehicle and on which the device must be housed

## Technical data

<b>DIMENSIONAL</b>	
Height	335 ± 1mm
Width	250 ± 1mm
Thickness	142 ± 1mm
Power supply cable length	2 m
Capacity	2 bags of liquid of a volume of 500ml each
Fastener dimensions	286x90x356mm (LxPxH)
<b>Weight</b>	
Weight	1 ± 0,1kg
Weight of the fastener	2,42 ± 0,1kg
<b>Power supply</b>	
Voltage	12-15 Vdc (-15% + 25%)
Power consumption	Max 45W
<b>FUSES</b>	
Cigarette lighter plug – 7.5A integrated fuse – non-replaceable Control electronics casing: ø5x20 250V F 5A	
<b>Sensors and alarms</b>	
Sensors	1 Temperature sensor and protection thermistor
Alarm	High priority – Visual + acoustic If the heated chamber exceeds a temperature of 42°C
Alarm sound level	48db
Failures warning/medium priority alarm <b>Note:</b> Failure means a temporary or permanent failure of an element	<ul style="list-style-type: none"> <li>- Temperature sensor failure</li> <li>- Heater failure</li> <li>- Microcontroller failure</li> <li>- Firmware failure</li> </ul>
<b>Other</b>	
Chamber heating time	5min @Tamb 22°C, T0 liquid 22°C
<b>Classifications</b>	
Classification according to Reg. 2017/745/UE	Ila
Insulation Class according to IEC 60601-1	Class II
Applied part (Infusion tube)	The device has no applied part – the infusion tube must be considered a type BF applied part.
Enclosure protection according to IEC 60529	IP33 First digit: Protection against the ingress of particles with a diameter greater than 2.5 mm Second digit: Protected against water splashes at an angle of 60°
Operation	Continuous

### Environmental conditions

Storage and Transport Conditions:

- Temperature: from -40°C to +70°C
- Humidity: from 15% to 90%
- Atmospheric Pressure: from 620 hPa to 1060 hPa

Transient and Operating Conditions:

- Temperature: from 0°C to +40°C
- Humidity: from 15% to 90%
- Atmospheric Pressure: from 620 hPa to 1060 hPa

Environmental conditions significantly affect the heating capacity of the device.

The optimal operating range, within which the expected bag temperature is guaranteed, is from 15°C to 40°C. Below 15°C, the bag temperature will be higher than the ambient temperature, but the target temperature will not be reached.

## 9. COMMISSIONING AND INSTALLATION

For first use, check that:

- Packaging is intact and has protected the device during transportation
- All parts included in the packing list are present.
- General functionality of the device
- Product cleanliness
- Do not modify the device or its parts for any reason as this could cause injury to the patient and/or rescuers.
- Failure to take the above measures will preclude safe use of the device, resulting in risk of injury to the patient, operators and damage to the device itself.

If the above conditions are met, the device may be considered ready for use; otherwise, you must immediately remove the device from service and contact the Manufacturer.

Do not alter or modify the device arbitrarily, as doing so could result in unpredictable operation and damage to the patient or rescuers and will void the warranty and release the Manufacturer from all liability.

### 9.1 INSTALLATION

The surfaces on which the device is installed must ensure secure fastening under all circumstances.

For the installation of the fastening system, holes must be drilled into the support surface of the vehicle to allow proper anchoring.

The installer must make the necessary assessments to ensure that future maintenance of the device is straightforward and that the fastening system is installed in compliance with the requirements of EN 1789. It is recommended to use the fastening system as a drilling template, making sure that the fixture is correctly leveled.

The holes to be used are shown in the figure.

Below are the requirements for the operating environment:


The emergency vehicle must be equipped with:

- A suitable compartment for positioning the device.
- A power supply source compliant with the requirements described in these instructions.

The device can be used for the infusion procedure while housed in the fastening system.

To allow fluid flow, the fastening system must therefore be positioned at a height above that of the patient under normal ambulance transport conditions.

Always check that the fastening system's belt buckle is securely fastened to prevent damage

- The device is intended for use in the electromagnetic environment specified below. The customer, installer, and user must ensure that these conditions are met at all times.

Guide and Manufacturer's Declaration		
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device must ensure that it is used in such an environment.		
EMISSION TESTS	Conformity	Guide to the electromagnetic environment
Emissions in RF CISPR 11	Group 1	The device uses RF energy only for its internal functions. Its RF emissions are therefore very low and unlikely to cause any interference with electronic equipment nearby.  The device is suitable for use in all environments including domestic as well as those directly connected to a low-voltage public network source supply of which supplies buildings used for domestic purposes.
Emissions in RF CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Compliant	
Emissions as a result of voltage fluctuations / flicker-IEC 6100 3-3		

Guide and Manufacturer's Declaration		
The device is intended for use in the electromagnetic environment specified below. The customer and / or user of the device must ensure that the device is used in such environment.		
IMMUNITY TEST	Conformity level	Guide to the electromagnetic environment
Electrostatic discharge (ESD) IEC 61000-4-2	± 8kV at contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV in air	Floors should be wood, concrete or ceramic tile. If the floors are covered with synthetic material, the relative humidity should not exceed most 30%
Electrical fast transient/burst IEC 61000-4-4	± 2kV power supply ± 1kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Over voltage IEC 61000-4-5	± 0.5, 1kV for line-to-line surge 0.5, 1, 2kV for line to ground surge	Mains power quality should be that of a typical commercial or hospital environment
Voltage dips IEC 61000-4-11	0 % UT; 0,5 cycles at 0°, 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°  0% UT 1 cycle and 70% UT 25/30 cycles (25 at 50Hz and 30 at 60Hz)  Single phase to 0°	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible

Voltage interruptions IEC 61000-4-11	0% UT; 250/300 cycles	power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment
<i>Note UT is the value of the tension of power source</i>		

<b>Guide and Manufacturer's Declaration</b>		
The device is intended for use in the electromagnetic environment specified below. The customer and / or end user of the device must ensure that the equipment is used in such environment.		
IMMUNITY TEST	Level of conformity	Guide to the electromagnetic environment
Conducted RF IEC 61000-4-6	6 V 150kHz to 80MHz in ISM bands and amateur radio bands	The equipment for communication in portable and mobile radio-frequency (RF) should not be placed near any part of the appliance, including cables etc. and should be kept at a distance which should never be less than the recommended and calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 0,583x\sqrt{P}$
Radiated immunity CEI EN 61000-4-3	80% AM a 1kHz 10 V/m 80MHz to 2.7 Ghz	$d = 1,2x\sqrt{P}$ from 80MHz to 800MHz $d = 2,3x\sqrt{P}$ from 800MHz to 2,7 GHz  where P is the maximum rated power output of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey should be at less than the compliance level in each frequency range b. Interference may occur in the vicinity of equipment marked with the following symbol: 
NOTE 1 At 80 MHz and 800 MHz the separation distance for the range of higher frequency is applied NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.		
<sup>a</sup> The ISM (industrial, scientific and medical) bands between 150 kHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz.		
<sup>b</sup> The compliance levels in the ISM frequency bands between 150 kHz and 80 MHz and in the frequency range 80 MHz to 2,5 GHz are intended to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas. For this reason, an additional factor of 10/3 has been incorporated into the formulae used in calculating the recommended separation distance for transmitters in these frequency ranges.		
<sup>c</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device.		
<sup>d</sup> Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 6 V/m.		

<b>Recommended separation distances between portable and mobile communications equipment and the ventilator 190.</b>			
The device is intended for use in an electromagnetic environment in which radiated RF disturbance is controlled. The customer or the user of the device may prevent electromagnetic interference by maintaining a minimum distance between the communications equipment radio frequency (RF) Portable and mobile equipment (transmitters) and the device, as described below and in accordance with the maximum output power of the communication device			
Maximum output power rating of the transmitter	Separation distances according to frequency of transmitter (m)		
(W)	From 150 kHz to 80 Mhz Inside and outside ISM bands $d = 0,583x\sqrt{P}$	From 80 MHz to 800 Mhz $d = 1,2x\sqrt{P}$	From 800 Mhz to 2,7 Ghz $d = 2,3x\sqrt{P}$
0,01	0.058	0,12	0.23
0,1	0.184	0,38	0.73

1	0.583	1,2	2,3
10	1.844	3,8	7,3
100	5.83	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in metres (m) can be determined using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 The ISM (industrial, scientific and medical) bands between 150 kHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz.

NOTE 3 An additional factor of 10/3 has been incorporated into the formulae used in calculating the recommended separation distance for transmitters in the ISM frequency bands between 150 kHz and 80 MHz and in the frequency range 80 MHz to 2,5 GHz to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas.

NOTE 4 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Guide and Manufacturer's Declaration		
Immunity to proximity fields from RF wireless communications equipment		
Test Frequency (MHz)	Modulation	Immunity level (V/m)
385	Pulse modulation <sup>(1)</sup> at 18Hz	27
450	FM <sup>(2)</sup> ±5Hz deviation 1kHz sine	28
710	Pulse modulation <sup>(1)</sup> at 217Hz	9
745	Pulse modulation <sup>(1)</sup> at 217Hz	9
780	Pulse modulation <sup>(1)</sup> at 217Hz	9
810	Pulse modulation <sup>(1)</sup> at 18Hz	28
870	Pulse modulation <sup>(1)</sup> at 18Hz	28
930	Pulse modulation <sup>(1)</sup> at 18Hz	28
1720	Pulse modulation <sup>(1)</sup> at 217Hz	28
1845	Pulse modulation <sup>(1)</sup> at 217Hz	28
1970	Pulse modulation <sup>(1)</sup> at 217Hz	28
2450	Pulse modulation <sup>(1)</sup> at 217Hz	28
5240	Pulse modulation <sup>(1)</sup> at 217Hz	9
5500	Pulse modulation <sup>(1)</sup> at 217Hz	9
5785	Pulse modulation <sup>(1)</sup> at 217Hz	9

The carrier shall be modulated using a 50 % duty cycle square wave signal.  
As an alternative to FM modulation, 50 % pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.

- If the environmental conditions are appropriate, the device can be installed while verifying that:
- If the optional fastener is used, the installation surface is level and strong enough to withstand accelerations and vibrations to which the device could be subjected during use on the medical vehicle. It's suggested the use of a backplate.
- If the optional fastener is not used, the device shall be stored in a suitable compartment with surfaces sturdy enough to withstand the forces they may be subjected to. Otherwise, the safety and functionality of the device will be compromised.
- The power supply system has been regularly maintained or, if it is being commissioned for the first time, periodic maintenance has been scheduled.
- Distances between other electrical and electronic devices are as described in this manual.
- The positioning of the device does not cause any type of obstruction inside the medical vehicle.
- The power supply has the characteristics described in this manual.

The device must be connected to electrical power supply having the following characteristics:

Power supply	
Voltage	From 12 to 15 V === (-15% + 25%)
Current	> 4A
Connection	DIN ISO 4165

After evaluating the layout inside the ambulance, it is possible to proceed with the installation according to one of the following methods:

• **1 – Installing the fixing bracket**

Install the mounting bracket according to the instructions.

Use of the bracket is recommended.

• **2 – Installation without fixing bracket**

The bag must be placed in a vertical position, at a height above the patient and at a distance such as to ensure that the infusion tube can reach the patient without being put under tension.

**9.2 START-UP**

For a proper and safe use of the product, proceed as follows:

- Ensure that the power supply complies with the specifications in this manual.
- Ensure that all components are properly connected

Spencer is available for consultations relating to the use and control activities of the device.


If the device works properly and the conditions stated are met, the device can be considered ready for use; otherwise the device must be immediately put out of service and the manufacturer contacted.

**Do not alter or modify the device arbitrarily; modifications may cause unpredictable operation and damage to the patient or rescuers and will also void the warranty, relieving the manufacturer from any liability.**

## 10. OPERATING CHARACTERISTICS

### 10.1 CONNECTION OF THE DEVICE


The device is supplied with the cigarette lighter plug coiled inside the front pocket and connected to the control electronics housing.

The device is not equipped with power buttons; the heating element is activated as soon as the cigarette lighter plug is inserted into the cigarette lighter socket.

To connect, proceed as follows:

- 1 – Open the zipper of the front pocket
- 2 – Take out the cable with the cigarette lighter plug and unroll it to its full length, making sure that the locking ring on the control electronics housing is firmly screwed on.
- 3 – Connect the cigarette lighter plug to the cigarette lighter socket inside the ambulance
- 4 – Check that the following indicator lights turn on

 	<p>Green light – Status indicator showing that the device is connected to the power source. Always on when power is available. There is an additional indicator on the cigarette lighter plug.</p>
  	<p>Flashing orange light – Heater active, chamber temperature not yet within the operating range. Steady orange light – Heater active, chamber temperature within the operating range. Orange light off – Heater not active. This may occur if the ambient temperature is very high which will prevent the heater from switching on.</p>

### 10.2 POSITIONING OF THE BAGS

The device can hold two 500 ml infusion bags.


Before placing the infusion bags inside the bag, make sure they are intact, ensuring that no liquid can leak inside the bag.

**Warning:** Leaks of liquid inside the bag may cause short circuits.

To place the infusion bags, open the zipper of the largest compartment of the bag and hang the bag on one of the two hooks located at the top. Repeat the same procedure for the second bag

### 10.3 USE



- Make sure the device is correctly positioned as required in section 9.1. During use, it must always be positioned vertically. Insert the infusion bags as described in section 10.2. Connect the device to the power source as described in section 10.1.
- Wait for the chamber to heat up. The chamber heating time is significantly affected by external environmental conditions.
- At an ambient temperature of 22°C, the chamber takes approximately 5 minutes to reach the target temperature (between 38 and 41°C). The liquid temperature is not directly monitored; therefore, it is recommended to place pre-heated bags inside the device and connect it to the power supply before arriving at the intervention site.
- To begin therapy, attach the tubing with drip set to the bag according to standard infusion procedures. Open the closure flaps of the holes located on the bottom side of the bag and pull the infusion tube out until it reaches the patient.

- The physician in charge of the emergency intervention is responsible for choosing the device to be used and for making the necessary clinical evaluations to ensure its correct use throughout the procedure.
- When used with the mounting bracket, the device must be securely placed into the dedicated housing of the bracket, then fastened using the closure strap with buckle. Make sure the buckle is properly locked.

#### 10.4 ALARMS AND FAILURE SIGNALS

- In addition to the signals described in paragraph 10.1 which identify the device's activity status, there are two other light and acoustic indicators as described below:

	<p>High priority alarm: Flashing red light accompanied by an acoustic alarm</p> <p><b>Meaning:</b> The heated chamber has exceeded the temperature of 42°C.</p> <p><b>Action:</b> Stop the infusion. Disconnect the device from the power supply, remove the bags from the device, and continue therapy only after removing the bags. Record the environmental conditions of use. If the alarm persists or occurs frequently, contact a service center.</p>
	<p>The signal can identify both medium and low priority alarm</p> <p>Medium priority alarm: Flashing yellow light accompanied by an acoustic alarm</p> <p><b>Meaning:</b> A fault is present in the temperature sensor or in the heater.</p> <p><b>Action:</b> Stop the infusion. Disconnect the device from the power supply. Check that the device has not been damaged. The activation of this warning may indicate that the heated chamber is unable to reach its operating temperature.</p> <p>Low priority alarm: Steady yellow light</p> <p><b>Meaning:</b> Inadequate power supply voltage (below 10V).</p> <p><b>Action:</b> Stop the infusion. Disconnect the device from the power supply. Check the condition of the power source. If the voltage drops further, the device will automatically shut down.</p>

## 11. CLEANING AND MAINTENANCE

### 11.1 CLEANING

Failure to carry out the correct cleaning routine could increase the risk of cross infection, due to presence of body fluids and/or residuals.

- The operator must always wear adequate personal protection such as gloves and mask, etc. during all checking and cleaning procedures
- The described operations must be performed after each use of the device. Failure to carry out the correct cleaning routine could increase the risk of cross infection, due to presence of body fluids and/or residuals.
  - The operator must always wear adequate personal protection such as gloves and mask etc. during all checking and cleaning procedures.
  - Disconnect the device from the power supply
  - Clean the exposed surfaces with neutral detergent and dry with a cloth.

### 11.2 ROUTINE MAINTENANCE

Establish a maintenance program and periodic testing routine, identifying an employee responsible for this. The person to whom the ordinary maintenance of the device is entrusted must ensure the basic requirements foreseen by the manufacturer in following paragraphs are checked. All maintenance and periodic servicing activities must be registered and kept together with the servicing reports. These documents have to be kept for a period of 10 years after the disposal of the device itself. This register will be made available to the competent authorities and/or manufacturer if requested.

Checks to be carried out before and after each use and at deadline indicated above, are as follows:

- General functionality of the device
- Cleanliness of the device (remember that the failure to clean may cause the risk of cross infections)
- Absence of alarm signals.
- No structural part is deformed or compromised. Damage to the compartments could expose electrical components which could then be damaged and/or cause harm to patients and operators.


- Visual check of cables integrity.
- Verification of the integrity of the bag. Damage to the bag could compromise the warming capabilities of the bag.

If the device does not power on even after checking the power source, it may be necessary to replace the 5A protection fuse.

To do this, open the front pocket of the bag, unscrew the fuse holder cartridge, and remove the blown fuse.

After replacing the fuse, before using the device again, check the vehicle's electrical system and make sure that no device alarms are present.

Use only original or Spencer Italia s.r.l. approved components/replacement parts and/or accessories to carry out any operation without causing alterations or modifications to the device. Failure to do so, we decline all responsibility regarding incorrect operation or any damage caused by the device to the patient or the operator, invalidating the warranty and invalidating compliance with EU Regulation 2017/745.

### 11.3 PERIODIC OVERHAUL

The device does not require periodic overhaul.

### 11.4 SPECIAL MAINTENANCE

Special maintenance can only be carried out by the Manufacturer, who uses internal and external technicians specialised and authorised by the Manufacturer itself.

Only maintenance activities carried out by specialised technicians authorised by the Manufacturer are considered valid by Spencer Italia s.r.l.

### 11.5 LIFE SPAN

The device, if used as indicated in the following instruction manual, has a life span of 5 years from the purchase date

The service life may be extended, at the sole discretion of the Manufacturer or an authorized service center, if the safety requirements of the device are still guaranteed, for a maximum of an additional 5 years.

Spencer Italia s.r.l. will accept no responsibility for incorrect operation or damage caused by the use of devices that have not been overhauled by the Manufacturer or authorized centre, or that have exceeded the maximum allowable life span.

## 12. TROUBLESHOOTING TABLE

PROBLEM	CAUSE	SOLUTION
The device does not turn on	Wrong electrical connection	Check that the power source meets the specifications described in this manual.
	Improper power supply	Ensure that the power cable connector is properly screwed into the control electronics housing.
	Blown fuse	Check the supply voltage. Voltages below 9V will prevent the device from powering on.  Check the integrity of the 5A fuse.
Activation of a failure alarm	Heater or other components and/or functions failure	Disconnect the device from the power supply and reconnect it. If the problem persists, take the device out of service and contact the service center.
Heater always active	Possible failures of internal components of the device	Disconnect the device from the power supply and reconnect it. If the problem persists, take the device out of service and contact the service center.
	Device not properly closed, causing heat loss or inadequate environmental conditions	Close the bag zipper. Ensure that the environmental conditions for use are met.


## 13. ACCESSORIES

Device	Description
IF03034	FASTENER ENERGY THERMOREGULATED BAG

## 14. SPARE PARTS

Device	Description
RIIF004	2m power supply cable w/cigarette lighter plug

The power cable can be replaced by opening the front pocket of the bag and unscrewing the connector nut attached to the control electronics housing. Attach the new power cable, ensuring that the connector is fully inserted, and tighten the nut manually.

When contacting the Manufacturer or an authorized service center the device code and serial number must always be supplied.

## 15. DISPOSAL

When devices and their accessories are no longer suitable for use, they can be disposed of as normal municipal solid waste if they have not been contaminated by special agents. Otherwise, follow the regulations in force regarding disposal.

### WARNINGS FOR THE CORRECT DISPOSAL OF THE PRODUCT IN ACCORDANCE WITH EC DIRECTIVE 2012/19/EU WEEE:

At the end of its life, the product must not be disposed as household waste. It can be taken to special recycling centers provided by local government, or on purchase of a new device of the same type and used for the same functions return it to the dealer. Disposal of the product separately avoids possible negative consequences for the environment and human health resulting from inappropriate disposal and allows recycling of the materials in order to obtain significant savings in energy and resources. The symbol on the label indicates separate collection of electrical and electronic equipment.

**Warning:** Incorrect disposal of electrical and electronic equipment could result in sanctions.

### Notice

The information in this manual is subject to change without notice.  
The images are included as examples and may vary slightly from the actual device.

© Copyright Spencer Italia s.r.l.

All rights reserved. No part of the document may be photocopied, reproduced or translated into another language without prior written consent from Spencer Italia s.r.l.