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
1. MODELS

The following basic models may be subject to implementation or change without notice.

100 SERIES	200 SERIES	300 SERIES	SLIM SERIES	STOCK SERIES	ELLA SERIES	PRIME SERIES
SPENCER 100	SPENCER 200	SPENCER 310	SPENCER 205	SPENCER 270	SPENCER 506	MSF
SPENCER 110	SPENCER 250			SPENCER 280	SPENCER 516	PROMYGES
SPENCER 122				SPENCER 281	ELLA SELF	TGS
SPENCER 128						

The above models are available with different coloured canvases and/or beds and can be supplied with a variable number of straps.

2. INTENDED USE

100/200/300/Slim/Stock series field stretchers:

Field stretchers are devices for lifting and carrying people which are also suitable for parking in field rescue conditions. If they are not equipped with supports, they can be parked using the Porta Stock stand.

Ellä series field stretchers:

Ellä series field stretchers are devices for lifting and carrying people for use in field rescue situations. The foldable, lockable legs allow for transferring either by means of wheels on the legs or by lifting when the legs are locked in the closed position. Parking is permitted with the stretcher alone or with a trolley with closed legs.

Prime series field stretchers:

PRIME series stretchers are devices for lifting and carrying people for use in field rescue situations. Parking is permitted preferably when using the parking stops installed in designated locations within the rescue area.

It is not foreseen that patients be able to intervene on devices.

TARGET PATIENTS

There are no particular indications related to the patient group.

The product configuration is able to accommodate any subject as long as he/she is within the maximum capacity and within the limits of the size of the device.

PATIENT SELECTION CRITERIA

The patients foreseen are those with injuries that prevent them from walking in a given rescue situation, or who are in a state of unconsciousness.

No particular contraindications or side effects are known with relation to use of the device, as long as it is used in accordance with the user manual.

CONTRAINDICATIONS AND UNWANTED SIDE EFFECTS

No particular contraindications or side effects are known with relation to use of the device, as long as it is used in accordance with the user manual.

USERS AND INSTALLERS

The intended users are field rescue workers, typically nurses, stretcher bearers and doctors present in the field.

These devices are not intended for lay people.

The device is intended for professional use only. All operators must be trained to transport patients safely and efficiently. Do not allow untrained persons to help while using the product, as they may cause injury to themselves or others.

Operators using the device should be physically able to use the device and have good muscle coordination, as well as strong back, arms, and legs for lifting and holding and be able to grasp the device firmly with both hands.

Operators must be able to provide the necessary patient care.

Users must be able to safely lift and handle the weight of the stretcher and patient together and any other equipment used with the device.

For patient loading techniques for particularly heavy patients, for operations on steep terrain or in special and unusual circumstances, the presence of multiple operators is recommended (not just 2 as expected in standard conditions).

The ability of each operator must be assessed before the definition of rescuer roles in use of the stretcher.

User training

Note: despite all efforts, laboratory tests, trials, and instructions for use, standards do not always reproduce practice, so the results obtained under actual conditions of product use in the natural environment may sometimes differ significantly.

The best instructions are the continuous practice of use under the supervision of competent and trained personnel.

• Regardless of your level of experience with similar devices in the past, you should carefully read and understand the contents of this manual before installing, operating, or servicing this product. In case of any questions, please contact Spencer Italia S.r.l. for the necessary clarifications.

• The product must be used only by personnel trained in the use of this product and not on other similar products.

• The suitability of the users for use of this product can be attested by the training registration, in which trained persons, trainers, date and place are specified. **This documentation must be kept for at least 10 years after the end of the product's life and must be made available to the competent authorities and/or the Manufacturer when requested.** In the absence of such documentation, the relevant bodies will apply any foreseen sanctions.

• Do not allow untrained persons to help while using the product, as they may cause injury to themselves or others.

• The product must be put into use only by personnel trained in the use of this product and not on other similar products.

Note: Spencer Italia S.r.l. is always available for training courses. Document user training on appropriate forms.

Installer training

If applicable, device installation must be carried out by qualified personnel who have been trained and certified by Spencer Italia S.r.l. Formazione installatore

3. REFERENCE STANDARDS

As Distributor or End-User of the products manufactured and/or marketed by Spencer Italia S.r.l., users are strictly required to be familiar with the legal provisions in force in the country of destination of the goods, applicable to the devices to be supplied (including regulations relating to technical specifications and/or safety requirements) and, therefore, to understand the requirements necessary to ensure compliance of the products themselves with all legal requirements of the territory.

REFERENCE	DOCUMENT TITLE
EU Regulation 2017/745	EU Regulation on Medical Devices

4. INTRODUCTION

4.1 USING THE MANUAL

The purpose of this manual is to provide healthcare professionals with the information necessary for safe and appropriate use and maintenance of the device.

Note: The Manual is an integral part of the device and therefore it must be kept for the entire life of the device and must accompany it in any changes of use or ownership. If any instructions for use for products other than the one received are present, please contact the Manufacturer immediately before use.

Spencer products User Manuals can be downloaded from the site <http://support.spencer.it> or by contacting the Manufacturer. Exceptions are those items whose essentiality and reasonable and predictable use are such that it is not necessary to draw up instructions, in addition to the following warnings and indications on the label.

IT

EN

DE

FR

ES

PT

CS

Regardless of your level of experience with similar devices in the past, it is advisable to carefully read and understand the contents of this manual before installing, operating, or servicing this product.

4.2 DEVICE LABELLING AND TRACEABILITY

Each device is provided with a label, placed on the device itself and/or on the packaging, which contains the Manufacturer's identification data, product, CE marking, serial number (SN) or lot number (LOT). This must never be removed or covered.

In the event of damage or removal, request a duplicate from the Manufacturer, or else the warranty will be void as the device can no longer be traced.

Regulation 2017/745/EU requires manufacturers and distributors of medical devices to keep track of their location. If the device is in a location other than the address to which it was shipped or sold, or if it was donated, lost, stolen, exported or destroyed, permanently removed from use, or if the device was not delivered directly from Spencer Italia S.r.l., please register the device at <http://service.spencer.it>, or inform Customer Service (see § 4.4).

4.3 SYMBOLS

Symbol	Meaning	Symbol	Meaning	Symbol	Meaning												
	Device in compliance with EU Regulation 2017/745		See the user manual		Product code												
	Medical device		Serial Number		Lot number												
	Manufacturer		Danger – Indicates a hazardous situation that may result in a situation directly related to serious injury or death.														
	Date of manufacture		Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner (only for USA Market)														
	Unique Device Identifier	Production identification Alphanumeric code that identifies the production units of the device, composed of:															
		<table><tr><td>(01)0805771123</td><td>company prefix</td></tr><tr><td>000</td><td>progressive GS1</td></tr><tr><td>6</td><td>control number</td></tr><tr><td>(11)200626</td><td>date of production (YYMMDD)</td></tr><tr><td>(21) 1234567890</td><td>SN</td></tr><tr><td>(10) 1234567890</td><td>Lot number</td></tr></table>				(01)0805771123	company prefix	000	progressive GS1	6	control number	(11)200626	date of production (YYMMDD)	(21) 1234567890	SN	(10) 1234567890	Lot number
(01)0805771123	company prefix																
000	progressive GS1																
6	control number																
(11)200626	date of production (YYMMDD)																
(21) 1234567890	SN																
(10) 1234567890	Lot number																
(01)08057711230006(11)200626(21)1234567890																	

4.4 WARRANTY AND SERVICE

Spencer Italia S.r.l. guarantees that products are free from defects for a period of **one year from the date of purchase**.

For information regarding correct interpretation of the instructions for use, maintenance, installation or return, please contact Spencer Customer Service tel. +39 0521 541111, fax +39 0521 541222, e-mail service@spencer.it.

To facilitate service, always indicate the lot number (LOT) or serial number (SN) on the label attached to the package or device itself.

Warranty and service conditions are available at <http://support.spencer.it>.

Note: Record and keep with these instructions: lot (LOT) or serial number (SN), if present, place and date of purchase, date of first use, date of checks, user name and comments.

To ensure the traceability of the products and protect maintenance and service procedures on your devices, Spencer has made the SPENCER SERVICE portal (<http://service.spencer.it/>) available to you. From this site, you can view the data of the products in your possession or placed on the market, monitor and update schedules for periodic checks and view and manage special maintenance.

5. WARNINGS

Warnings, notes, and other important safety information are provided in this section and are clearly visible throughout the manual.

 At least every 6 months, it is important to check for updated instructions and any changes involving your product. This information is freely available on the website www.spencer.it on the specific product page.

Product features

Use of the product for any purpose other than that described in the User Manual is prohibited.

- Before each use, always check the conditions of the product, as specified in the User Manual. In the event of faults/damage that could compromise its functionality/safety, immediately remove it from service and contact the Manufacturer.
- If the product is found to be malfunctioning, immediately use a similar device to ensure continuity of ongoing operations. Non-compliant devices must be taken out of service.
- The product must not be tampered with or modified without the manufacturer's authorisation (modification, tweaking, additions, repair, use of non-approved accessories), as they may constitute imminent danger of injury to persons and material damage. Should these operations be performed, we decline any responsibility for incorrect operation or any damage caused by the product itself; moreover, the CE marking and the product warranty shall be null and void.
- Anyone who modifies or causes products manufactured by Spencer Italia S.r.l. to be modified or re-processed in such a way that they no longer serve the intended purpose or no longer provide the intended performance shall comply with the conditions valid for the first placing on the market.
- When using the devices, position and adjust them in such a way that they do not hinder operator works or the use of any other equipment.
- Be sure to take every precaution to avoid hazards from contact with blood or body secretions, if applicable.
- Avoid contact with sharp objects.
- Device installation must be carried out by qualified personnel trained and enabled by Spencer Italia S.r.l. The times and methods of such shall be agreed upon between the customer and our Sales Offices.
- Operating temperature: from -10°C to + 50°C.

Storage

- The product must not be exposed or come into contact with thermal sources of combustion or flammable agents, but must instead be stored in a dry, cool place, away from light and sun.
- Do not store the product under other more or less heavy materials that may damage the structure of the product.
- Store and transport the product with its original packaging, otherwise the warranty shall be invalidated.
- Storage temperature: -20°C to +60°C.

Regulatory requirements:

As Distributor or End-User of the products manufactured and/or marketed by Spencer Italia S.r.l., users are strictly required to be familiar with the legal provisions in force in the country of destination of the goods, applicable to the devices to be supplied (including regulations relating to technical specifications and/or safety requirements) and, therefore, to understand the requirements necessary to ensure compliance of the products themselves with all legal requirements of the territory.

- Promptly and in detail notify Spencer Italia S.r.l. (already in the quotation request phase) about possible fulfillments by the Manufacturer necessary for the compliance of products with specific legal requirements of the territory (including those deriving from regulations and/or regulatory provisions of another nature).

- Act with due care and diligence to help ensure compliance with the general safety requirements of the devices placed on the market, providing end-users with all the information necessary to carry out periodic revisions on the supplied devices, exactly as indicated in the User Manual.
- **Participate in safety checks on products** placed on the market, transmitting information regarding product risks to the Manufacturer as well as to the Competent Authorities for their respective actions.
- Without prejudice to the above, the Distributor or End-User shall assume wider liability related to non-compliance with non-fulfilment of the above-mentioned obligations, with consequent obligation to indemnify and/or hold Spencer Italia S.r.l. harmless from any possible injurious effect.
- With reference to EU Regulation 2017/745, please note that public or private operators who, when exercising their activity, detect an incident involving a medical product are required to notify the Ministry of Health, within the terms and in the manner established by one or more ministerial decrees, and notify the Manufacturer. Public or private health care professionals are required to notify the Manufacturer of any other incident that may allow the adoption of measures to ensure the protection and health of patients and users.

General warnings for medical devices

The user must carefully read the following in addition to the general warnings.

- Application of the device should not last longer than the time strictly necessary for pre-hospital treatment.
- Qualified personnel and at least two operators must be present during use of the device.
- Follow the internal procedures and protocols approved by your organisation.
- Disinfection operations must be carried out in accordance with the validated cycle parameters, as stated in the specific technical standards.

6. SPECIFIC WARNINGS

To use the product, you must also have read, understood and carefully follow all the instructions in the user manual.



- Follow approved Emergency Medical Service procedures for patient immobilization and transportation.
- Follow approved Emergency Medical Service procedures for patient positioning and transportation.
- Do not use if the device or parts of it are punctured, torn, frayed, or excessively worn.
- Before any type of handling, make sure that operators have a secure grip on the device.
- In case of exposed and/or injured skin, cover the surfaces in contact with the patient with a surgical canvas that respects bio-compatibility regulations to protect the patient's health.
- Do not drag the device on uneven surfaces.
- Do not lift with a crane or other mechanical lift.
- Do not use drying machines.
- The device is a transport device and **cannot be used as a stationing device**.
- Do not use with devices other than those expressly approved by the manufacturer.
- Practice with a device without a patient to make sure you are familiar with the manoeuvres.
- For patient loading techniques for particularly heavy patients, for operations on steep terrain or in special and unusual circumstances, the presence of multiple operators is recommended (not just 2 as expected in standard conditions).
- Make sure that the patient is adequately immobilized before loading him/her on the stretcher. A failure to immobilize can cause serious damage to the patient.
- Make sure that the sheet does not interfere with any stretcher handling or control mechanisms.
- Do not move the device if the weight is not properly distributed.
- Always use straps that are anchored to the stretcher frame to ensure patient safety.
- Use only the perimeter frame or the handles to move the stretcher and not other points not intended for this purpose.
- **Hold the device firmly if the patient is lying on it.**
- **Where present, parking brakes are aids for the operator but they do not in any way replace his or her supervision.**
- Pay the utmost attention to any obstacles (water, ice, debris, etc.) present on the route, as they may cause the operator to lose balance and compromise proper functioning of the device. If you cannot clear the way, choose an alternative route.
- Condensation, water, ice and accumulations of dust can affect correct operation of the device, making it unpredictable and causing a sudden change in the weight the operators have to support.
- When closing stretchers, do not place your fingers near the moving parts of the frame as there is a finger crushing hazard.
- Always comply with the maximum capacity indicated in the User Manual. Maximum load capacity means the total weight distributed according to human anatomy. When determining the total weight load on the product, the operator should consider the weight of the patient, equipment and accessories. In addition, the operator should assess whether the overall size of the patient reduces the functionality of the product.
- Make sure that operators are in good physical condition before lifting, as listed in the User Manual.
- **The maximum weight, which weighs on each operator, must comply with local health and safety requirements.**

7. RESIDUAL RISK

No residual risks, or rather risks that could arise despite compliance with all warnings in this user manual, have been identified.

8. TECHNICAL DATA AND COMPONENTS

Note: Spencer Italia S.r.l. reserves the right to make changes to specifications without notice.

PARTS		DESCRIPTION AND MATERIALS
 100 series	1	Lifting handles Made of rubber, these are the gripping points provided for lifting the device.
	2	Support canvas Made of PVC, this is the part used to support the patient.
	3	Joint points These are used to fold the stretcher. Models 110/122/128 also have hinge points underneath the PVC canvas to allow folding.
	4	Support feet Made of aluminium and attached to the frame, these keep the stretcher raised from the ground.
PARTS		DESCRIPTION AND MATERIALS
 200 series	1	Support feet Made of steel and attached to the frame, these keep the stretcher raised from the ground. Present on model 200 only; 250 is not equipped with them.
	2	Support canvas Made of PVC, this is the part used to support the patient.
	3	Joint points These are used to fold the stretcher.
	4	Supports on wheels Made of steel and attached to the frame, these are complete with wheels to facilitate small movements on the ground. These keep the stretcher raised from the ground together with the support feet. Present on model 200 only; 250 is not equipped with them.

IT

EN

DE


FR


ES

PT


CS

IT

PARTS		DESCRIPTION AND MATERIALS
300 series		Support feet 1 Made of steel and attached to the frame, these keep the stretcher raised from the ground.
		Support canvas 2 Made of PVC, this is the part used to support the patient.
		Joint points 3 These are used to fold the stretcher.
		Backrest adjustment rack 4 Used to adjust the backrest to several tilt positions.
		Supports on wheels 5 Made of steel and attached to the frame, these are complete with wheels to facilitate small movements on the ground. These keep the stretcher raised from the ground together with the support feet.

PARTS		DESCRIPTION AND MATERIALS
Slim series		Supports on wheels 1 Made of steel and attached to the frame, these are complete with wheels to facilitate small movements on the ground. These keep the stretcher raised from the ground together with the support feet. These are foldable on Slim series stretchers.
		Support canvas 2 Made of PVC, this is the part used to support the patient.
		Joint points 3 These are used to fold the stretcher.
		Support feet 4 Made of steel and attached to the frame, these keep the stretcher raised from the ground. These are foldable on this type of stretcher.

EN

PARTS		DESCRIPTION AND MATERIALS
Stock series		Support canvas 1 Made of PVC, this is the part used to support the patient.
		Supports 2 Made of steel and attached to the frame through welding, these keep the stretcher raised from the ground. They also act as a guide when you want to stack several stretchers.
		Backrest adjustment rack 3 Used to adjust the backrest to several tilt positions. Not present on model 270.


DE

FR


ES

PT

CS


PARTS		DESCRIPTION AND MATERIALS
516		
	1	Reclining sidebars Made of aluminium, steel and nylon and installed on the 506 stretcher, these allow the patient to be held laterally.
	2	Trendelenburg/Fowler (or Treendelenburg) lever Made of steel, positions the patient bed in the respective position.
	3	Feet side patient bed Made of polyethylene, equipped with footrests.
	4	Telescopic stretcher handle Made of steel and equipped with a rubber knob, these are used to grip the stretcher for lifting.
	5	Stretcher release control Made of steel, if pulled this releases the stretcher from the trolley.
	6	Rear leg unlocking lever Made of plastic-coated steel, unlocks rear leg movement.
	7	Leg Lock Control Made of steel, when operated with the legs in the closed position, it keeps the legs in this position.
	11	Support wheels Made of TPU and Nylon, these allow the trolley to be slid with the legs closed onto the board on which it is loaded.
	12	Front wheels Made of polyurethane, these allow the trolley to slide over the ground.
	13	Trolley knobs Made of steel and fitted with a rubber knob, these are used to grip the trolley for lifting.
	14	Leg locking hooks Made of steel, they are part of the locking system described in point 7.
	15	Loading trolley Made of steel and adjustable, it rests on the loading surface and supports the stretcher when the front leg is closed.
	16	Stretcher wheels Made of TPU and Nylon, these allow the stretcher to slide on the ground and are the part that rests on the trolley when the trolley and stretcher are coupled.
	17	Stretcher crossbars Made of steel, it allows coupling with the trolley.

8	Front leg unlocking lever Made of plastic-coated steel, it unlocks front leg movement.	18	Telescopic handle button Made of nylon, it unlocks telescopic handle sliding when pressed.
9	Rear leg Made of steel, this is the mobile part of the trolley that allows the equipment to be moved and loaded onto any support surfaces.	19	Backrest unlocking knob Made of nylon, it unlocks movement of the backrest when pulled, to return it to a lower position than it is in.
10	Bracket with wheel Made of chrome-plated iron, rotating and with a brake, it allows the stretcher to be moved and driven.	Note	If only the detachable stretcher (506) is used, please only consider the applicable callouts.

Ella Self is equipped with the same components as 516 except for those used for attaching the two parts.
There is in fact no detachable stretcher on Ella Self, but the patient surface is integrated in the trolley

PRIME SERIES

MSF

PROMYGES

TGS

Wall-mounted anchorage		Front wall-mounted lock	
1	Made of steel, this must be attached to the wall for stretcher anchorage.	7	Made of aluminium, allows the stretcher to be anchored by means of the front wheel axis.
Stretcher stop		Rear stop arch	
2	Made of steel, it must contain a stretcher handle in order to keep the stretcher in position.	8	Made of steel, it allows the stretcher to be anchored by means of the rear stop (accessory).
Wheel with brake		Patient bed	
3	Made of rubber and nylon, allows the stretcher to slide on the ground.	9	Made of polyethylene, it is the patient's support surface.
Stabiliser belt		Telescopic handles	
4	Helps further stabilise the stretcher by anchoring it to any elements attached to the ground.	10	Made of steel, these are unlocked by pressing the corresponding nylon button to facilitate operator lifting.
Front stop arch		Wall-mounted backrest	
5	Made of steel, it allows the stretcher to be anchored by means of the front stop (accessory).	11	Made of PVC, this provides back support if the patient sits on the stretcher rather than lying on it.
Arch for rear stop			
6	Made of steel, it allows the stretcher to be anchored by means of the rear stop (accessory).		

TECHNICAL DATA

Characteristics	Spencer 100	Spencer 110	Spencer 122	Spencer 128
Length when open (mm)	2010 ± 10 mm	2010 ± 10 mm	2010 ± 10 mm	2095 ± 5 mm
Length when folded (mm)	1000 ± 10 mm	2010 ± 10 mm	1000 ± 10 mm	Nd
Width when open (mm)	480 ± 5 mm	550 ± 5 mm	550 ± 5 mm	540 ± 5 mm
Width when folded (mm)	–	110	400	Nd
Height when open (mm)	140 ± 5 mm	140	140	Nd
Height when folded (mm)	170 ± 5 mm	140 ± 5 mm	170 ± 5 mm	nd
Height from the ground	nd	nd	nd	70 ± 5 mm
Weight (kg)	5,5 ± 0,3 kg	6 ± 0,3 kg	6,5 ± 0,3 kg	6 ± 0,3 kg
Capacity (kg)	170	170	170	159
Gripping points (no.)	4	4	4	4
Support feet (no.)	4	4	4	4
Dimensions when folded	nd	nd	nd	580x210x120 ± 10mm

Characteristics	Spencer 200	Spencer 250	Spencer 310	Spencer 205
Length when open (mm)	1885 ± 5 mm	1890 ± 10 mm	1870 ± 10 mm	1890 ± 10 mm
Length when folded (mm)	925 ± 5 mm	950 ± 10 mm	945 ± 10 mm	950 ± 10 mm

IT

EN

DE

FR

ES

PT

CS

IT

Width (mm)	490 ± 5 mm	480 ± 5 mm	500 ± 5 mm	480 ± 5 mm
Thickness (mm)	–	30 mm	–	–
Thickness when folded (mm)	165 ± 5 mm	60 ± 5mm	160 ± 5 mm	80 ± 5mm
Height from the ground	–	–	160 ± 5 mm	180 ± 5 mm
Weight (kg)	5 ± 0,2 kg	5,5 ± 0,3 kg	6,5 ± 0,3 kg	6,2 ± 0,2 kg
Maximum load capacity (kg)	170	170	170	170
Gripping points (no.)	–	–	–	8
Wheels	2 each 100mm	–	2 each 100mm	2
Support feet	2	–	2	2

Characteristics	Spencer 270	Spencer 280	Spencer 281
Length when open (mm)	1885 ± 10 mm	1885 ± 10 mm	1885 ± 10 mm
Width (mm)	480 ± 5 mm	480 ± 5 mm	480 ± 5 mm
Height (mm)	175 ± 10 mm	175 ± 10 mm	175 ± 10 mm
Weight (kg)	5 ± 0,5 kg	6 ± 0,5 kg	6 ± 0,5 kg
Capacity (kg)	170	170	170
Gripping points (no.)	8	8	6
Support feet (no.)	4	4	4
Approximate overall dimensions 10 stacked stretchers	610 ± 30 mm	nd	nd

EN

Characteristics	MSF	PROMYGES	TGS
Stretcher length	1895	1820	1915
Stretcher width	575	575	575
Possible anchorage points	3	2 with Front and rear stop accessory	2 with Rear stop accessory

DE

Characteristics	Stretcher 506	Trolley	ELLA SELF
Length w/ handles closed	1930 ± 10 mm	1980 ±10 mm	1980 ±10 mm
Length w/ handles open	2340 ± 10 mm	–	–
Width	570 ± 10 mm	590 ± 20 mm	570 ± 10 mm
Wheels	ø 100 mm	ø200 mm	ø200 mm
Dimensions sidebars	680x200 mm	–	680x200 mm
Materials	Al, PE, steel	Steel/Nylon/Al	Steel/Nylon/Al/PE
Weight	22 ± 1 kg	34 ± 1 kg	40 ± 1 kg
Capacity	170 kg	170 kg	170 kg
Loading trolley underwheel height	–	–	Min 44 ± 1cm • Max 55 ± 1 cm

FR

9. COMMISSIONING

For first use, check that:

- Packaging is intact and has protected the device during transportation
- Check that all parts included in the packing list are present
- General functionality of the device
- If present, check that the patient bed adjustment systems lock in place properly
- If present, check that all controls are working properly
- If present, check that all unlocking mechanisms are working properly

ES

The procedures for carrying out the above checks are described in par.11.

Do not modify the stretcher structural, leverage or traction parts in any way as this could cause damage to the patient and/or the rescuers.

 Failure to take the above measures will preclude safe use of the device, resulting in risk of damage to the patient, operators and the device itself.

For subsequent use, perform the operations specified in paragraph 10.

If the above conditions are met, the device may be considered ready for use; otherwise, you must immediately remove the device from service and contact the Manufacturer.



Do not alter or modify the device arbitrarily, as doing so could result in unpredictable operation and damage to the patient or rescuers and will void the warranty and release the Manufacturer from all liability.

PT


CS

10. PROPER USE


Primary medical evaluations must be carried out before the patient is moved, lifted or transported. Once the diagnosis has been made, it is preferable to advise the patient to actively contribute to the transition from bed to stretcher/chair, while at the same time informing them of the risks they may be running. Before loading the patient, bring the device as close as possible to the patient.

Opening and closing widthwise stretchers 110, 122 and 128	Opening and closing lengthwise (100, 122, 128, 200, 205 and 310)
 Opening Open the stretcher by spreading the handles as far apart as possible. Press the joint in the direction of the handles until it locks. Perform the same operation on the opposite side of the stretcher. Closing Pull the joint in the opposite direction to the handles until it is released. Perform the same operation on the opposite side of the stretcher. Bring the handles closer together Caution: Danger of pinching, only press on the rigid parts, keeping your fingers away from the joints.	 To fold the stretcher, simply bring the ends of the frame together as shown in the image. On model 128, the two halves of the stretcher can be folded into two more parts, further reducing the bulk. To unfold the stretcher, perform the reverse unfolding operation.


IT

Opening and closing stretchers 250
 The safety sleeves must be unlocked in order to fold stretchers 250. Rotate the sleeve to release the locking pin from its seat. Once unlocked, slide the sleeve along the frame until the joints are fully exposed. Perform the same operation on the other side of the stretcher. The stretcher can now be folded as illustrated above for the other models. Perform the above steps in reverse order to open the stretcher. When locking the sleeve, make sure that the locking pin engages properly when you turn the sleeve. Do not lift the patient without making sure that the safety mechanism is properly engaged.

EN

Folding feet model 205
 Model 205 is equipped with folding feet. To fold, press the unlocking buttons located at the base of the brackets. After unlocking the brackets, rotate them inwards 90° until they are in the closed position. The stretcher can now be folded as illustrated above for the other models. Perform the above procedures in reverse order to open the stretcher.


DE

Using the adjustable backrest (281, 310)
 To tilt the backrest, lift it to approximately the desired position. Then find the nearest point of engagement on the rack, which will hold the backrest in position. Before releasing the grip on the backrest, ensure that both pins of the backrest support frame are correctly engaged in the rack. To return the backrest to the horizontal position, lift the backrest slightly to pull the pins out of the rack and, while keeping the backrest support frame slightly raised, guide the backrest back to the horizontal position.

FR

Stacking Stock stretchers (270, 281)
Stock series stretchers can be stacked in order to optimise storage space. To stack the stretchers, simply set them in the same direction and stack them one on top of the other. Do not stack more than 10 stretchers, as stacking more than that will reduce stability due to the height and may cause the stack to tip over, posing a risk to those nearby.

ES

Using telescopic handles (Ella and Prime series)
 Telescopic handles make lifting easier and less tiring. To open them, press the buttons near the handles, then pull the handles outwards. To return them to their original position, press the buttons and push the handles inwards until they lock. Telescopic handles without an unlocking button are released by pressing down on the knob. After applying this pressure, pull the handles outwards. To return them to their original position, push them inwards until they lock.


PT

CS


Ella series controls and operation

The 516 and Ella Self are equipped with an adjustable loading trolley that enables the loading height of the stretcher to be adjusted to the surface on which it is to be loaded. The possible range of adjustment is as follows:

Height under wheel (cm)		VLF height (loading board)	
Min	Max	Min	Max
49	62	55	68

B

These heights can be obtained by adjusting the fastening screws on the loading trolley and positioning the trolley in a suitable position for the surface on which the stretcher is to be loaded.

Image A shows how to adjust the loading trolley in order to obtain the minimum height position above the ground.


Image B shows how to adjust the loading carriage in order to obtain the maximum height position above the ground.

Caution: It is not possible to use the leg locking system in the closed position when the loading trolley is adjusted to the minimum height positions.

IT

EN

DE

Adjustments to the trolley have a direct effect on the positioning of the rear fastening system installed on the board. Adjustment should be made taking into account that, in simulated loading conditions and with the loading trolley on the board, the wheels of the front legs should be 5/6 cm off the ground.

The rear stop can be installed after adjusting the height of the loading trolley (if not carried out by the manufacturer according to the specifications provided by the customer).


To correctly install the rear stop, bring the stretcher into contact with the front stop, pull the lever out of the rear stop as far as it will go and hold it in a vertical position.

Position the brake assembly so that it is perfectly centred on the stretcher and the upper end of the lever is approximately 1cm from the stretcher arch.

Then mark the position of the holes and drill the holes using the clamp assembly as a template.

After screwing the assembly to the table, check that the lever is pressed against the stretcher stop and that there is no more than 5mm of a gap.

Folding sidebars

The sidebars are designed to be used with two hands to prevent accidental movement.

To unlock the sidebars from the closed position, use both hands to pull outwards on the drives at the bottom as shown in the figure. The sidebar will automatically lower, turning approximately 180°.

To close the sidebars back up, bring them back to the vertical position and press simultaneously with both hands on the side profiles until the locking system is engaged. Check that the sidebars are correctly locked by making a few repeated pulling movements.

FR

ES

Lowering the stretcher


Make sure that the sidebar is in the closed position; otherwise the stretcher is likely to become damaged. Lower the stretcher without the patient as follows:

- Grasp the trolley handles on 516 from the "foot" (rear) side.
- Lift the foot end (rear) until the front loading trolley wheels touch the ground.
- Engage both trolley closing control levers.
- Lower the rear part down to the ground, supporting the weight.

Before loading the patient, while the stretcher is on the ground, make sure that the ground underneath is flat, as non-horizontal surfaces can affect the static balance of the stretcher. Spencer 516 is easier to manoeuvre in the open position; therefore, whenever possible, the stretcher should be brought to the desired point before lowering it.

PT

Lifting the stretcher with the patient


Secure the patient on the stretcher with the straps provided.

The operators should position themselves at each end of the stretcher (one at the feet, the other at the head) and check that the leg locks are disengaged.

Using the appropriate lifting technique to avoid fatigue, the operators should grasp the stretcher frame at each end and lift the stretcher until the legs have opened and check to make sure that the legs are secure and stable. If the entire stretcher is to be lifted, the legs can be locked in the closed position by pulling the appropriate locking lever (image opposite).

CS

Loading the stretcher on the board

Make sure that the folding sidebars are in the closed position; otherwise the stretcher can become damaged. When loading the stretcher:

- Push the stretcher in the direction of the loading board.
- Set the loading wheels on the board and move forward until the front legs have touched the edge of the surface.
- The operator must make sure that both wheels of the loading trolley are on the board (safety position) and simultaneously assess the distance of the front wheels from the ground (at least 5/6 cm).
- Once the front legs touch the board, the operator, without pushing the stretcher, must engage the right lever (1) (close front legs).
- After the lever has been engaged, push the stretcher, allowing the front legs to close and fold. Avoid extreme bumping of the front legs against the edge of the board. The legs and their mechanisms may become damaged.
- Push the stretcher until the rear legs touch the board; then raise the stretcher and pull up the left lever (2) (rear leg lock); then complete stretcher loading.
- Make sure that the wheels of the loading trolley fit into the front stop (if supplied).
- Secure the stretcher with the rear stop (if supplied).

We recommend the use of the original Spencer front and rear stops for anchoring the stretcher as they guarantee optimal stability in all directions. If you use fastening systems other than those listed above, make sure that they work properly.

Unloading the stretcher

Unlock the stretcher from its fastening systems.

Grasp the stretcher by the rear handles on the frame and pull it towards you without operating any controls, taking care to keep the stretcher in a horizontal position until the rear leg is fully open, where it will lock automatically. After checking that the rear leg has been opened correctly, proceed with unloading until the front leg has been fully opened, where the same check must be carried out.



Make sure that the ground underneath is safe and stable before extracting the front loading trolley wheels.

Brakes

The brake can be engaged by pulling the lever on the wheel. The brake system not only locks wheel rotation also stops complete wheel unit rotation with respect to its vertical axis, creating adequate immobilization of the stretcher.

Press on the same from the opposite side to disengage the brake.

Securing the stretcher to the trolley


Engage the parking brakes on the lower trolley.

Place the stretcher on top of the trolley and lay it on top of the trolley so that the sliding guides engage with each other.

Slowly push Stretcher 506 towards the front of the trolley by operating the hook lock on the trolley, making them a stable and secure system.


Removing the stretcher from the trolley

Pull the knob on the rear of the trolley with your right hand and grasp the left hand knob of the stretcher with your left hand and simultaneously pull it towards you.

The two operators, one at the front and one at the rear, must grasp the knobs of stretcher 506 and lift it up.

Stretcher board adjustments

Always inform the patient before making adjustments.

Changing the backrest inclination from a horizontal to vertical position:

Lift the backrest until the first position is reached, where it will automatically lock. Continue in this way to reach subsequent positions (seven positions are possible in addition to horizontal), always verifying correct activation of the locking system.

Changing the backrest inclination from a vertical to horizontal position:

Support the backrest with one hand (to prevent unexpected movements) and take away the weight from the backrest by moving it upward.

At the same time, pull the plastic knob on the backrest piston upward until the safety mechanism unlocks. Then, supporting it with the other hand, lower the backrest toward the horizontal position until you reach the desired height and then release the knob to the rest position. To reach subsequent positions, lift and support the backrest movement (seven positions are possible in addition to the horizontal position). Always lift the weight of the backrest before pulling the knob. Incorrect procedure could permanently damage the backrest piston.

IT

EN

DE


FR

ES

PT

CS

Trendelenburg/Fowler position

The stretcher boards can be moved in Trendelenburg and Fowler positions, which raises the patient's lower limbs.

The Trendelenburg/Fowler system can be put into three different positions:

Position A of the system allows the two boards to be tilted so as to obtain an angle with the vertex pointing upwards.

Position B of the system allows the boards to be raised in relation to the horizontal axis of the stretcher.

The middle position between the two described above keeps the board horizontal.

Restraint straps

To apply the two-piece restraint straps, identify a portion of the frame at chest height and one at the patient's leg level. The restraint straps should be applied in a position that allows for adequate immobilization of the patient.

Secure to the frame by means of a slipknot.

The restraint straps supplied with the stretcher must be fastened immediately after positioning the patient.

To do this, insert the male buckle into the female buckle.

11. CLEANING AND MAINTENANCE

Spencer Italia S.r.l. declines all responsibility for any direct or indirect damage which is the consequence of improper use of the product and spare parts and/or in any case of any repair carried out by a person other than the Manufacturer, who uses internal and external technicians authorised to do so; moreover, doing so will invalidate the warranty.

- The operator must wear suitable personal protective equipment, such as gloves, goggles, etc. during all checking, maintenance and cleaning operations.
- Establish a maintenance schedule, periodic inspections and extend the average life span, if foreseen by the Manufacturer in the User Manual, identifying a reference person who meets the basic requirements set forth in the User Manual.
- **The frequency of inspections is determined by factors such as legal requirements, type of use, frequency of use, and environmental conditions during use and storage.**
- Repairs of products manufactured by Spencer Italia S.r.l. must be carried out by the Manufacturer, who shall make use of specialised internal or external technicians who, using original spare parts, shall provide quality repair service in strict compliance with the technical specifications indicated by the Manufacturer. Spencer Italia S.r.l. declines any responsibility for any direct or indirect damage which is a consequence of improper use of spare parts and/or any repair work carried out by unauthorised parties.
- Use only original or Spencer Italia S.r.l. approved components/spare parts and/or accessories in order to carry out all operations without causing alterations or modifications to the product.
- All maintenance and overhaul activities must be recorded and documented with the relevant technical operation reports. This documentation must be kept for at least 10 years after the end of the product's life and must be made available to the competent authorities and/or the Manufacturer when requested.
- Cleaning, provided for reusable products, must be carried out in accordance with the Manufacturer's instructions in the User Manual in order to avoid the risk of cross-infection due to the presence of body fluids and/or residues.
- The product and all its components must be washed and left to dry completely before storage.
- If the product requires lubrication, this must be done after cleaning and complete drying.

11.1 CLEANING

Failure to carry out the correct cleaning operations could increase the risk of cross-infection due to presence of body fluids and/or residues.

The operator must wear suitable personal protective equipment, such as gloves, goggles, etc. during all checking and cleaning operations.

Metal parts exposed to external agents undergo surface treatments and/or coating in order to obtain better resistance. Clean the exposed parts with water and delicate soap.

Never use solvents or stain removers.

Do not use cleaning agents containing chlorine / sodium hypochlorite as corrosion of components may occur.

Rinse thoroughly with lukewarm water, making sure you have removed all traces of soap, which may deteriorate or compromise conditions and durability. **Avoid using high pressure water**, as it penetrates the joints and removes lubricant, creating the risk of corrosion on components. Let dry completely before storing. Drying after washing or after use in a wet environment must be natural and not forced. Do not use flames or other direct heat sources.

If **disinfecting**, use products that do not have a solvent or corrosive action on materials constituting the device. Be sure to take all precautions to ensure that there is no risk of cross-infection or contamination of patients and operators.

11.2 ROUTINE MAINTENANCE

Establish a maintenance program and periodic checks, identifying a designated reference person. The person entrusted with device maintenance must ensure the basic requirements envisaged by the manufacturer in the following paragraphs.

All routine and special maintenance activities and all general overhauls must be recorded and documented with the relevant technical intervention reports. This documentation must be kept for at least 10 years after the end of the device's life and must be made available to the competent authorities and/or the Manufacturer when requested.

Routine device maintenance must be carried out by operators with specific qualifications, training and instruction in the use and maintenance of the device.

The operator must wear suitable personal protective equipment, such as gloves, goggles, etc. during all checking, maintenance and cleaning operations.

Checks to be carried out before and after each use and in any case every 3 months are as follows:

- General functionality of the device
- Cleanliness of the device (remember that the failure of cleaning may cause the risk of cross infections)
- Correct tightening of all nuts, bolts and screws
- Absence of cuts, holes, tears on the structure, including the straps
- None of the tubes or metal sheets present bends or cracks
- Welded areas are intact, without any cracks or breaks
- Moving parts, wheels, levers, handles are intact and work properly
- Lubrication of moving parts
- Conditions of use of wheels and brake system
- The wheels are correctly fixed, are stable and turn properly
- The wheels are free of debris
- The device opens and locks properly

- The device opens and closes properly
- Functioning of springs

The frequency of inspections is determined by factors such as legal requirements, type of use, frequency of use, and environmental conditions during use and storage.

Please note that cleaning must be carried out as described in this manual and that functionality must be verified before and after each use. Spencer Italia S.r.l. declines any responsibility for improper functioning or damages caused to the patient or user by the use of devices not subject to routine maintenance and will void the warranty and compliance to Medical Device Directive 2017/745/EU.

Use only original or Spencer Italia S.r.l. approved components/replacement parts and/or accessories to carry out any operation without causing alterations or modifications to the device. Otherwise, we decline all responsibility regarding incorrect operation or any damage caused by the device to the patient or the operator, invalidating the warranty and invalidating compliance with Regulation 2017/745/EU Medical Devices.

11.3 PERIODIC OVERHAUL

Periodic overhaul is only foreseen for the Ella series, which must be overhauled every year by the Manufacturer, who uses internal and external technicians specialised and authorized by the Manufacturer itself.

In the absence of an above-mentioned overhaul, the device must be **DECOMMISSIONED**, as compliance with Regulation 2017/745/EU will cease and, despite the CE marking, the device will no longer meet the safety requirements guaranteed by the Manufacturer at the time of supply.

Spencer Italia srl declines any responsibility for incorrect operation or damage caused by use of devices that have not been regularly serviced.

Only overhaul activities carried out by specialised technicians authorised by the Manufacturer are considered valid by Spencer Italia S.r.l.

11.4 SPECIAL MAINTENANCE

Special maintenance can only be carried out by the Manufacturer, who uses internal and external technicians specialised and authorised by the Manufacturer itself.

Only maintenance activities carried out by specialised technicians authorised by the Manufacturer are considered valid by Spencer Italia S.r.l..

The end-user can replace only the spare parts indicated in § 15.

11.5 LIFE SPAN

The device, if used as described in the following instructions, has a life span of 5 years from the date of purchase.

This life span can be extended following annual overhauls, for the Ella series only.

The mattress and straps shall be replaced every two years.

Overhauls must be carried out by the Manufacturer, who uses internal and external technicians specialised and authorised by the Manufacturer. **In the absence of such annual overhauls, the device must be DISPOSED OF ACCORDING TO THE INFORMATION IN PARAGRAPH 16 AND NOTIFICATION MUST BE GIVEN TO THE MANUFACTURER.**

The life span may be extended at the sole discretion of the Manufacturer or authorised centre if the safety requirements of the device are still met.

Spencer Italia S.r.l. will accept no responsibility for incorrect operation or damage caused by the use of devices that have not been overhauled by the Manufacturer or authorized centre, or that have exceeded the maximum allowable life span.

12. TROUBLESHOOTING TABLE

PROBLEM	CAUSE	REMEDY
Canvas not sufficiently taut	Torn or broken canvas or bent frame	Immediately remove the device from service and contact customer service.
Locks defective	Damaged hinges	Immediately remove the device from service and contact customer service
Stretcher unstable	Damage to the frame	Immediately remove the device from service and contact customer service.
The Ella series leg unlocking mechanisms do not work or are difficult to engage.	The moving mechanisms are damaged.	Immediately remove the device from service and contact customer service.
	Means of connection between components have been lost.	Immediately remove the device from service and contact customer service.
Damage to the structure	Improper use	Immediately remove the device from service and contact customer service.
The front legs on the Ella series do not lock.	The moving mechanisms are damaged.	Immediately remove the device from service and contact customer service.
	The height of the board is not suited to the device, the safety height is not being respected.	Adjust the loading board and/or device to meet the requirements described in this manual. If the problem is not solved by adjustments, take the device out of service immediately and contact service.
The attachment of the stretcher to the trolley is not secure.	Wear or damage of the components that constitute the locking mechanisms.	Immediately remove the device from service and contact customer service.
The telescopic handles are not working.	Failure to engage the unlocking system or injury to the release mechanism	Make sure that you have activated the unlocking control (button or press down on the handle). If the problem persists, immediately remove the device from service and contact customer service.

13. ACCESSORIES

ST42800	FRONT FIELD STRETCHER ANCHORAGE	ST42802	REAR FIELD STRETCHER 506 ANCHORAGE
ST42801	REAR FIELD STRETCHER ANCHORAGE		

14. SPARE PARTS

There are no spare parts available for these devices.

15. DISPOSAL

When devices and their accessories are no longer suitable for use, they can be disposed of as normal municipal solid waste if they have not been contaminated by special agents. Otherwise, follow the regulations in force regarding disposal.

Warning

The information contained in this document is subject to change without notice and is to be intended as a commitment by Spencer Italia S.r.l. subject to change.
The images are included as examples and may vary slightly from the actual device.

© Copyright Spencer Italia S.r.l.

All rights reserved. No part of the document may be photocopied, reproduced or translated into another language without prior written consent from Spencer Italia S.r.l.

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© Copyright Spenser Italia S.r.l.