

CINCO MAS

10g Self loading stretcher with 5 heights



Spencer Cinco Mas is a self-loading stretcher with 5 adjustable heights for patient transport. It is a compact and light stretcher, usable in rescue operations on streets and rescue vehicles. The legs fold away independently, activating the dedicated commands during the unloading procedure, they open automatically. 10G certified if used with the fastener *Sensor Lock*.

Specific features

Polyethylene board easy to clean
Adjustable backrest
Antishock position Trendelenburg
Extreme maneuverability thanks to the Twist system which, when operated, makes pivoting all 4 wheels
Brakes on rear wheels
Foldable side bars <i>Safe bar</i>
5 Adjustable heights system
Integrated SHARP SENSOR LOCK system on the leg unlocking handle, which allows to unlock the stretcher from the fastener without acting on the fastener
Adjustable heights unlocking system logically placed to improve usability

Needed optional accessories

CR90013B	Loading trolley extender for Cinco Mas
----------	--

Optional accessories

ST42706	10G fastener sensor lock (Mandatory for 10g versions)
EN90003C	END-T Table for stretchers 10g cert.
ST00491A	STX 90 Telescopic headrest for stretchers
ST00499B	STX 499 – 4 pcs adjustable thorax belt
IF01047C	Track 4-30 telescopic IV holder
IF01049B	Track 5 telescopic IV holder for Carrera
CB09025C	Tanker oxygen tank holder for stretcher side bar
CB09028C	Tanker support for Safe Bar

Standard Equipment

ST70000A	QMX 777 BLACK ANATOMIC MATTRESS,
ST00497B	DNA STRAP Belt w/integrated re-winding system
ST00498B	DNA STRAP Th.Belt w/integr.re-winding system

Technical data

Length	1970 mm
Width	570 mm
Weight	42 Kg
Loading capacity	250 kg
Wheels	ø200mm ± 5%
Height loading carriage wheel	65 cm
Frame material	Stainless steel
Side bars length	680 mm
Side bars Height	200 mm

Class I Med.Device according to Regulation 2017/745/UE

CINCO MAS 65 T RED BOARD 10G
CERTIFIED

CND Classification
V0899



Rev.1 (26/05/2020)

Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it