

D-L-Handle

Adult fibre optic laryngoscope handle



Adult handle to be used with fibre optic blades.

Specific features

Non-slip surface treatment
High-brightness white led light
Needs 2 batteries type C not included

Blades connection standard

UNI EN ISO 7376

Technical data

Materials	Chromed brass
Handle diameter	30 mm
Overall lenght	160 mm
Handle lenght	130 mm
Weight (without batteries)	230 g

Class I Med.Dev. according to Dir. 93/42/CEE

D-L-HANDLE FIBRE-OPTIC
ADULT BLADE
DG03721B

CND Classification



Z12021003

Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222
www.spencer.it