

D-L-Handle

Pediatric fibre optic laryngoscope handle



Pediatric laryngoscope handle to be used with fibre optic blades.

Specific features

- Non-slip surface treatment
- High-brightness white led light
- Needs 2 batteries type AA not included

Technical data

Materials	Alluminum
Handle diameter	19 mm
Overall length	160 mm
Handle length	130 mm
Weight (without batteries)	160 g

Blades connection standard

Class I Med.Dev. according to Dir. 93/42/CEE
UNI EN ISO 7376

D-L-HANDLE FIBRE-OPTIC
PAEDIATRIC

DG03726

CND Classification



Z12021003

www.spencer.it

© Copyright - This document is property of SPENCER ITALIA S.r.l and can't be reproduced, used and disclosed without the written approval of SPENCER ITALIA S.r.l.