

SELF-LOADING STRETCHERS COMPARISON TABLE



SPENCER SELF-LOADING STRETCHERS	1 CINCO MAS	2 CROSS	3 CROSSOVER	4 CROSS UP 8409	5 CROSSOVER CHAIR	6 CARRERA TEC	7 CARRERA ACTIV	8 CARRERA PRO	9 CARRERA XL T TEC	10 CARRERA XL T TEC MAX
N° OF LEVELS	5	3	3	3	3	2	2	2	3	3
WEIGHT (kg)	42	45	45	52	59	37	36	36	37	38
LOAD CAPACITY (kg)	250	170	250	250	250	170	170	170	170	250
LOADING HEIGHT ¹ (cm)	65-75	From 52 to 70	From 47 to 70	65	65	50-60-63-68-73-75	60-63-64-68	50-60	55-58-60-66-68-75	50-60-63-68-75
SHORTENING BOARD LENGTH		●	●				●			
4 SWIVELLING WHEELS	●		●		●					
DETACHABLE TROLLEY/FIELD STRETCHER				●	●					
HOUSING FOR SPINE BOARD								●		

Dimensions and weights subject to tolerances of ± 10 mm and ± 1 kg respectively.
(1) Height from the ground as measured from right under the loading wheels.