



## EMERGENCY BACKPACK

**PAX WASSERKUPPE L****MAIN FEATURES**

- Durable and comfortable thanks to padded shoulder straps and back
- Made of durable and easy-care materials
- Equipped with an extended front compartment and a zipped side compartment
- 2 handles allow it to be carried both sideways and lengthways
- Riveted feet allow the backpack to remain stable on different surfaces
- Reflective elements increase the visibility of the backpack
- Customisable with options for inner bag sets, divider panels, materials, various colors and sizes
- Customised inscription available on request.



MANUFACTURED BY X-CEN-TEK

The **PAX WASSERKUPPE L** is specifically designed for everyday use in rescue services. It allows for customizable configurations with options for internal pouches sets, different materials and colors. Constructed for longevity and comfort, it is reinforced with padding and riveted feet to ensure stability on various terrains.

Two handles enable both longitudinal and transverse transport, while two-way zippers on each pocket, along with adjustable, removable shoulder straps, provide maximum convenience. Reflective stripes on all sides enhance the backpack's visibility.

The main compartment, which can be fully opened and is equipped with Velcro® straps, can accommodate color-coded PAX inner pouches, divider panels or modules and features additional inner pockets. An extended front compartment can accommodate essential equipment such as a collar, while a zipped side document compartment provides easy access to important documents.

42 x 57 x 23 cm  
16.5 x 22.4 x 9 in.



DIMENSIONS

2,50 kg  
5.51 lb



WEIGHT

38,88 L  
10.27 gal








CAPACITY



EASY TO CLEAN

## TECHNICAL DATA

Code	200050301	200050307	200050308	200054501	200054507
Dimension	42 x 57 x 23 cm 16.5 x 22.4 x 9 in.	42 x 57 x 23 cm 16.5 x 22.4 x 9 in.	42 x 57 x 23 cm 16.5 x 22.4 x 9 in.	42 x 57 x 23 cm 16.5 x 22.4 x 9 in.	42 x 57 x 23 cm 16.5 x 22.4 x 9 in.
Capacity	38,88 L 10.27 gal	38,88 L 10.27 gal	38,88 L 10.27 gal	38,88 L 10.27 gal	38,88 L 10.27 gal
Weight without accessories	2,50 kg 5.51 lb	2,50 kg 5.51 lb	2,50 kg 5.51 lb	2,50 kg 5.51 lb	2,50 kg 5.51 lb
Material	PAX-Tec	PAX-Tec	PAX-Tec	PAX-Rip Tec	PAX-Rip Tec
Color					



## AVAILABLE ACCESSORIES



### PAX INNER POUCH S - WELDED

S-size inner bag with Velcro®, for customising PAX bags/backpacks.

102093201 ●

102093203 ●

102093204 ●

102093205 ●

102093207 ●

102093211 ●



### PAX INNER POUCH M - WELDED

M-size inner bag with Velcro®, for customising PAX bags/backpacks.

102103201 ●

102103203 ●

102103204 ●

102103205 ●

102103207 ●

102103211 ●



### PAX INNER POUCH L - WELDED

L-size inner pouches with Velcro®, for customising PAX bags/backpacks.

102113201 ●

102113204 ●

102113205 ●

102113207 ●

102113211 ●



### PAX INNER POUCH FT - WELDED

FT-size inner bag with Velcro®, for customising PAX bags/backpacks.

102123201 ●

102123204 ●

102123205 ●

102123207 ●

102123211 ●



### PAX DIVIDER WASSERKUPPE L

Divider with compartments and Velcro® to customise the Wasserkuppe backpack.

165820310



### PAX INNER POUCH SET CHILD EMERGENCY - WELDED

Set of 5 inner pouches with Velcro®, for customising PAX bags/backpacks: 4 colors and 2 sizes.

130383200

## AVAILABLE ACCESSORIES



### **PAX INNER POUCH SET 1 - WELDED**

Set of 4 inner pouches with Velcro®, for customising PAX bags/backpacks: 4 colors and 2 sizes.

130373200



### **PAX INNER POUCH SET 2 - WELDED**

Set of 5 inner pouches with Velcro®, for customising PAX bags/backpacks: 4 colors and 3 sizes.

130373200



### **PAX SET INNER POUCH SET 3 - WELDED**

Set of 9 inner pouches with Velcro®, for customising PAX bags/backpacks: 6 colors and 4 sizes.

133743200



### **PAX WASSERKUPPE CENTRAL PANEL INFUSION/INTUBATION**

Central panel for infusion/intubation equipments and transparent pockets, to customise the PAX Wasserkuppe backpack.

206260310



### **PAX WASSERKUPPE CENTER PANEL MODULE**

Central divider panel with Velcro® for internal PAX pouches and transparent pockets, to customise the PAX Wasserkuppe backpack.

202180310

202200210



### **PAX CENTER PANEL WK AMPOULE KIT - PCI**

Central divider panel with PAX clip-in system installed in the middle segment and Velcro® to customise the PAX Wasserkuppe backpack. Accessories such as syringes and cannulae or aerosols and ampoules can be stored in the transparent pockets.

162490310

Rev.0 (23/10/2024) - UNCERTIFIED COPY  
DISTRIBUTED BY SPENCER ITALIA S.R.L.

©Copyright - This document is property of SPENCER ITALIA S.r.l and can't be reproduced, used and disclosed without the written approval of SPENCER ITALIA SRL.