



MULTI-LEVEL STRETCHER

AGILE**MAIN FEATURES**

- Stainless steel stretcher with micro shot-peening for enhanced longevity
- 6 adjustable heights
- Easy to access handle to change stretcher's height
- Assisted lowering system
- Shortening frame for enhanced mobility in tight spaces
- Trendelenburg and Fowler positions
- 4 pivoting wheels with brake
- Twist system to unlock front wheels
- Available in 3 models with different loading heights.

Tested according to ☒ UNI EN 1865-3 ☒ EN 1789

Class I Medical Device compliant with MDR (EU) 2017/745

PRODUCT DESCRIPTION: AGILE is a self-loading stretcher with 6 adjustable heights and 260 kg (573 lbs) load capacity, which can be used in confined spaces thanks to its shortening frame and 4 swivel wheels that facilitate the stretcher's movement. Color-coded handles assist with loading the stretcher into an ambulance. Available in three models with different loading heights, the stretcher can be adjusted to both Trendelenburg and Fowler positions.

INTENDED USE: Ambulance stretchers like AGILE are designed primarily for the safe transport of sick or injured individuals in a lying position. They are not intended for the patient to interact with the device, nor for prolonged use or as a substitute for a hospital bed. When used with the appropriate fixing systems, they are certified to withstand 10G forces.





573 lb
260 kg
**LOAD
CAPACITY**

**FOLDABLE
FRAME**

**ASSISTED
LOWERING**

6
**ADJUSTABLE
HEIGHTS**

TECHNICAL DATA

Models	AGILE 60	AGILE 68	AGILE 75
Codes	AG00001-DO AG00001-ZO	AG00002-DO AG00002-ZO AG00002-D1	AG00003-DO AG00003-ZO
Length	197,5 cm / 77.7 in	197,5 cm / 77.7 in	197,5 cm / 77.7 in
Shortened length	170 cm / 67 in	170 cm / 67 in	170 cm / 67 in
Width	57,5 cm / 22.6 in	57,5 cm / 22.6 in	57,5 cm / 22.6 in
Weight (without accessories) (1)	56,5 kg / 124.6 lb	56,5 kg / 124.6 lb	56,5 kg / 124.6 lb
Side bar dimension	68 x 19 cm / 26.7 x 7.48 in	68 x 19 cm / 26.7 x 7.48 in	68 x 19 cm / 26.7 x 7.48 in
Wheels diameter	20 cm / 7.8 in	20 cm / 7.8 in	20 cm / 7.8 in
Load capacity (1)	260 kg / 573 lbs	260 kg / 573 lbs	260 kg / 573 lbs
Loading height (2)	60 cm / 23.6 in	68 cm / 26.7 in	75 cm / 29.5 in
Number of brakes	4	4	4
Twist system	✓	✓	✓
6 adjustable heights	✓	✓	✓
Shortening frame	✓	✓	✓
Trendelenburg	✓	✓	✓
Fowler	✓	✓	✓
DNA Restraint belts (thoracic, femur/pelvis, ankles)	✓	✓	✓
QMX 777 Mattress	Included in AG00001-DO	Included in AG00002-DO and AG00002-D1	Included in AG00003-DO
Stretcher advancement extension	–	Included in AG00002-D1	–
Materials	Stainless steel, Al, PE	Stainless steel, Al, PE	Stainless steel, Al, PE
Frame color			
Board color			

(1) The weight is subject to a tolerance of ± 1 kg, depending on the loading height of the model of the stretcher.

(2) Height measured from the ground to the loading wheels.

RECOMMENDED FASTENING SYSTEM



ST42720

S-MAX PRO

10G Fastening system
for stretchers



AVAILABLE ACCESSORIES



EN90003

END-T
10G Instruments table
for stretchers



ST00491

STX 90
Telescopic headrest
for stretchers



ST00480

FIXO KID
Stretcher restraint
system for paediatric
patients



IF01047

TRACK 4
Adjustable IV pole



CB09012

OXYPACK
Oxygen tank carrier
for stretchers - 2 L



ST00046

S-WINGS
Bariatric enlargement
supports for AGILE

Rev.0 (06/12/2024)
UNCERTIFIED COPY

©Copyright - This document is property of SPENCER ITALIA S.r.l and can't be reproduced, used and disclosed without the written approval of SPENCER ITALIA SRL.

Spencer Italia s.r.l. a socio unico soggetta a direzione e coordinamento da Protect Medical Holding GmbH
Via Provinciale, n° 12 - 43038 Sala Baganza (PR) Italy - +39 0521 541111 - info@spencer.it - www.spencer.it