



10G FASTENING SYSTEM FOR SELF-LOADING STRETCHERS

S-MAX PRO**MAIN FEATURES**

- S-MAX PRO engaging and release mechanisms are easy to locate and activate
- It is easy to install thanks to its one-piece structure
- Its flat silhouette shape avoids unnecessary bulk when the stretcher is not loaded on the ambulance
- The position of the release mechanism minimize the risk of accidental injuries when the stretcher is pulled out
- The fixation system acts simultaneously both on the front and the rear of the stretcher.



Tested according to EN 1789

Class I Medical Device compliant with MDR (EU) 2017/745



S-MAX PRO is a 10G fastener for self-loading stretchers easy to install either on ambulance floors or on top of stretcher support systems. Its advanced anchoring system acts both on the front and the rear of the stretcher, ensuring a stable and safe stretcher fixation.

TECHNICAL DATA

Code	ST42720
Dimension	179 x 19 x 8 cm / 70.4 x 7.4 x 3 in
Metal profile thickness	3 cm / 1.81 in
Weight	5,2 kg / 11.4 lb
Materials	Al, Nylon

**ANTI-RUST
SURFACE**

5,2 kg / 11.4 lb

**DEVICE
WEIGHT****QUICK TO
ENGAGE**Product technical sheet **S-MAX PRO** – Code **ST42720** – CDN Classification **V9099** – Registration N° **2227449**

Rev.2 (13/02/2024)

UNCERTIFIED COPY

©Copyright - This document is property of SPENCER ITALIA S.r.l and can't be reproduced, used and disclosed without the written approval of SPENCER ITALIA SRL.

Spencer Italia s.r.l. a socio unico soggetta a direzione e coordinamento da Protect Medical Holding GmbH
Via Provinciale, n° 12 - 43038 Sala Baganza (PR) Italy - +39 0521 541111 - info@spencer.it - www.spencer.it