



MULTILEVEL SELF-LOADING STRETCHER

CROSS**MAIN FEATURES**

- Intermediate level
- Shortenable frame for improved mobility in confined spaces
- The loading height is adjustable during installation, offering greater flexibility in the choice of destination vehicle
- Wheels with 20 cm (7.8 in) of diameter
- Trendelenburg position.



CE **MD** Compliant with ☒ EN ISO 1865-1 ☒ EN 1789 ☒ EN ISO 1865-3

Class I Medical Device compliant with MDR (EU) 2017/745

PRODUCT DESCRIPTION: The CROSS is a self-loading stretcher with an intermediate-level position feature, supporting loads of up to 170 kg (375 lb). Its shortening frame allows it to be maneuvered easily in confined spaces, while color-coded handles ensure straightforward loading into ambulances. The 20 cm (7.8 in) diameter wheels enable smooth movement on flat surfaces, enhancing mobility and efficiency during transport.







INTENDED USE: Ambulance stretchers like CROSS are designed primarily for the safe transport of sick or injured individuals in a lying position. They are not intended for the patient to interact with the device, nor for prolonged use or as a substitute for a hospital bed. When used with the appropriate fixing systems, they are certified to withstand 10G forces.

375 lb
170 kg
**LOAD
CAPACITY**

**SHORTENING
FRAME**

3
**ADJUSTABLE
LEVELS**

TECHNICAL DATA

Codes	CR72010-D0 CR72010-Z0	CR72020-D0 CR72020-D1 CR72020-Z0 CR72020-Z1	CR72030-D0 CR72030-Z0
Length	197 cm / 77.5 in	197 cm / 77.5 in	197 cm / 77.5 in
Width	57 cm / 22.4 in	57 cm / 22.4 in	57 cm / 22.4 in
Weight (without accessories)	48 kg / 105.8 lb	48 kg / 105.8 lb	48 kg / 105.8 lb
Wheels diameter	20 cm / 7.8 in	20 cm / 7.8 in	20 cm / 7.8 in
Load capacity	170 kg / 375 lb	170 kg / 375 lb	170 kg / 375 lb
Loading height (1)	From 58 to 67 cm From 22.8 to 26.3 in	From 58 to 67 cm From 22.8 to 26.3 in	From 58 to 67 cm From 22.8 to 26.3 in
Intermediate height	✓	✓	✓
Shortening frame	✓	✓	✓
Trendelenburg	✓	✓	✓
Fowler	–	–	✓
Restraint belts (thoracic, femur/pelvis, ankles)	✓	✓	✓
QMX 777 Mattress	Included in CR72010-D0	Included in CR72020-D0 and CR72020-D1	Included in CR72030-D0
Stretcher advancement extension	–	Included in CR72020-D1 and CR72020-Z1	–
Materials	Stainless steel, Al, PE	Stainless steel, Al, PE	Stainless steel, Al, PE
Frame color			
Board color			

Note: Dimensions and weights subject to tolerances of ± 10 mm (0.4 in) and ± 1 kg (2.2 lb) respectively.
(1) Height from the ground as measured from right under the loading wheels.

SISTEMA DI FISSAGGIO RICHIESTO



ST42720

S-MAX PRO
10G Fastening system
for stretchers



AVAILABLE ACCESSORIES



ST00497

DNA STRAPS
Belts with integrated
re-winding system



ST00498

DNA STRAPS
Chest belts with
integrated retractor
system



ST00480

FIXO KID
Stretcher restraint
system for paediatric
patients



ST00491

STX 90
Telescopic headrest
for stretchers



EN90003

END-T
10G Instruments table
for stretchers



CB09012

OXYPACK
Oxygen tank carrier
for stretchers – 2 L



IF01047

TRACK 4
Telescopic IV pole



ST00043

S-WINGS
Enlargment platform
for stretchers

Rev.0 (10/06/2025)
UNCERTIFIED COPY

©Copyright - This document is property of SPENCER ITALIA S.r.l and can't be reproduced, used and disclosed without the written approval of SPENCER ITALIA SRL.