



EMERGENCY BACKPACK

LIFEBAGS L**MAIN FEATURES**

- Backpack can be fully opened via double zip
- Large front pocket for additional material and smaller front pocket with organized compartments for small items
- 2 side pockets for quick access to small items
- A MOLLE system on both sides to accommodate additional pockets or equipment
- Reinforced top handle and removable shoulder straps attached to stable D-rings with a buckle
- Reflective elements on three sides for optimal visibility even in low-light conditions
- 5 different-colored inner bags included.



LIFE BAGS L is an emergency backpack designed to ensure excellent space organization, quick access to tools and documents, ergonomic comfort, and maximum visibility in any situation. The main compartment equipped with dividers, elastic loops and Velcro® strips, ensures organised and functional management of objects and additional material, and allows easy attachment of inner bags.

The removable shoulder straps provide comfort and balanced load distribution, while the reflective elements ensure visibility even in low-light conditions. On the front, a 23 x 5 cm Velcro® surface provides space for attaching identification patches, while a small transparent front pocket enables fast and visible labeling.



TECHNICAL DATA

Codes	0406605	0406606
Size	L	L
Dimensions	37 x 47 x 22 cm	37 x 47 x 22 cm
Weight without accessories	2,27 kg	2,27 kg
Inner bags green & orange	5 x 30 x 11 cm	5 x 30 x 11 cm
Inner bags yellow & red	8 x 22 x 10 cm	8 x 22 x 10 cm
Inner bag blue	8 x 22 x 22 cm	8 x 22 x 22 cm
Color		
Material	PES with PVC coating	PES with PVC coating

Note: Dimensions and weight are subject to tolerances of ± 2 cm and ± 0.2 kg respectively.



THE PRODUCT INCLUDES

- 1 green inner bag
- 1 orange inner bag
- 1 red inner bag
- 1 yellow inner bag
- 1 blue inner bag.



Rev.0 (16/07/2025)
UNCERTIFIED COPY

©Copyright - This document is property of SPENCER ITALIA S.r.l and can't be reproduced, used and disclosed without the written approval of SPENCER ITALIA SRL.

Spencer Italia s.r.l. a socio unico soggetta a direzione e coordinamento da Protect Medical Holding GmbH
Via Provinciale, n° 12 - 43038 Sala Baganza (PR) Italy - +39 0521 541111 - info@spencer.it - www.spencer.it