



EMERGENCY BACKPACK

LIFEBAGS XL**MAIN FEATURES**

- Backpack can be fully opened via double zip
- Two large front pockets: the lower section is equipped with elastic loops, while the upper section features Velcro® surfaces for attaching inner bags
- Three side pockets for quick access to small items
- Large main compartment with multiple loops and Velcro® strips to fix the inner pouches
- Reinforced top and side handles
- Removable shoulder straps with chest strap attached to stable D-rings with a buckle
- Reflective elements on three sides for optimal visibility even in low-light conditions
- 6 colored inner bags included
- Integrated 2 L oxygen bottle holder
- Integrated center wall with elastic loops and transparent pockets.



LIFEBAGS XL is a spacious emergency backpack designed to ensure excellent space organization, quick access to tools and documents, ergonomic comfort, and maximum visibility in any situation.

A side-access document compartment allows for quick retrieval of important papers without opening the main compartment. The main compartment equipped with elastic loops and Velcro® strips, ensures organised and functional management of objects and additional material, and allows easy attachment of inner bags. The central divider is equipped with numerous elastic loops in various sizes and transparent compartments to help organize tools and supplies effectively.

The removable shoulder straps with chest strap provide comfort and balanced load distribution, while the reflective elements ensure visibility even in low-light conditions. On the lower front pocket, a 26 x 5 cm Velcro® surface provides space for attaching identification patches, while a small transparent front pocket enables fast and visible labeling.

TECHNICAL DATA

Codes	0406610	0406611
Size	XL	XL
Dimensions	44 x 52 x 25 cm	44 x 52 x 25 cm
Weight without accessories	4,7 kg	4,7 kg
Inside pockets green & orange	5 x 30 x 11 cm	5 x 30 x 11 cm
Inside pockets yellow & red	8 x 22 x 10 cm	8 x 22 x 10 cm
Inside pocket blue	8 x 22 x 22 cm	8 x 22 x 22 cm
Color		
Material	PES with PVC coating	PES with PVC coating

Note: Dimensions and weight are subject to tolerances of ± 2 cm and ± 0.2 kg respectively.



THE PRODUCT INCLUDES

- 1 green inner bag
- 1 orange inner bag
- 1 red inner bag
- 1 yellow inner bag
- 2 blue inner bags.



Rev.0 (16/07/2025)
UNCERTIFIED COPY

©Copyright - This document is property of SPENCER ITALIA S.r.l and can't be reproduced, used and disclosed without the written approval of SPENCER ITALIA SRL.

Spencer Italia s.r.l. a socio unico soggetta a direzione e coordinamento da Protect Medical Holding GmbH
Via Provinciale, n° 12 - 43038 Sala Baganza (PR) Italy - +39 0521 541111 - info@spencer.it - www.spencer.it