



EMERGENCY BAG



LIFEBAGS M**MAIN FEATURES**

- Made of durable and waterproof material, available in two colours
- Three side pockets for quick access to small items
- A flat transparent pocket on top
- Large main compartment with multiple loops and internal dividers for smaller sections
- An additional zipped pocket over the main compartment
- MOLLE system on the side pockets for attaching additional pouches or equipment
- Velcro® surface on the back (32 x 4 cm) for attaching identification patches
- Small transparent front insert for quick labeling
- Reflective elements around the entire perimeter ensure visibility in low-light environments
- Removable shoulder strap and two handles with Velcro® closure for comfortable, versatile transport.



LIFEBAGS M is a spacious, water-repellent emergency bag specifically designed to hold a large amount of rescue equipment, available in different colours. The large main compartment is divided into smaller sections by multiple internal dividers. The front pocket features a small transparent insert for easy labelling, while the two side pockets are equipped with the MOLLE system for attaching additional pouches or equipment. Inside the lid is a document compartment with elastic loops to secure smaller accessories. On the outside of the lid is a transparent zippered pocket. The M version includes an additional zipped pocket above the main compartment.

TECHNICAL DATA

Codes	0405050	0405051
Dimensions	49 x 30 x 28 cm	49 x 30 x 28 cm
Adjustable shoulder strap	103-133 cm	103-133 cm
Weight without accessories	1,6 kg	1,6 kg
Color		
Material	PES with PVC coating	PES with PVC coating

Note: Dimensions and weight are subject to tolerances of ± 2 cm and ± 0.2 kg respectively.

THE PRODUCT INCLUDES

- Removable shoulder strap.

Product technical sheet **LIFEBAGS S** – Codes **0405053** - **0405054** – **LIFEBAGS M** – Codes **0405050** - **0405051**

Rev.0 (02/10/2025)
UNCERTIFIED COPY

©Copyright - This document is property of SPENCER ITALIA S.r.l and can't be reproduced, used and disclosed without the written approval of SPENCER ITALIA SRL.

Spencer Italia s.r.l. a socio unico soggetta a direzione e coordinamento da Protect Medical Holding GmbH
Via Provinciale, n° 12 - 43038 Sala Baganza (PR) Italy - +39 0521 541111 - info@spencer.it - www.spencer.it