



COMPACT UNIVERSAL HEAD IMMOBILIZER

SUPER

MAIN FEATURES

- Cushions material helps to absorb bumps.
- Large hoop and loop surface on two sides offer optimal cushions adherence
- Forehead and chin straps are adjustable and padded with an high-adhesion fabric to optimize the stabilization of the head.
- Two holes on the sides allow for ear inspection.
- The padding of the base increases patient comfort.
- Pillows are made of non-absorbent material for an easy sanitation.
- Super allows for preliminary diagnostic examinations with X-rays.



Class I Medical Device compliant with MDR (EU) 2017/745

SUPER is a compact head immobilizer that can be used with spine boards. Made of high density foam, it is waterproof and provides high levels of durability and resistance. It doesn't contain ferromagnetic materials that are known to interfere with the resonance environment.

TECHNICAL DATA

Code	SH00110	SH00111
Pillow dimensions (1 piece)	25,5 x 8,5 x 16 cm 10 x 3.34 x 6.3 in	25,5 x 8,5 x 16 cm 10 x 3.34 x 6.3 in
Base dimensions (Buttonhole excluded)	36,5 x 26 x 2 cm 14.3 x 10.2 x 0.78 in	36,5 x 26 x 2 cm 14.3 x 10.2 x 0.78 in
Hole for ear inspection	Ø 7 cm / 2.75 in	Ø 7 cm / 2.75 in
Base weight	360 g / 0.79 lb	360 g / 0.79 lb
Pillow weight	500 g / 1.10 lb	500 g / 1.10 lb
Overall weight	1,46 kg / 3.21 lb	1,46 kg / 3.21 lb
Materials	PE, PVC, Nylon	PE, PVC, Nylon
Color	●	●

Product technical sheet **SUPER** – CND Classification Y060312
Code **SH00110** – Registration number 2418457 • Code **SH00111** – Registration number 161554

Rev.0 (28/11/2025)
UNCERTIFIED COPY

©Copyright - This document is property of SPENCER ITALIA S.r.l and can't be reproduced, used and disclosed without the written approval of SPENCER ITALIA SRL.