

4BELL STAIR+ MAX

10G Fixation System



Chair mounting system designed for installation inside emergency vehicles, intended for use with the 4BELL STAIR+ chair. The device must be securely fixed to the vehicle's wall, with careful consideration of the vehicle's structural design and the positioning of its onboard equipment.

Main features

The one-piece structure allows for easy installation and cleaning.
The chair can be secured effortlessly.

Technical Data

Dimensions	490x143x780 ± 5mm
Materiali	Al, PP, Acciaio

Reference standard

EN 1789

Medical Device class I in compliance with EU Reg. 2017/745

4BELL STAIR+ MAX SISTEMA DI
FIXATION SYSTEM 10G
ST10455B



CND Classification

V9099

Rev.0 (01/09/2022)
UNCHECKED COPY

Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it