




SPINAL BOARD WITH PINS
ROCK PIN

MAIN FEATURES

- Rock Pin's robust and lightweight structure ensures impact resistance
- It can float and used for water rescue
- The number, positioning and shape of the handles simplify rescue operations
- Tapered shape is storage-friendly
- Holes and pins offer multiple strapping and fixation points for spider straps and head immobilizers
- The surface in polyethylene is easy to clean
- Radiolucent spinal board allowing diagnostic imaging, including X-ray (XR), MRI, and CT scans.
- With the Spencer Mark Up system you can tag the spinal board and minimize the risk of losing it.



Class I Medical Device compliant with MDR (EU) 2017/745

ROCK PIN is a spinal board designed to immobilize and carry patients with suspected spinal injuries. Made of a single piece of high-density polyethylene injected with expanded polyurethane, it offers great resistance, is waterproof and suitable for water rescue. It can carry up to 350 kg / 772 lb. It is compatible with any standard spider straps system or head immobilizer.





772lb
350 kg
LOAD CAPACITY

7,3 kg / 16 lb
DEVICE WEIGHT

MARK-UP SYSTEM

EASY TO STORE

TECHNICAL DATA

Codes	ST02010	ST02012	ST02013	ST02014
Length	186 cm / 73.2 in	186 cm / 73.2 in	186 cm / 73.2 in	186 cm / 73.2 in
Width	44,5 cm / 17.5 in	44,5 cm / 17.5 in	44,5 cm / 17.5 in	44,5 cm / 17.5 in
Thickness	5 cm / 1.97 in	5 cm / 1.97 in	5 cm / 1.97 in	5 cm / 1.97 in
Weight	7,3 ± 0,3 kg 16.1 ± 0.7 lb	7,3 ± 0,3 kg 16.1 ± 0.7 lb	7,3 ± 0,3 kg 16.1 ± 0.7 lb	7,3 ± 0,3 kg 16.1 ± 0.7 lb
Handles	16	16	16	16
Pins	14	14	14	14
Load capacity	350 kg / 772 lb	350 kg / 772 lb	350 kg / 772 lb	350 kg / 772 lb
Material	PE	PE	PE	PE
Color				



RECOMMENDED ANCHORING SYSTEMS




ST02105

FIX BOARD
 10G vertical fastener
 for spinal board



ST02108

FIX BOARD
 10G horizontal fastener
 for spinal board




ACCESSORIES




ST02015

PIN STRAPS
Strap system with hooks



ST02022

ECS STRAPS
Spinal board/mattress strap system



ST02035

ROCK STRAPS
Universal spinal board belt system



ST02039

T-STRAPS
Universal spinal board strap system



SH00201

SPENCER CONTOUR
Anatomic universal head immobilizer



SH00203

SPENCER CONTOUR HP
Anatomic universal head immobilizer



SH00260

PEDI FIX
Paediatric head immobilizer



SH00110

SUPER
Blue compact universal head immobilizer



SH00162

SUPER BLUE HP
Universal head immobilizer

ACCESSORIES



ST02106

MARK UP
Permanent
personalization system



SH00240

TANGO FIX
Adult/paediatric head
immobilizer



SH00246

TANGO FIX HP
Adult/paediatric head
immobilizer



SH00151

FXA PRO
Anatomic universal
ankle immobilizer



ST00494

STX 494
Belt with metal buckle
with carabiner – 2
pieces



ST02102

SPINE PACK
Empty spinal board
transport bag