




SPINAL BOARD WITH PINS  
**ROCK PIN**

**MAIN FEATURES**

- Rock Pin's robust and lightweight structure ensures impact resistance
- It can float and used for water rescue
- The number, positioning and shape of the handles simplify rescue operations
- Tapered shape is storage-friendly
- Holes and pins offer multiple strapping and fixation points for spider straps and head immobilizers
- The surface in polyethylene is easy to clean
- Radiolucent spinal board allowing diagnostic imaging, including X-ray (XR), MRI, and CT scans.
- With the Spencer Mark Up system you can tag the spinal board and minimize the risk of losing it.



Class I Medical Device compliant with MDR (EU) 2017/745

**ROCK PIN** is a spinal board designed to immobilize and carry patients with suspected spinal injuries. Made of a single piece of high-density polyethylene injected with expanded polyurethane, it offers great resistance, is waterproof and suitable for water rescue. It can carry up to 350 kg / 772 lb. It is compatible with any standard spider straps system or head immobilizer.





772lb  
350 kg  
**LOAD CAPACITY**

6,7 kg / 14,7 lb  
**DEVICE WEIGHT**

**MARK-UP SYSTEM**

**EASY TO STORE**

**TECHNICAL DATA**

Codes	ST02010	ST02012	ST02013	ST02014
Length	186 cm / 73.2 in	186 cm / 73.2 in	186 cm / 73.2 in	186 cm / 73.2 in
Width	44,5 cm / 17.5 in	44,5 cm / 17.5 in	44,5 cm / 17.5 in	44,5 cm / 17.5 in
Thickness	5 cm / 1.97 in	5 cm / 1.97 in	5 cm / 1.97 in	5 cm / 1.97 in
Weight	6,7 ± 0,3 kg 14,7 ± 0.7 lb	6,7 ± 0,3 kg 14,7 ± 0.7 lb	6,7 ± 0,3 kg 14,7 ± 0.7 lb	6,7 ± 0,3 kg 14,7 ± 0.7 lb
Handles	16	16	16	16
Pins	14	14	14	14
Load capacity	350 kg / 772 lb	350 kg / 772 lb	350 kg / 772 lb	350 kg / 772 lb
Material	PE	PE	PE	PE
Color				



**RECOMMENDED ANCHORING SYSTEMS**




ST02105

**FIX BOARD**  
 10G vertical fastener  
 for spinal board



ST02108

**FIX BOARD**  
 10G horizontal fastener  
 for spinal board




**ACCESSORIES**




**ST02015**

**PIN STRAPS**  
Strap system with hooks



**ST02022**

**ECS STRAPS**  
Spinal board/mattress strap system



**ST02035**

**ROCK STRAPS**  
Universal spinal board belt system



**ST02039**

**T-STRAPS**  
Universal spinal board strap system



**SH00201**

**SPENCER CONTOUR**  
Anatomic universal head immobilizer



**SH00203**

**SPENCER CONTOUR HP**  
Anatomic universal head immobilizer



**SH00260**

**PEDI FIX**  
Paediatric head immobilizer



**SH00110**

**SUPER**  
Blue compact universal head immobilizer



**SH00162**

**SUPER BLUE HP**  
Universal head immobilizer

## ACCESSORIES



ST02106

**MARK UP**  
Permanent  
personalization system



SH00240

**TANGO FIX**  
Adult/paediatric head  
immobilizer



SH00246

**TANGO FIX HP**  
Adult/paediatric head  
immobilizer



SH00151

**FXA PRO**  
Anatomic universal  
ankle immobilizer



ST00494

**STX 494**  
Belt with metal buckle  
with carabiner – 2  
pieces



ST02102

**SPINE PACK**  
Empty spinal board  
transport bag