



SPINE BOARD WITH PINS

B-BAK PIN MAX

MAIN FEATURES

- B-Bak Pin is only 45 mm thick and very easy to store
- With a weight of about 6.7 kg, it can support up to 454 kg
- 14 ergonomic handholes offer a high number of gripping points.
- 8 molded pins around the perimeter offer additional fixation points for spider straps
- Large oval space at the bottom facilitate immobilization of the lower limbs without interfering with the handles
- X-ray translucent
- If you use the Spencer Mark Up system you can tag the B-Bak Pin and not lose it.



Class I Medical Device compliant with MDR (EU) 2017/745

B-BAK PIN MAX is a lightweight, yet sturdy spine board designed for immobilizing patients with suspected spinal injuries. Weighing only 6,7 kg, it can support up to 454 kg (1000 lb) and features numerous handles and pins for secure grip and accessory attachment.

Its design facilitates easy transport and storage, even in compact spaces like ambulance stretchers. Additionally, its X-ray translucent nature ensures it can be used without interfering with diagnostic imaging, making it an essential tool for pre-hospital and in-hospital emergency care.




1000 lb
↔
454 kg
LOAD CAPACITY

MARK-UP SYSTEM

EASY TO STORE

EASY TO CLEAN

TECHNICAL DATA

Codes	ST02067	ST02047	ST02068
Length	184 cm	184 cm	184 cm
Width	40,5 cm	40,5 cm	40,5 cm
Thickness	4,5 cm	4,5 cm	4,5 cm
Handholes	14	14	14
Molded pins	8	8	8
Weight	6.7 kg ± 0,3kg	6.7 kg ± 0,3kg	6.7 kg ± 0,3kg
Load capacity	454 kg	454 kg	454 kg
Material	PE	PE	PE
Color			



RECOMMENDED WALL MOUNT SYSTEMS



ST02105A

FIX BOARD
 10G vertical fastener
 for spine board



ST02108A

FIX BOARD
 10G horizontal fastener
 for spine board




ACCESSORIES



ST02015A

PIN STRAPS
 Strap system with hooks



ST02022A

ECS STRAPS
 Spine board/mattress strap system



ST02035A

ROCK STRAPS
 Universal spine board belt system



ST02039A

T-STRAPS
 Universal spine board strap system



SH00201A

SPENCER CONTOUR
 Anatomic universal head immobilizer



SH00203A

SPENCER CONTOUR HP
 Anatomic universal head immobilizer



SH00162C

SUPER BLUE HP
 Universal head immobilizer



SH00110C

SUPER
 Blue compact universal head immobilizer



SH00240A

TANGO FIX
 Adult/paediatric head immobilizer

ACCESSORIES

SH00260A

PEDI FIX
Paediatric head
immobilizer

SH00246A

TANGO FIX HP
Adult/paediatric head
immobilizer

SH00151A

FXA PRO
Anatomic universal
ankle immobilizer

ST02106A

MARK UP
Permanent
personalization
system

ST00494B

STX 494
Belt with metal buckle
with carabiner – 2
pieces

ST02102B

SPINE PACK
Empty spine board
transport bag