




SPINE BOARD WITH PINS
B-BAK PIN

MAIN FEATURES

- B-Bak Pin is only 4,5 cm thick and very easy to store
- With a weight of only 6,7 kg, it can support up to 180 kg
- 14 ergonomic handholes offer a high number of gripping points
- 8 molded pins around the perimeter offer additional fixation points for spider straps
- Large oval space at the bottom facilitate immobilization of the lower limbs without interfering with the handles
- X-ray translucent
- If you use the Spencer Mark Up system you can tag the B-Bak Pin and not lose it.






Class I Medical Device compliant with MDR (EU) 2017/745

B-BAK PIN is a lightweight, yet sturdy spine board designed for immobilizing patients with suspected spinal injuries. It features numerous handles and pins for secure grip and accessory attachment. Weighing about 6,7 kg, it can support up to 180 kg and features numerous handles and pins for secure grip and accessory attachment.

Its design facilitates easy transport and storage, even in compact spaces like ambulance stretchers. Additionally, its X-ray translucent nature ensures it can be used without interfering with diagnostic imaging, making it an essential tool for pre-hospital and in-hospital emergency care.



TECHNICAL DATA

Codes	ST02061	ST02046	ST02069
Length	184 cm	184 cm	184 cm
Width	40,5 cm	40,5 cm	40,5 cm
Thickness	4,5 cm	4,5 cm	4,5 cm
Handholes	14	14	14
Molded pins	8	8	8
Weight	6.7 kg ± 0,3kg	6.7 kg ± 0,3kg	6.7 kg ± 0,3kg
Load capacity	180 kg	180 kg	180 kg
Material	PE	PE	PE
Color			



RECOMMENDED WALL MOUNT SYSTEMS



ST02105A

FIX BOARD
 10G vertical fastener
 for spine board



ST02108A

FIX BOARD
 10G horizontal fastener
 for spine board




ACCESSORIES



ST02015A

PIN STRAPS

Strap system with hooks



ST02022A

ECS STRAPS

Spine board/mattress strap system



ST02035A

ROCK STRAPS

Universal spine board belt system



ST02039A

T-STRAPS

Universal spine board strap system



SH00201A

SPENCER CONTOUR

Anatomic universal head immobilizer



SH00203A

SPENCER CONTOUR HP

Anatomic universal head immobilizer



SH00162C

SUPER BLUE HP

Universal head immobilizer



SH00110C

SUPER

Blue compact universal head immobilizer



SH00240A

TANGO FIX

Adult/paediatric head immobilizer

ACCESSORIES



SH00260A

PEDI FIX
Paediatric head
immobilizer



SH00246A

TANGO FIX HP
Adult/paediatric head
immobilizer



SH00151A

FXA PRO
Anatomic universal
ankle immobilizer



ST02106A

MARK UP
Permanent
personalization
system



ST00494B

STX 494
Belt with metal buckle
with carabiner – 2
pieces



ST02102B

SPINE PACK
Empty spine board
transport bag