




MILITARY SPINE BOARD WITH PINS

B-BAK PIN

MAIN FEATURES

- B-Bak Pin is only 45 mm thick and very easy to store
- It can support up to 180 kg (in the Max version up to 454 kg)
- 14 ergonomic handholes offer a high number of gripping points.
- 8 molded pins around the perimeter offer additional fixation points for spider straps
- Large oval space at the bottom facilitate immobilization of the lower limbs without interfering with the handles
- X-ray translucent
- If you use the Spencer Mark Up system you can tag the B-Bak Pin and not lose it.



Class I Medical Device compliant with MDR (EU) 2017/745

B-BAK PIN is a lightweight, yet sturdy spine board designed for immobilizing patients with suspected spinal injuries. Weighing about 6,7 kg, it can support up to 180 kg / 397 lb (in the Max version up to 454 kg / 1000 lb) and features numerous handles and pins for secure gripping and attachment of accessories.

Its design facilitates transport and storage, even in compact spaces such as ambulance stretchers. In addition, its X-ray translucency ensures that it can be used without interfering with imaging, making it an essential tool for pre-hospital and hospital care.



MARK-UP
SYSTEM



EASY TO
STORE



EASY TO
CLEAN

TECHNICAL DATA

Codes	ST02069	ST02068
Models	B-BAK PIN	B-BAK PIN MAX
Length	184 cm	184 cm
Width	40,5 cm	40,5 cm
Thickness	4,5 cm	4,5 cm
Handholes	14	14
Molded pins	8	8
Weight	6.7 kg ± 0,3kg	6.7 kg ± 0,3kg
Load capacity	180 kg	454 kg
Material	PE	PE
Color	●	●



RECOMMENDED WALL MOUNT SYSTEMS



ST02105A

FIX BOARD
 10G vertical fastener
 for spine board



ST02108A


FIX BOARD
 10G horizontal fastener
 for spine board

ACCESSORIES




ST02015A

PIN STRAPS
 Strap system with hooks



ST02022A

ECS STRAPS
 Spine board/mattress strap system



ST02035A

ROCK STRAPS
 Universal spine board belt system



ST02039A

T-STRAPS
 Universal spine board strap system



SH00201A

SPENCER CONTOUR
 Anatomic universal head immobilizer



SH00203A

SPENCER CONTOUR HP
 Anatomic universal head immobilizer



SH00162C

SUPER BLUE HP
 Universal head immobilizer



SH00110C

SUPER
 Blue compact universal head immobilizer



SH00240A

TANGO FIX
 Adult/paediatric head immobilizer

ACCESSORIES



SH00260A

PEDI FIX
Paediatric head
immobilizer



SH00246A

TANGO FIX HP
Adult/paediatric head
immobilizer



SH00151A

FXA PRO
Anatomic universal
ankle immobilizer



ST02106A

MARK UP
Permanent
personalization
system



ST00494B

STX 494
Belt with metal buckle
with carabiner – 2
pieces



ST02102B

SPINE PACK
Empty spine board
transport bag