



BIVALVE RIGID CERVICAL COLLAR

NEC LOC**MAIN FEATURES**

- Ergonomic two-piece design with a bivalve structure and securing strap to maintain the collar position during application.
- Dedicated sizes for adult and pediatric patients, developed according to patient morphology, facilitate correct size selection.
- Easy and quick to apply.
- Featuring a tracheostomy opening and a shaped, comfortable padding.
- Radiolucent for X-ray and CT imaging.
- MRI compatible under specific conditions.
- Also available in kits (3–6 sizes), including a carrying bag.



Class I medical device compliant with EU Regulation 2017/745 (MDR)

INTENDED USE: NEC LOC is a rigid bivalve cervical collar designed to help maintain neutral alignment of the cervical spine during first aid and emergency situations. Its use helps reduce the risk of paralysis and further spinal injury in trauma patients during transport.

PRODUCT DESCRIPTION: The NEC LOC bivalve cervical collar is available in kits and in different sizes for both adult and pediatric patients. The sizes are designed according to patient morphology, facilitating correct size selection.

PEDIATRIC SIZES**TECHNICAL DATA**

Codes	JM10007	JM10008	JM10009	JM10021
Sizes	1 (Pediatric, 0–2 years)	2 (Pediatric, 2–6 years)	3 (6–12 years)	4 (XS)
Material	PE/Plastizote	PE/Plastizote	PE/Plastizote	PE/Plastizote
Color				

NEC LOC Bivalve Cervical Collar – Product Technical Data Sheet: Code **JM10007** SIZE 1 (0-2 YEARS) • Code **JM10008** SIZE2 (2-6 YEARS)
 Code **JM10009** SIZE3 (6-12 YEARS) • Code **JM10021** SIZE 4 (XS) • Code **JM10030** SIZE 5 (STOUT) • Code **JM10040** SIZE 6 (SMALL)
 Code **JM10050** SIZE 7 (MEDIUM) • Code **JM10060** SIZE 8 (LARGE)
NEC LOC 3-Size Collar kit (4-6-7) with bag – Code JM10100 • NEC LOC 6-Size Collar kit (2-4-5-6-7-8) with bag – Code JM10101

ADULT SIZES



TECHNICAL DATA

Codes	JM10030	JM10040	JM10050	JM10060
Sizes	5 (Stout)	6 (Small)	7 (Medium)	8 (Large)
Material	PE/Plastizote	PE/Plastizote	PE/Plastizote	PE/Plastizote
Color				

KIT



TECHNICAL DATA

Codes	JM10100	JM10101
Models	3-Size Collar kit	6-Size Collar kit
Sizes included in the kit	4 (XS), 6 (Small), 7 (Medium)	2 (Pediatric, 2-6 years), 4 (XS), 5 (Stout) 6 (Small), 7 (Medium), 8 (Large)
Material	PE/Plastizote	PE/Plastizote

THE KIT INCLUDES

- Carrying bag (Dimensions: 42 x 34 cm).